

DENDRITIC POLYESTER SCAFFOLDS: FUNCTIONAL AND BIOCOMPATIBLE PRECISION POLYMERS FOR DRUG DELIVERY APPLICATIONS



Sandra García-Gallego,¹ Michael Malkoch^{2*}

¹Post-doctoral Researcher, Division of Coating Technology, Fibre and Polymer Technology, KTH Royal Institute of Technology

²Associate Professor, Division of Coating Technology, Fibre and Polymer Technology, KTH Royal Institute of Technology, and Chief Executive Officer, Polymer Factory, Sweden AB

*Email: malkoch@kth.se

Introduction

Dendrimers and dendrons are highly branched polymer structures with a precise distribution of functional groups. In contrast to their linear analogs, they are monodisperse in nature and are extremely consistent from batch-to-batch. A broad range of sophisticated applications have capitalized on these dendritic materials including catalysis, optics, and nanomedicine.¹ In nanomedicine, great efforts have been devoted toward the design of selective drug delivery systems and efficient imaging probes.² Thus, dendrimers present significant advantages over linear polymers, including the following:

- They are precision nanomedicines carrying an exact number of drugs, all arising from their controlled synthesis and monodisperse nature.
- Their properties can be directly correlated to their structure and present predictable and repeatable pharmacokinetics.

- They can be tuned to present optimal properties such as biodegradability, low toxicity, and water solubility.
- They provide an increased loading capacity, achieved by encapsulation or direct conjugation of drugs with different hydrophilic/hydrophobic properties.
- The drug release profile can selectively be tuned and controlled.
- They are intrinsically highly functional and can be designed to present multipurpose actions (e.g., therapeutic, diagnosis, targeting).

Non-toxic Biodegradable bis-MPA Dendritic Scaffolds

A promising family of dendritic macromolecules is based on the 2,2-bis(methylol)propionic acid (bis-MPA) building block.³⁻⁴ These scaffolds exhibit highly valuable properties for biomedical applications including low or no *in vitro* toxicity, no specific *in vivo* organ accumulation, no immunogenic profile, and biodegradability under physiological conditions.⁵⁻⁶ Moreover, the use of both the bis-MPA monomer and dendritic macromolecules has increased due to the commercial availability of a wide range of hyperbranched polymers, dendrimers, and dendrons (Figure 1). Herein, the most advanced applications of bis-MPA based dendritic materials as drug delivery agents will be summarized.




Figure 1. Examples of bis-MPA based dendritic structures of generation 4.