

TSK-GEL® Nonporous Resin Columns for Rapid Analysis of Proteins, Peptides, and Nucleic Acids

TSK-GEL nonporous resin columns are especially suited for microscale analysis of high molecular weight proteins, polypeptides, and nucleic acids. TSK-GEL columns are available for gel filtration, ion exchange, reversed phase, hydrophobic interaction, and affinity chromatography.

Key Words

● nonporous resin ● proteins ● peptides ● nucleic acids

Proteins, peptides, and nucleic acids are successfully analyzed by HPLC with the use of microparticulate, nonporous resin (NPR) packings. These chemically and mechanically stable nonporous resins enable large, high molecular weight proteins and polypeptides to be analyzed rapidly, with high resolution and excellent recovery.

Porous resin HPLC columns have limited capability for protein and peptide analyses. Large molecules, such as proteins, diffuse slowly through the porous matrix, slowing the separation process. If proteins are trapped in the pores, diffusion basically stops and the molecule is essentially lost — making sample recovery impossible. Protein diffusion slows the kinetics of separation and can cause peak broadening, loss in resolution, and poor sample recovery.

High resolution is attained with TSK-GEL NPR columns because their phases are covalently bonded to spherical, nonporous beads of a very small diameter — 2.5µm. Since there are no pores to trap molecules or slow the separation process, resolution and recovery are maximized.

Interaction between sample and phase is limited to the surface of these resin beads, making interaction kinetics extremely fast. Fast kinetics, combined with excellent resolution and recovery, make TSK-GEL NPR columns particularly useful for both quantitative recovery of submicrogram levels of protein and microscale purification.

Applications

Especially suited for microscale analysis — generally less than 100µg of sample — TSK-GEL NPR columns are available for ion exchange, reversed phase, and hydrophobic interaction chromatography.

DEAE-NPR and SP-NPR columns are packed with anion and cation exchangers prepared by introducing diethylaminoethyl and sulfopropyl groups into the 2.5µm diameter spherical nonporous hydrophilic polymer. Both DEAE-NPR and SP-NPR columns are well suited for protein ion exchange analyses. Oligonucleotides, RNA, and DNA fragments are best separated on DEAE-NPR.

Exposed silanols of silica-based columns can exhibit secondary binding with your sample. On the reversed phase Octadecyl-NPR column, octadecyl (C18) groups are covalently bonded to a hydrophilic polymeric backbone that eliminates this secondary binding. Reversed phase analyses of high molecular weight peptides, oligonucleotides, and proteins can be performed on this column.

Butyl groups bonded to nonporous resin beads form the packing of the TSKgel Butyl-NPR column. Proteins are separated by hydrophobic interaction chromatography, using small sample loads, with high resolution and excellent recovery of enzymatic activity.

Compared to porous columns, TSK-GEL NPR columns separate proteins and peptides more rapidly — without sacrificing resolution. A mixture of five proteins is separated at least 35 minutes faster on a Butyl-NPR column than on conventional porous hydrophobic interaction columns (Figure A). A mixture of ovalbumin and trypsin inhibitor is separated five times faster on a nonporous anion exchange column than on a conventional porous anion exchange column (Figure B). By cation exchange analysis, another protein mixture is separated in less than 1/6 the time required by a porous column (Figure C). In all cases, not only is the analysis more rapid, but excellent resolution is obtained.

TSK-GEL columns also separate large DNA restriction fragments that were traditionally performed by agarose or polyacrylamide gel electrophoresis. A plasmid DNA, pBR322, is digested into fragments with the restriction enzyme BstNI. TSKgel DEAE-NPR columns separate these DNA fragments without shearing (Figure D). DNA fragments of up to 50,000 base pairs can be separated — in 30 minutes or less — in predictable elution orders with high recovery. Separation is obtained without the frequent band broadening caused by sample diffusion through porous packing. Nonporous, hydrophilic, weak anion-exchange resins also eliminate size range limitations characteristic of most porous resin ion exchange gels.

A fast, precise, quantitative assay for PCR product, which separates specific and non-specific products, is required for optimization of PCR yield. PCR analysis traditionally utilizes gel electrophoresis, fluorescence assay, and liquid- or solid-phase hybridization to a labeled probe (1). None of these techniques provides a fast, simple, specific quantitation procedure. Figure E shows efficient separation of a *Hae III* digest of pBR322 DNA obtained on a TSKgel DEAE-NPR column with a 5mm guard column.

TSK-GEL NPR columns are especially effective for micropreparative separations, accurate analyses of rare proteins, or instances where samples are very small (low nanogram levels). In these cases, optimum recovery is important. TSK-GEL NPR columns provide excellent recovery of a wide range of proteins, even though their loading capacity is very low (Table 1).

Figure A. TSK-GEL Nonporous Resin Columns Offer Higher Resolution and More Rapid Analyses than Porous Resin Columns


Mobile Phase: **A-1** (NH₄)₂SO₄, 12 min linear gradient elution, 2.3M-0 in 0.1M phosphate buffer (pH 7.0)
A-2, A-3 (NH₄)₂SO₄, 60 min linear gradient elution, 1.8M-0 in 0.1M phosphate buffer (pH 7)
 Flow Rate: 1.0mL/min
 Det.: 280nm UV
 Temp.: 25°C

- Mixture**
1. Myoglobin
 2. Ribonuclease
 3. Lysozyme
 4. α-Chymotrypsin
 5. α-Chymotrypsinogen

Hydrophobic Interaction Separation

A-1 — Butyl-NPR


Column: TSKgel Butyl-NPR, 3.5cm x 4.6mm ID, 2.5µm particles
 Cat. No.: 8-14947



794-0014

A-2 — Phenyl-5PW


Column: TSKgel Phenyl-5PW, 7.5cm x 7.5mm ID, 10µm particles
 Cat. No.: 8-07573



794-0015

A-3 — Ether-5PW

Column: TSKgel Ether-5PW, 7.5cm x 7.5mm ID, 10µm particles
 Cat. No.: 808641



794-0016

Figure B. Higher Resolution by Anion Exchange


B-1

Column: TSKgel DEAE-NPR, 3.5cm x 4.6mm ID, 2.5µm particles
 Cat. No.: 813075


B-2

Column: TSKgel DEAE-5PW, 7.5cm x 7.5mm ID, 10µm particles
 Cat. No.: 807164
 Mobile Phase: NaCl linear gradient 0-0.5M in 20mM Tris-HCl buffer (pH 8.0), **B-1** 10 min; **B-2** 60 min
 Flow Rate: **B-1** 1.5mL/min, **B-2** 1.0mL/min
 Det.: 280nm UV
 Temp.: 25°C

B-1 — DEAE-NPR



B-2 — DEAE-5PW



- Mixture**
1. Ovalbumin
 2. Trypsin inhibitor

794-0017,0018

Figure C. Cation Exchange Analysis on Nonporous Resin Saves Time


C-1

Column: TSKgel SP-NPR, 3.5cm x 4.6mm ID, 2.5µm particles
 Cat. No.: 813076


C-2

Column: TSKgel SP-5PW, 7.5cm x 7.5mm ID, 10µm particles
 Cat. No.: 807161
 Mobile Phase: NaCl linear gradient 0-0.5M in 20mM MES buffer (pH 6.0) (z[N-morpholine]ethanesulfonic acid-NaOH), **C-1** 10 min, **C-2** 60 min
 Flow Rate: **C-1** 1.5mL/min, **C-2** 1.0mL/min
 Det.: 280nm UV

C-1 — SP-NPR



C-2 — SP-5PW



- Mixture**
1. α-Chymotrypsinogen
 2. Ribonuclease

794-0019, 0020

Figure D. Effective Separation of pBR 322 DNA, Digested with BstNI




Figure E. pBR322 DNA-HAE III Digest Using DEAE-NPR Column




Table 1. Protein Recovery from TSK-GEL NPR Columns

Protein	Recovery (%)			
	DEAE-NPR [■]	SP-NPR [▲]	Octadecyl-NPR [▼]	Butyl-NPR [♦]
α-Chymotrypsin	---	100	---	91
α-Chymotrypsinogen	---	---	---	93
α-Chymotrypsinogen A	---	95	---	---
β-Lactoglobulin	101	---	---	---
Bovine serum albumin	102	---	96	---
Cytochrome C	---	93	95	---
Ferritin	99	---	---	---
δ-Globulin	104	---	---	---
Hemoglobin	91	88	---	---
Insulin	---	---	102	---
Lysozyme	---	96	104	101
Myoglobin	---	---	97	95
Ovalbumin	103	---	74	90
Ribonuclease	---	95	94	89
Thyroglobulin	100	---	---	---
Transferrin	---	---	101	---
Trypsin inhibitor	98	---	---	84
Trypsinogen	---	87	---	---

All Columns: 3.5cm x 4.6mm ID.

■ DEAE-NPR: 5µg protein, 10 min., linear gradient NaCl 0-0.5M in 20mM Tris-HCl buffer (pH 8.0), Flow Rate: 1.5mL/min., Det.: 280nm.

▲ SP-NPR: 5µg protein, 10 min., linear gradient NaCl 0-0.5M in 20mM phosphate buffer (pH 7.0), Flow Rate: 1.5mL/min., Det.: 280nm.

▼ Octadecyl-NPR: 0.5µg protein, 8 min., linear gradient acetonitrile 20-80% in 0.05% TFA, Flow Rate: 1.5mL/min., Det.: 220nm.

♦ Butyl-NPR: 2µg protein, 12 min., linear gradient (NH₄)₂SO₄ from 2.3M-0 in 0.1M phosphate buffer (pH 7.0), Flow Rate: 1.0mL/min., Det.: 280nm.

Using TSK-GEL NPR Columns

Preparing Your HPLC System for Column Installation

TSK-GEL NPR columns are provided with zero dead volume compression fittings and include nuts and ferrules. A TSK-GEL NPR column requires 1/16" fittings for connection to your HPLC system. If your system has other fittings, adapters are available from our current Supelco Chromatography Supplies Catalog.

Since most HPLC column end fittings differ from one another, we recommend you make up two short pieces of 1/16" OD x 0.010" (0.25mm) ID capillary tubing to connect the column to the injector (or better, in-line frit filter) and the detector. Use enclosed 1/16" nuts and ferrules, plus the appropriate fittings for connecting the injector and detector, to prepare these connectors. Or you can use our 1/16" handtight fittings, Cat. No. 58478 or 58462.

In either case, make sure you insert the capillary tubing fully into the column end fitting before setting the ferrule. Once the tubing is fully inserted, set the ferrule by tightening the nut. In the same manner, make connections at the other end of the 1/16" capillary tubing.

Before installing your new TSK-GEL NPR column, flush all nonwater-miscible solvents from the system. Be sure to flush injector loops, valves, and other portions of the system that could trap previously used mobile phases. Flush the entire system with at least 100mL of distilled water to remove the last traces of organic solvents, then flush with 50mL of each of the concentrated salt buffers you intend to use.

Preparing Buffers and Solvents

Prevent air bubbles from entering the column during installation, use, and storage. Air bubbles may cause degradation of column performance due to formation of channels in the packed bed. Mobile phases must be thoroughly degassed by filtration or helium sparging. Filtration will also remove small particles that could clog the column frit.

Only high quality reagents, water, and solvents should be used for preparing buffers (e.g., super distilled or HPLC grade distilled water, and special grade reagents). Resin fouling, leading to a loss in resolution and/or efficiency, occurs faster in TSK-GEL NPR columns due to the small surface area of nonporous resin particles.

Use 0.20 or 0.45µm microporous membranes to filter all buffers, salt solutions, aqueous and aqueous/organic mobile phases (e.g., a 0.45µm Nylon-66 membrane).

Be careful not to mix, or use in sequence, solutions that might precipitate or gel in the column or system. When replacing an aqueous buffer with an organic solvent, flush the column with distilled water as an intermediate step.

Column Equilibration

In ion exchange chromatography, counterions can have a significant influence on a separation. Cationic SP-NPR columns are shipped in the Na⁺ form, and anionic DEAE-NPR columns are shipped in the Cl⁻ form. Both columns are shipped in an aqueous form, and should be equilibrated by using a buffered mobile phase containing appropriate ions specific to your analyses. Equilibration typically requires ~40-50 column volumes and is indicated by a stable baseline.

Solvent changes are often needed to fine-tune a specific separation. Gradient changes made to achieve a new mobile phase composition must be done *slowly* and *gradually*. We recommend a gradient change of not more than 2%/minute and a flow rate of not more than 0.5mL/minute. After the new composition is reached, equilibrate the column with ~20-30 column volumes of mobile phase. The system will then be ready for use. To ensure a stable baseline, run the selected gradient profile once or twice without injecting the sample.

Salt Solutions

Continuous, indiscriminate use of highly concentrated salt solutions can cause instrument and analysis problems. The following precautions are recommended during analyses utilizing salt solutions.

Be sure your system has no leaks. Salt deposits at a fitting (buffer creep) reveal a leak that can cause poor reproducibility, and, eventually, instrument downtime. Wash the area thoroughly with distilled water. Remove the fitting and wash the inside with distilled water. *Do not try to compensate for the leak by overtightening the fitting, or you might damage it.*

Do not shut off the pump for more than a few minutes at a time, while the system contains concentrated salts. When you do shut down the system, turn off the pump and remove the column, then thoroughly flush the entire system with distilled water. (Note: the column should not be stored for any length of time in a concentrated salt solution. See Table 2 for recommended storage solutions.)

Periodically — or when you want to use an organic mobile phase in the system (e.g., a reversed phase analysis on another column) — completely flush the system, *but not the column*, with a 10% nitric acid solution. Then rinse it with distilled water. Refer to your instrument manuals or contact the manufacturer, however, before using nitric acid.

Sample Preparation

If possible, always dissolve your sample in mobile phase, or in the starting mobile phase when operating under gradient conditions. Try to match the pH, salt concentration, and organic solvent of the sample with those of the mobile phase and run a test to ensure that no precipitate, suspension, or flocculate is formed. Finally, before making an injection, filter the sample through a 0.20µm porosity filter.

Table 2. TSK-GEL NPR Column Operating Conditions and Specifications

Condition/Specification	DEAE-NPR	SP-NPR	Octadecyl-NPR	Butyl-NPR
Column Dimensions	3.5cm x 4.6mm	3.5cm x 4.6mm	3.5cm x 4.6mm	3.5cm x 4.6mm
Shipping and Storage Solvent	Distilled Water	Distilled Water	70% Methanol/30% Water	Distilled Water
Max. Flow Rate	1.6mL/min.	1.6mL/min.	1.6mL/min.	1.2mL/min.
Standard Flow Rate	1.0-1.5mL/min	1.0-1.5mL/min	1.0-1.5mL/min	0.5-1.0mL/min
Max Pressure	200kg/cm ² (3000psi)	200kg/cm ² (3000psi)	200kg/cm ² (3000psi)	200kg/cm ² (3000psi)
pH Range [□]	2-12	2-12	2-12	2-12
Salt Conc.	≤ 1 Molar	≤ 1 Molar	≤ 0.5 Molar	< 4 Molar
Organic Conc.	≤ 20%	≤ 20%	0-100%	≤ 20% (Confirm salt does not precipitate)
Sample Loading Capacity [▲]	10µg (proteins) 1-5µg (oligonucleotides)	10µg	0.5µg	2µg
Temperature	≥ 0°C	≥ 0°C	≤ 40°C	≥ 0°C
Resolution [▼]	≥ 6.0 (a)	≥10.0 (b)		≥ 3.0 (c)

□ pH above 12 or below 2 can be used only for a short time.

▲ Defined as the amount of pure protein that can be injected onto the column before peak broadening, due to overload, occurs. The loading capacity for crude samples can be as much as 100 times greater since it is the sum of the loading capacity of the proteins in the (crude) mixture.

▼ Resolution:
$$(RX1, X2) = \frac{2(V2-V1)}{1.7(W2 + W1)}$$

where V1, V2 = elution volumes

W1, W2 = widths of peaks X1 and X2, at half height

(a) = Resolution for trypsin inhibitor and ovalbumin

(b) = Resolution for α-trypsinogen and trypsinogen

(c) = Resolution for ovalbumin and lysozyme

Column Flow Direction

The recommended flow direction through a TSK-GEL NPR column is indicated by the arrow on the tag (FLOW →). Operating the column with flow in the reverse direction is recommended only as a last resort when removing particulates from a clogged frit, or as part of the following column installation procedure. (Refer to the *Column Instruction Manual* for the specific procedure to clean a clogged frit.)

Column Installation

Your column was tested by our quality assurance department after manufacture, then flushed with storage solvent and closed with caps to prevent solvent evaporation. During installation, use, and storage, you should take precautions to prevent air from entering the column.

Perform the following checks before using a new column:

- Remove the cap from the column *inlet* side. Solvent should be visible at the inlet fitting. If the fitting appears dry, see step (f).
- Connect the column inlet to the injector outlet tubing.
- Remove the cap from the column outlet.
- Start mobile phase flow using a slow rate (0.2-0.3mL/min).
- After liquid flows from column exit, connect the column to the detector inlet line. The column is now ready for use.

If the column inlet fitting appears dry, we recommend the following:

- Disconnect the column outlet cap. Hook the column outlet to the injector outlet tubing.
- Slowly start flow (0.1mL/min) in this reversed flow direction until a few drops of mobile phase exit the column.
- Turn off flow. Let the pressure go to zero and disconnect the column from the system.
- Turn the column around, connecting it so flow is now in the direction of arrow indicated on the tag (FLOW →).
- Now you can hook the column to the detector inlet, and slowly increase the flow rate to its desired setting.

Column Operating Conditions

Table 2 summarizes operating conditions and specifications for each TSK-GEL NPR column. Additional information can also be found in the *Operating Conditions and Specifications* sheet, and the *Column Instruction Manual*, included with each TSK-GEL column you receive.

Connecting Columns in Series

For complex samples, or those containing components with widely differing molecular weights, you can use up to TSK-GEL NPR columns in series to maximize resolution. Equilibrate all columns with the same mobile phase *before* connecting them. Operate the system within recommended flow rates for each column (Table 2) to avoid exceeding maximum back pressure. As in other types of HPLC, doubling column length (e.g., by coupling two columns in series) increases resolution by $\sqrt{2}$.

To minimize dead volume, use small bore connecting tubing. Insert it fully into the end fitting before tightening. Keep air out of the columns and connecting tubing to avoid channeling in the packing bed. Make connections while maintaining a mobile phase flow (0.3mL/minute), but be sure to follow safety precautions (hand protection, etc.) appropriate for the mobile phase you are using.

Column Protection

Guard columns are not available — and are usually not needed — for these very short TSK-GEL columns. However, we recommend using a 0.45 μ m in-line frit filter (Supelco Cat. No. 59124) to prevent particulates from plugging the column frit. For PCR applications, it is important to include a DEAE-NPR guard column.

Column Cleaning

Occasionally — usually after long-term use — elution times may differ significantly from original analyses. To remedy, follow cleaning procedures, using recommended solvents for each column (Table 3).

Column Storage

The column does not need to be removed from the system if used daily. The mobile phase may be left in the column overnight if the buffer salt has good solubility.

When storing the column, flush it with the shipping solvent, remove it from the instrument and seal the ends with the screws provided. All columns should be stored at room temperature.

During all the above steps, prevent air from entering the column.


Column Longevity

Column life depends on mobile phases, pH, sample matrix, contaminants, and other factors. You can make more than 1000 injections onto a properly protected column. Longer column life can be anticipated when you use simple, clean samples and mild mobile phases.

Table 3. Cleaning Solvents and Procedures for TSK-GEL NPR Columns

	DEAE-NPR & SP-NPR	Octadecyl-NPR	Butyl-NPR
Solvents	<ol style="list-style-type: none"> 0.1-0.2M NaOH 20-40% Acetic Acid, aqueous Urea or nonionic surfactant in buffer Aq. buffer in 20% organic solvent 	<ol style="list-style-type: none"> Acetonitrile or methanol Aq. buffer in organic solvent 0.1-0.2M NaOH 20-40% Acetic Acid, aqueous 60% Acetonitrile/40% 0.2M NaOH 	<ol style="list-style-type: none"> 0.1-0.2M NaOH 20% Acetic Acid, aqueous
Procedures	Clean column regularly by injecting up to one column volume 0.1-0.2M NaOH with <ul style="list-style-type: none"> 100-250μL increments on DEAE-NPR 250μL-2mL increments on SP-NPR 	Clean column regularly by injecting up to one column volume cleaning solution in 100-250 μ L increments.	Clean column daily after use, or regenerate column, by injecting 100-250 μ L NaOH. If repeated injections are not effective, use repeated 100-250 μ L injections of acetic acid.

Ordering Information:



910-0274

TSK-GEL NPR Ion Exchange Columns

3.5cm x 4.6mm ID

Description	Cat. No.	Exchange Mechanism	Functional Group/Counterion
SP-NPR	813076	Strong Cation	$-\text{SO}_3^-/\text{Na}^+$
DEAE-NPR	813075	Weak Anion	$-\text{CH}_2\text{CH}_2\text{N}^+(\text{Et})_2/\text{Cl}^-$

Description	Cat. No.
Guard Column	
DEAE-NPR, 0.5cm x 4.6mm, 5 μm particles	817088
Reversed Phase Column	
3.5cm x 4.6mm ID, Octadecyl-NPR	814005
Hydrophobic Interaction Column	
3.5cm x 4.6mm ID, Butyl-NPR	814947




910-0142

Column Inlet Filter

Prevents particulates from clogging column and frit. Minimizes dead volume and band broadening.

Description	Cat. No.
Frit, 3mm, 0.5 μm pores	59124
Replacement Frits, 3mm, 0.5 μm pores, pk. of 5	59126




910-0146

High Pressure Fingertight Upchurch Fittings

Use these stainless steel fingertight fittings with all common 1/16" external fittings used on HPLC columns. Will hold to 6000psi (420kg/cm²) when finger tightened, or to 10,000psi (700kg/cm²) with a slight turn of a wrench.

Description	Cat. No.
Fingertight Fittings, ea.	58478-U
Ferrules, 1/16", pk. of 5	58479




913-0340

Dynaseal Handtight Fittings

Will seal against pressures to over 7000psi (490kg/cm²). Each fittings consists of a reusable nut, collet, and ferrule. The ferrule will not deform stainless steel seats and can be reused.

The fittings are compatible with all 1/16" female compression-type fittings now used on HPLC columns. Spacers are helpful for leak-free sealing with some column end fittings where thread depth varies.

Description	Cat. No.
Handtight Fittings, pk. of 2	58462
Replacement Ferrules, pk. of 5	58463
Collets, pk. of 2	58464
Ferrules, KEL-F [®] polymer, pk. of 10	58468
Spacers, pk. of 2	58679



910-0160

Stainless Steel Capillary Tubing

Capillary tubing for use between pump and injector.
1/16" OD x 0.010" ID. One piece each.

Description	Cat. No.
5cm lengths	56713
10cm lengths	56714
20cm lengths	56715-U
30cm lengths	56716
0.5m lengths	56717
1.0m lengths	56718-U

Peek Tubing

Polyetheretherketone, 1/16" OD, 10-foot lengths.

ID (in.)	Color	Cat. No.
0.005	Red	Z227307
0.007	Yellow	Z226688
0.01	Blue	Z226661
0.02	Orange	Z227293
0.03	Green	Z226955



910-0354

Replaceable Frit and Screen Filters

Replaceable 1/8" frit has 0.5µm pores to protect 3µm or larger column packings from particles. Use with 1/16" tubing.

Description	Cat. No.
Frit Filter	58420-U
Replacement Frits (pk. of 10)	
0.5µm pores	59037
2.0µm pores	59129

Use this efficient, low dead volume in-line filter between your pump and injector to trap particles. The replaceable screen has 2µm pores. Use with 1/16" tubing.

Screen Filter	58279-U
Replacement Screen (pk. of 10)	58284

TSK-GEL Columns:

Affinity Chromatography

Column	Particle Size (µm)	Column Length (cm)	Dimensions Diameter (mm)	Cat. No.
TSKgel ABA-5PW	10	7.5	7.5	813067
TSKgel ABA-5PW guard column kit	20	1	6.0	813127
TSKgel ABA-5PW packing (5mL)	20	—	—	813128
TSKgel Boronate-5PW	10	7.5	7.5	813066
TSKgel Boronate-5PW packing (5mL)*	20	—	—	813126
TSKgel Chelate-5PW	10	7.5	7.5	808645
TSKgel Chelate-5PW guard column kit	20	1	6.0	808647
TSKgel Chelate-5PW packing (5mL)	20	—	—	808648
TSKgel Heparin-5PW	10	7.5	7.5	813064
TSKgel Heparin-5PW guard column kit	20	1	6.0	813121
TSKgel Heparin-5PW packing (5mL)	20	—	—	813122
TSKgel Tresyl-5PW	10	4	6.0	814455
TSKgel Tresyl-5PW	10	7.5	7.5	814456

*Requires a TSKgel Guardgel holder, 6mm x 1cm (803432). Replacement frits, 2µm, pk. of 10 (803430).

Gel Filtration Chromatography — Polymer-Based

Column	Particle Size (µm)	Column Length (cm)	Dimensions Diameter (mm)	Cat. No.
TSKgel G-Oligo-PW	6	30	7.8	808031
TSKgel Oligo guard column	12	4	6.0	808034
TSKgel G2500PWXL	6	30	7.8	808020
TSKgel G3000PWXL	6	30	7.8	808021
TSKgel G4000PWXL	10	30	7.8	808022
TSKgel G5000PWXL	10	30	7.8	808023
TSKgel G6000PWXL	13	30	7.8	808024
TSKgel GMPWXL	13	30	7.8	808025
TSKgel G-DNA-PW	10	30	7.8	808032
TSKgel PWXL guard column	12	4	6.0	808033
TSKgel PWXL packing (1g)	10	—	—	808035
TSKgel G1000PW**	10	30	7.5	805760
TSKgel G2000PW**	10	30	7.5	805761
TSKgel G2500PW	10	30	7.5	808028
TSKgel G3000PW	10	30	7.5	805762
TSKgel G4000PW	17	30	7.5	805763
TSKgel G5000PW	17	30	7.5	805764
TSKgel G6000PW	17	30	7.5	805765
TSKgel GMPW	17	30	7.5	808026
TSKgel PWH guard column	17	7.5	7.5	806732

**Use Oligo guard column (808034) for G1000PW and G2000PW; PWH guard column (806732) for G2500PW through GMPW.

Gel Filtration Chromatography — Silica-Based

Column	Particle Size (µm)	Column Length (cm)	Dimensions Diameter (mm)	Cat. No.
TSKgel G2000SWxL	5	30	7.8	808540
TSKgel G3000SWxL	5	30	7.8	808541
TSKgel G4000SWxL	8	30	7.8	808542
TSKgel SWxL guard column	7	4	6.0	808543
TSKgel SWxL packing (1g)	5	—	—	808544
TSKgel G2000SW	10	30	7.5	805788
TSKgel G2000SW	10	60	7.5	805102
TSKgel G3000SW	10	30	7.5	805789
TSKgel G3000SW	10	60	7.5	805103
TSKgel G4000SW	13	30	7.5	805790
TSKgel G4000SW	13	60	7.5	805104
TSKgel SW guard column	10	7.5	7.5	805371
TSKgel SW packing (1g)	10	—	—	806819
TSKgel QC-PAK GFC 200 (stainless steel)	5	15	7.8	816215
TSKgel QC-PAK GFC 200GL (glass)	5	15	8.0	816214
TSKgel QC-PAK GFC 300 (stainless steel)	5	15	7.8	816049
TSKgel QC-PAK GFC 300GL (glass)	5	15	8.0	816216

Guard column kits include 1 cm x 6.0mm column, column holder, 5mL packing, 2 frits, 2 nuts and ferrules for 1/16" tubing, 5cm of 1/16" tubing.

Hydrophobic Interaction Chromatography

Column	Particle Size (µm)	Column Length (cm)	Dimensions Diameter (mm)	Cat. No.
TSKgel Butyl-NPR	2.5	3.5	4.6	814947
TSKgel Ether-5PW	10	7.5	7.5	808641
TSKgel Ether-5PW guard column kit	20	—	—	808643
TSKgel Ether-5PW packing (5mL)	20	—	—	808644
TSKgel Phenyl-5PW	10	7.5	7.5	807573
TSKgel Phenyl-5PW guard column kit	20	—	—	807652
TSKgel Phenyl-5PW packing (5mL)	20	—	—	807651

Hydroxyapatite Chromatography

Column	Particle Size (µm)	Column Length (cm)	Dimensions Diameter (mm)	Cat. No.
TSKgel HA-1000 packing (5mL)	10	—	—	813130

Trademarks

KEL-F — 3M Co., Chemical Div.
TSK-GEL — Tosoh Corp.

Reference

1. *Amplifications*, Perkin-Elmer Corp., June 1992, Issue 8, pp 10-13, Elena Katz, Will Block, and John Wages.
Reference not available from Supelco.

Contact our Technical Service Department
(phone 800-359-3041 or 814-359-3041, FAX 814-359-5468)
for expert answers to your questions.

Ion Exchange Chromatography — Anion Exchange

Column	Particle Size (µm)	Column Length (cm)	Dimensions Diameter (mm)	Cat. No.
TSKgel DEAE-NPR	2.5	3.5	4.6	813075
TSKgel DEAE-NPR guard column	5	0.5	4.6	817088
TSKgel DEAE-5PW	10	7.5	7.5	807164
TSKgel DEAE-5PW guard column kit	20	—	—	807210
TSKgel DEAE-5PW packing (5mL)	20	—	—	807207
TSKgel DEAE-2SW	5	25	4.6	807168
TSKgel DEAE-3SW	10	7.5	7.5	807163
TSKgel DEAE-SW guard column kit	10	—	—	807648
TSKgel DEAE-SW packing (5mL)	10	—	—	807647

Ion Exchange Chromatography — Cation Exchange

Column	Particle Size (µm)	Column Length (cm)	Dimensions Diameter (mm)	Cat. No.
TSKgel CM-5PW	10	7.5	7.5	813068
TSKgel CM-5PW guard column kit	20	—	—	813069
TSKgel CM-5PW packing (5mL)	20	—	—	813070
TSKgel SP-NPR	2.5	3.5	4.6	813076
TSKgel SP-5PW	10	7.5	7.5	807161
TSKgel SP-5PW guard column kit	20	—	—	807211
TSKgel SP-5PW packing (5mL)	20	—	—	807208
TSKgel CM-2SW	5	25	4.6	807167
TSKgel CM-3SW	10	7.5	7.5	807162
TSKgel CM-SW guard column kit	10	—	—	807650
TSKgel CM-SW packing (5mL)	10	—	—	807649

Reversed Phase Chromatography

Column	Particle Size (µm)	Column Length (cm)	Dimensions Diameter (mm)	Cat. No.
Silica-Based Columns				
TSKgel Oligo-DNA RP	5	15	4.6	813352
TSKgel Oligo-DNA RP	5	15	7.8	813353
TSKgel Super-ODS	2	5	4.6	818154
TSKgel Super-ODS	2	10	4.6	818197
TSKgel G Filter** (pk. of 3)	—	—	—	818207
Resin-Based Columns				
TSKgel Octadecyl-NPR	2.5	3.5	4.6	814005
TSKgel Octadecyl-4PW	7	15	4.6	813351
TSKgel Phenyl-5PW RP	10	7.5	4.6	808043
TSKgel Phenyl-5PW RP packing (5mL)*	20	—	—	814019

Resin-Based Columns

TSKgel Octadecyl-NPR	2.5	3.5	4.6	814005
TSKgel Octadecyl-4PW	7	15	4.6	813351
TSKgel Phenyl-5PW RP	10	7.5	4.6	808043
TSKgel Phenyl-5PW RP packing (5mL)*	20	—	—	814019

*Requires a TSKgel Guardgel holder, 6mm x 1cm (**803432**); Replacement frits, 2µm, pk. of 10 (**803430**)

Requires a guard holder, 4mm x 4cm (818206**)

Guard column kits include 1 cm x 6.0mm column, column holder, 5mL packing, 2 frits, 2 nuts and ferrules for 1/16" tubing, 5cm of 1/16" tubing.

BULLETIN 862

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA · Sigma-Aldrich de Argentina, S.A. · Buenos Aires 1119 AUSTRALIA · Sigma-Aldrich Pty. Ltd. · Castle Hill NSW 2154 AUSTRIA · Sigma-Aldrich Handels GmbH · A-1110 Wien
BELGIUM · Sigma-Aldrich N.V./S.A. · B-2880 Bornem BRAZIL · Sigma-Aldrich Quimica Brasil Ltda. · 01239-010 São Paulo, SP CANADA · Sigma-Aldrich Canada, Ltd. · 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC · Sigma-Aldrich s.r.o. · 186 00 Praha 8 DENMARK · Sigma-Aldrich Denmark A/S · DK-2665 Vallensbaek Strand FINLAND · Sigma-Aldrich Finland/YA-Kemia Oy · FIN-00700 Helsinki
FRANCE · Sigma-Aldrich Chimie · 38297 Saint-Quentin-Fallavier Cedex GERMANY · Sigma-Aldrich Chemie GmbH · D-82041 Deisenhofen GREECE · Sigma-Aldrich (o.m.) Ltd. · Ilioupoli 16346, Athens
HUNGARY · Sigma-Aldrich Kft. · H-1067 Budapest INDIA · Sigma-Aldrich Co. · Bangalore 560 048 IRELAND · Sigma-Aldrich Ireland Ltd. · Dublin 24 ISRAEL · Sigma Israel Chemicals Ltd. · Rehovot 76100
ITALY · Sigma-Aldrich s.r.l. · 20151 Milano JAPAN · Sigma-Aldrich Japan K.K. · Chuo-ku, Tokyo 103 KOREA · Sigma-Aldrich Korea · Seoul MALAYSIA · Sigma-Aldrich (M) Sdn. Bhd. · 58200 Kuala Lumpur
MEXICO · Sigma-Aldrich Quimica S.A. de C.V. · 50200 Toluca NETHERLANDS · Sigma-Aldrich Chemie BV · 3330 AA Zwijndrecht NORWAY · Sigma-Aldrich Norway · Torshov · N-0401 Oslo
POLAND · Sigma-Aldrich Sp. z o.o. · 61-663 Poznań PORTUGAL · Sigma-Aldrich Quimica, S.A. · Sintra 2710 RUSSIA · Sigma-Aldrich Russia · Moscow 103062 SINGAPORE · Sigma-Aldrich Pte. Ltd.
SOUTH AFRICA · Sigma-Aldrich (pty) Ltd. · Jet Park 1459 SPAIN · Sigma-Aldrich Quimica, S.A. · 28100 Alcobendas, Madrid SWEDEN · Sigma-Aldrich Sweden AB · 135 70 Stockholm
SWITZERLAND · Supelco · CH-9471 Buchs UNITED KINGDOM · Sigma-Aldrich Company Ltd. · Poole, Dorset BH12 4QH
UNITED STATES · Supelco · Supelco Park · Bellefonte, PA 16823-0048 · Phone 800-247-6628 or 814-359-3441 · Fax 800-447-3044 or 814-359-3044 · email: supelco@sial.com

H

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.