


# CHOZN® & UCOE® Combined Platform

Accelerate Biopharmaceutical Development with an  
Integrated Upstream Solution



## Isolate more high-producing clones with fewer resources

The CHOZN® & UCOE® Combined Platform's streamlined processes offer numerous efficiencies for cell line development. While maintaining all the benefits of the CHOZN® GS expression system, through the CHOZN® & UCOE® Combined Platform you can also:

- Accelerate your cell line development by increased efficiency of isolating high producing clones. Eight-fold more clone candidate pools and twice the number of clones expressing more than 3 g/L were obtained with the combined platform.
- Increase cell line development success with hard-to-produce molecules


- Achieve high and stable titers to expedite your cGMP manufacturing with fewer resources


For nearly a decade, the CHOZN® GS platform has been used to develop high expressing, stable biomanufacturing clones. Consisting of a ZFN-modified CHO-K1 line, paired cell culture media and feeds, an expression vector, documentation and protocols, the platform is optimized to consistently isolate the high-expressing clones within a transfected pool. However, cell line development can be a lengthy and resource-intensive task because these high-expressing clones are relatively infrequent within the pool.

The UCOE® technology offers a path to increase the frequency of these high-expressing clones, thereby reducing the time and resources needed for cell line development.

The UCOE® (Ubiquitous Chromatin Opening Elements) technology uses naturally-occurring sequences capable of creating and maintaining open chromatin to enable high, stable transgene expression, regardless of the integration site. When included in a plasmid alongside the transgene of interest, these UCOE® elements greatly increase the number of high expressing clones that can be isolated within a single cell line development process.

To demonstrate the ability of the CHOZN® & UCOE® Combined Platform to increase the number of high-expressing clones, two concurrent IgG1 cell line development procedures were initiated. CHOZN® GS<sup>-/-</sup> cells were transfected with either the standard CHOZN® GS expression vector containing gene cassettes






**Figure 1.** Standard CHOZN<sup>®</sup> expression vector (left) and CHOZN<sup>®</sup> expression vector with the addition of UCOE<sup>®</sup> sequences (right).

for an IgG1 heavy and light chain (Control), or the CHOZN<sup>®</sup> expression vector with the addition of UCOE<sup>®</sup> sequences (UCOE), pictured in **Figure 1**. The CHOZN<sup>®</sup> GS<sup>-/-</sup> protocols were followed. At each screening step, the UCOE conditions yielded multi-fold more high producing pools and clones (Figure 2). This increased efficiency in generating high producing pools and clones in the CHOZN<sup>®</sup> & UCOE<sup>®</sup> Combined Platform can be leveraged to:

- Reduce screening and scale-up efforts over the course of cell line development without compromising productivity goals
- Help ensure cell line development success with hard-to-produce molecules
- Identify more pools and clones that match both titer and protein quality criteria



**Figure 2.** At each screening step, the UCOE conditions yielded multi-fold more high producing pools and clones.

## CHOZN® GS<sup>-/-</sup> Cell line

- ECACC CHO K1 origin, adapted to suspension growth in chemically defined, animal component-free media
- GS<sup>-/-</sup> phenotype knock out by targeted mutagenesis using the non-virus ZFN gene editing technology
- Enables rapid cell line development by metabolic selection without the addition of MSX (GS inhibitor)
- cGMP-banked and fully traceable cell lines
- Paired, batch-to-batch consistent and chemically defined cell culture media and feeds using highest grade of critical raw materials

## UCOE® Technology

- Highly-characterized Ubiquitous Chromatin Opening Elements create and maintain open chromatin
- Offers major improvements in gene expression in stably transfected mammalian cells

## CHOZN® & UCOE® Combined Platform

- Clear IP path for the complete platform
- Complete documentation and cell line history
- Comprehensive user protocols guide customers from transfection through small scale bioreactors
- Optimized UCOE® GS vector for expression of antibodies, Fc-fusion or recombinant proteins
- Simple and accessible combined platform licensing structure
- Technical support and consultation

To place an order or receive technical assistance in the U.S. and Canada, call toll-free 1-800-645-5476  
For other countries across Europe and the world, please visit: [EMDMillipore.com/offices](https://www.emdmillipore.com/offices)  
For Technical Service, please visit: [EMDMillipore.com/techservice](https://www.emdmillipore.com/techservice)

MilliporeSigma  
400 Summit Drive  
Burlington, MA 01803

[EMDMillipore.com](https://www.emdmillipore.com)

© 2020 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. MilliporeSigma, the vibrant M, Millipore, SAFC, CHOZN and UCOE are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

