User Guide

Millex[®]-VV/GV/HV Filter Unit (Sterile) with Durapore[®] Membrane

REF	SLVV033RS	SLGV033RS*	SLHV033RS
		SLGV033RB*	SLHV033RB


Available languages

English	3	Lietuviškai (Lithuanian).	33
Français (French).	6	Magyar (Hungarian).	36
Italiano (Italian)	9	Český (Czech)	39
Deutsch (German)	12	Polski (Polish)	42
Español (Spanish)	15	日本語 (Japanese).	45
Português (Portuguese)	18	繁體中文 (Traditional Chinese).	47
Nederlands (Dutch)	21	Русский (Russian)	49
Dansk (Danish)	24	Norsk (Norwegian).	52
Svenska (Swedish).	27	Slovenščine (Slovenian)	55
Latviski (Latvian).	30		

*Brazil Product Registration
Number 81167700001.
 See page 2.

PR05131w
 Rev 06/23
 Date of issue: 2023-06

Definition of Symbols Used

Symbol	Definition	Symbol	Definition
	Sterilized using irradiation		Temperature limit
	Non-pyrogenic		Use-by date
	Consult instructions for use		Date of manufacture
	Do not use if package is damaged and consult instructions for use		Manufacturer
	Do not re-sterilize		Caution
	Do not re-use		CE conformity marking
	Catalogue number		Swiss authorized representative
	Batch code		Unique device identifier

Made in Ireland
 Merck Millipore Ltd.
 Tullagreen, Carrigtwohill,
 Co. Cork, Ireland
 T45 KD29
 an affiliate of Merck KGaA,
 Darmstadt, Germany

MDSS CH GmbH
 Laurenzenvorstadt 61
 5000 Aarau
 Switzerland

Brazilian Registration Holder (BRH)
 Merck S/A
 Rua Torre Eifel, nº 100
 Lote 3ª, Gleba A
 Parque Rincão
 Cotia/SP
 an affiliate of Merck KGaA,
 Darmstadt, Germany
 *Brazil Product Registration
 Number 81167700001

Millex®-VV/GV/HV Filter Unit (Sterile) with Durapore® Membrane

Millex®-VV Filter	SLVV033RS	Qty: 50 pack
Millex®-GV Filter	SLGV033RS	Qty: 50 pack
	SLGV033RB	Qty: 250 pack
Millex®-HV Filter	SLHV033RS	Qty: 50 pack
	SLHV033RB	Qty: 250 pack

- 33 mm
- Sterile
- Single use only
- Non-pyrogenic

Intended Purpose

The Millex®-VV/GV/HV Filter Unit (Sterile) are medical devices intended for use as a syringe filter to sterilize and/or clarify low volume solutions in direct patient care and pharmacy admixture applications. Devices are intended for medical and/or laboratory professional use only.

Introduction

This document provides compatibility information, operating steps, and specifications for the Durapore® hydrophilic polyvinylidene fluoride (PVDF) family of sterile filters. The bidirectional support of this filter membrane enables users to filter aqueous solutions in either direction; forward (from the syringe into the container) or backward (from the container into the syringe). These filters remove microorganisms, particles, precipitates, and undissolved powders larger than the membrane's rated pore size. These single-use filters are non-pyrogenic and non-toxic. They consist of a membrane filter sealed in a modified acrylic copolymer (MMA) housing.

Applications

Typical laboratory applications include the sterile filtration (VV/GV) and/or clarification (VV/GV/HV) of protein solutions, tissue culture media, additives, buffers, and water. Clinical benefits include the removal of microorganisms, particles, precipitates, and undissolved powders larger than the membrane's rated pore size from solutions for clinical use. In the hospital pharmacy, filters with Durapore® membranes can be used for the sterile filtration (VV/GV) and/or clarification (VV/GV/HV) of small volumes of protein pharmaceuticals, diagnostic imaging agents, chemotherapeutics, aqueous solutions, or water during admixture preparation. Direct patient care applications include sterilization (VV/GV) and/or particulate removal (VV/GV/HV) from epidural and other liquid anaesthetics as well as from irrigation solutions used in ophthalmic, otic, and other surgical procedures.

Chemical Compatibility




These filters are compatible with most aqueous solutions. Compatibility of solution should be tested before use. Visit [SigmaAldrich.com/MedMillex](https://www.sigmaaldrich.com/medmillex) for more information.

WARNINGS

- Filter solutions with a temperature limit of 45 °C (113 °F).
- Syringe filters are for manual use only; do not use on automated systems.
- Not suitable for filtering high viscosity solutions or blood.
- Do not filter solutions containing 5 mg or less of active drug materials unless binding studies have been performed.
- To ensure sterility, do not use this product if the package is damaged.
- Use aseptic technique to maintain sterility removing product from primary package.
- Do not use with syringes smaller than 10 mL because pressures in excess of the maximum pressure rating may be reached, potentially causing damage to the filter and/or personal injury.
- Do not use for direct neuraxial applications which may come in contact with cerebrospinal fluid.
- Do not use the same syringe filter to filter solutions in both directions.
- Sudden loss of pressure could indicate failure of the filter.
- Single use only; do not re-use or re-sterilize.
- Make sure to wet the filter membrane thoroughly before injecting the solution; improperly wetted filters can become airlocked.
- Not suitable for use with medical pressure infusion equipment or gravity fed intravenous (IV) administration sets.
- Discard appropriately after single use. See "Disposal" below.

Directions for Use

1. Wearing gloves using aseptic technique: fill syringe with solution to be filtered.
2. Remove cover from package.
3. Attach suitable syringe to filter inlet and remove assembly from package.
4. Hold syringe with filter pointing up and top off by pushing a few drops through to wet filter.
⚠ Excess solution may be hazardous and should be disposed with care.
5. Depress syringe plunger to filter solution into suitable container or compatible attachment on outlet.

Specifications and Performance Characteristics

Materials

Membrane	Hydrophilic Durapore® polyvinylidene fluoride (PVDF)
Pore size	Millex®-VV filter: 0.10 µm Millex®-GV filter: 0.22 µm Millex®-HV filter: 0.45 µm
Housing	Modified acrylic copolymer (MMA)

Dimensions

Inlet to outlet	27 mm (1.06 in.)
Diameter	33 mm (1.30 in.)
Filtration area	4.52 cm ² (0.70 in ²)

Temperature limit	45 °C (113 °F)
--------------------------	----------------

Pressure limit at 21 °C	10.3 bar (150 psi) inlet maximum
--------------------------------	----------------------------------

Filtration volume	10 mL to 100 mL
--------------------------	-----------------

Hold-up volume	≤0.1 mL after air purge
-----------------------	-------------------------

Sterilization method	Gamma irradiation
-----------------------------	-------------------

Connections	Female Luer-Lok™ inlet; male Luer-slip outlet
--------------------	---

Flow rate at 2.1 bar (30 psi), 21 °C	Millex®-VV filter: ≥ 25 mL/min Millex®-GV filter: ≥ 40 mL/min Millex®-HV filter: ≥ 300 mL/min
---	---

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://www.sigmaaldrich.com/offices).

Disposal

Follow precautions for disposal of items contaminated with potentially infectious or hazardous biological material according to all applicable international, federal, state, and local regulations.

Technical Assistance

Visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Any serious incident of this device should be reported to manufacturer and competent authority of country where user is established.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at [SigmaAldrich.com/terms](https://www.sigmaaldrich.com/terms).

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MilliporeSigma, Millipore, Millex, Durapore, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.