

Product Information

PCR Low Ladder Marker Set

Product Code **D 7808**

Store at 2 to 8 °C

Product Description

The PCR Low Ladder Marker Set has been especially designed for size determination of PCR generated DNA fragments. Mixing the two markers together will produce a band every 20 bp, with every fifth band (100 bp increments) more intense. These markers can be used with either agarose or polyacrylamide gels.

Fragment Sizes : base pairs (bp)

20 bp Ladder

1,000	700	400	100
980	680	380	80
960	660	360	60
940	640	340	40
920	620	320	20
900	600	300	
880	580	280	
860	560	260	
840	540	240	
820	520	220	
800	500	200	
780	480	180	
760	460	160	
740	440	140	
720	420	120	

100 bp Ladder

1,000
900
800
700
600
500
400
300
200
100

Components

Suitable for 50 applications of the 20 bp ladder and approx. 150 applications of the 100 bp ladder

- 25 µg vial of 20 bp ladder, Product Code P 2849, containing 50 bands ranging from 20 to 1,000 bp in exact 20 bp increments. (Approx. 200 µg/ml)
- 25 µg vial of 100 bp ladder, Product Code P 9364, containing 10 bands ranging from 100 to 1,000 bp in exact 100 bp increments. (Approx. 100 µg/ml)
- 1 ml vial of 6X Sample Loading Buffer, Product Code B 6924, containing 15% (w/v) Ficoll, 0.06% (w/v) bromophenol blue, 0.06% (w/v) xylene cyanol FF, 30 mM EDTA.

Storage/Stability

Markers are stable for 12 months when stored at 2-8 °C. This product is shipped at ambient temperature. Routinely used marker should be stored at 2-8 °C. For infrequently used markers, storage at -20 °C may increase shelf life.

Storage Buffer for Ladder Components

10 mM Tris-HCl, pH 7.5-8.0
 1.0 mM EDTA

Product Summary

Suitable for use as electrophoresis marker for DNA. Typically, 1.5-2.0 µl of DNA ladder component, 2.5 µl of 6X loading buffer plus 10.5 µl of water is used. This translates to 15-20 ng/band, giving sharp, easily visible and photographed bands. Adjustments may be made for different well sizes and individual preferences.

Suitability Assay

Each DNA Ladder was prepared for electrophoresis as follows:

20 bp DNA Ladder

- 2.5 μ l 20 bp DNA Ladder
- 5.8 μ l water
- 1.7 μ l 6X Loading Buffer (Product Code B 6924)

100 bp DNA ladder

- 1.6 μ l 100 bp DNA Ladder
- 6.7 μ l water
- 1.7 μ l 6X Loading Buffer (Product No. B 6924)

10 μ l of each of the above solutions were loaded on a 20 cm, 3% agarose gel (agarose, wide range:routine, 3:1, Product Code A 7431). The gel was run with the appropriate DNA standards at 4 °C in 1X TBE (Product Code T 9525) for approximately 16 hours, until the bromophenol blue tracking dye reaches the bottom of the gel. **Note:** The 20 bp band of the 20 bp ladder may run off the end of the gel. After staining for 15-20 minutes in 5 μ g/ml ethidium bromide and destaining with water for 15-20 minutes, the resulting banding pattern for each ladder was consistent with the indicated DNA sizes.

Factors Influencing Resolution and Accurate Sizing of DNA Fragments

1. DNA loading: For accurate size determinations, load the smallest practical amount of sample DNA. A single application of 10-20 ng is readily visible and will yield a sharp, accurately sized band.
2. Well thickness: For best resolution of DNA bands, use only properly formed sample wells, ≤ 1 mm in thickness.
3. Salt concentration: It is important to accurately match the salt concentration of the ladder mixture to that of the DNA being evaluated in order to obtain the best size determinations. One useful technique for very precise sizing of sample-fragments, which eliminates concerns over matching salt concentrations, is to co-electrophorese the sample and the ladder in the same well. Ladder-only and sample-only lanes should be run to aid in interpretation of electrophoresis patterns.
4. Anomalies: The 100 bp ladder may show a double- or triple-banding pattern in some types of polyacrylamide gels, particularly under higher run temperatures.

7/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.