

Multiplex detection of four Alzheimer's Disease biomarkers in human CSF using a novel MILLIPLEX® assay panel

Introduction

Progressive neurodegenerative disorders, such as Alzheimer's Disease (AD), affect millions worldwide and are becoming more prevalent as our population ages. Two key neuropathological features that exemplify AD are extracellular Amyloid β ($A\beta$) plaques and intracellular neurofibrillary tangles, which are composed of the abnormally hyperphosphorylated protein Tau. Biochemical changes in $A\beta$ and Tau reflect AD pathologic processes in the brain.

Monitoring $A\beta$, Tau and other protein biomarkers in cerebrospinal fluid (CSF) of patients with these neurological disorders may be highly beneficial to understanding the pathological processes involved. Previous studies in which several protein biomarkers were monitored using a multiplexed panel showed

that analyzing combinations of two or more CSF biomarkers (such as phosphorylated Tau [Thr181] in combination with $A\beta$ 1-42) provided a more accurate diagnosis of AD than any single CSF biomarker^{1,2}.

We developed a new multiplexed assay panel, the MILLIPLEX® MAP Human Amyloid Beta and Tau Magnetic Bead Panel (Cat. No. HNABTMAG-68K), to simultaneously quantitate levels of $A\beta$ 1-40, $A\beta$ 1-42, total Tau and phosphorylated Tau (Thr181) in human CSF. Compared to other commercially available panels that do not include $A\beta$ 1-40, this 4-plex panel enables determination of the $A\beta$ 1-40: $A\beta$ 1-42 ratio, which has also been suggested as crucial for AD pathogenesis³. In fact, many mutations associated with familial AD increase the $A\beta$ 1-40: $A\beta$ 1-42 ratio.

Materials and Methods

Human CSF samples from undiagnosed (normal) and AD individuals were acquired from Discovery™ Life Sciences and Precision Medicine, respectively.

The multiplex assays for AD biomarkers were performed in 96-well plates according to product instructions supplied for the MILLIPLEX® MAP Human Amyloid Beta and Tau Magnetic Bead Panel (Cat. No. HNABTMAG-68K).

Results

The specificity of this kit was tested and showed no cross-reactivity nor any significant difference in analyte concentrations, regardless of whether each analyte was measured in a single-plex or multiplex (full 4-analyte

kit) format (Table 1). Furthermore, sensitivity, intra- and inter-assay precision, linearity of dilution and spike recovery exhibited excellent analytical performance (Table 2).

Plex	4-plex	1-plex	4-plex	1-plex	4-plex	1-plex	4-plex	1-plex
Analyte	A β 40	A β 40	A β 42	A β 42	tTau	tTau	pTau (T181)	pTau (T181)
CSF1	1993	2060	189	172	843	805	21	20
CSF1	2136	2094	193	190	867	886	21	22
CSF2	3050	3090	400	392	915	898	28	29
CSF2	2900	3068	412	381	857	898	29	25
CSF3	1160	1178	95	93	547	508	29	28
CSF3	1240	1191	102	100	544	573	30	31
CSF4	2174	2208	260	252	919	859	27	27
CSF4	2170	2284	254	249	929	879	28	28

Table 1.

Comparison of CSF sample concentration in single-plex vs. multiplex. Measurement of human CSF sample value in single-plex bead and detection vs. 4-plex beads and detection in duplicates.

	A β 1-40	A β 1-42	Total Tau	Phosphorylated Tau (Thr181)
Intra-assay (%CV)	3.8	2.6	3.9	2.9
Inter-assay (%CV)	7.9	9.0	14.2	4.8
Spike recovery in buffer (%)	93	91	92	95
Spike recovery in CSF (%)	101	109	98	91
CSF dilution linearity (%)	96	101	105	112
Cross-reactivity (%)	A β 40	< 2	0	< 0.8
	A β 42	0.4	100	0
	Total Tau	0	0	100
	Phosphorylated Tau (Thr181)	0	0	0
MinDC +2 Std Dev	10.2	5.4	14.2	1.5

Table 2.

Assay performance characteristics. Inter-assay %CVs and intra-assay %CVs were calculated using positive controls from different assays. Percent spike recovery was calculated from spiked recombinant proteins in assay buffer and CSF samples. The average percent linearity was calculated from three diluted human CSF samples. Percent cross-reactivity was performed using single, purified recombinant proteins with multiplexed beads and detection cocktails. Minimum detectable concentration (MinDC) was calculated using MILLIPLEX® Analyst 5.1 Software (Cat. No. 40-087).

Given the physicochemical characteristics of A β peptides, we investigated whether the presence of A β peptides would interfere with the detection or quantitation of the other protein biomarkers when measured simultaneously. We performed an interference test (Figure 1) in which

increasing concentrations of A β 40 or A β 42 recombinant proteins were spiked into human CSF samples. The measured concentrations of the analytes did not change, indicating that there was no detectable interference from the A β peptides.




Figure 1.

MILLIPLEX® MAP Human Amyloid Beta and Tau Magnetic Bead Panel (Cat. No. HNABTMAG-68K) interference test. Increasing concentrations of Aβ40 or Aβ42 recombinant proteins were spiked into human CSF samples (n=4). No significant differences in sample concentrations were detected.




Figure 2.

CSF sample values plotted on standard curves using the MILLIPLEX® MAP Human Amyloid Beta and Tau Panel. Aβ40, Aβ42, total Tau and phosphorylated Tau (Thr181) concentrations from human CSF samples are plotted on the standard curves using the MILLIPLEX® MAP Human Amyloid Beta and Tau Magnetic Bead Panel (Cat. No. HNABTMAG-68K).

Results (continued)

Using the MILLIPLEX® MAP multiplex immunoassay kit, A β 1-40, A β 1-42, total Tau and phosphorylated Tau (Thr181) were detected simultaneously in age-matched AD and non-AD human CSF samples (Figure 3), minimizing the use of valuable CSF samples (12.5 μ L/well). Phosphorylated Tau (Thr181) and A β 1-42 correlated significantly with the presence of an AD diagnosis (p-value < 0.01 and p-value < 0.05,

respectively). Neither total Tau nor A β 1-40 showed a significant correlation with AD diagnosis in the CSF samples used in this study.

The ratio of phosphorylated Tau (Thr181) over A β 1-42 improved the p-value < 0.001. This result was consistent with published data indicating that the ratio of phosphorylated Tau (Thr181) to A β 1-42 correlates with cognitive decline⁴.




Figure 3.

Multiplex measurement of human non-AD-CSF and AD-CSF samples. Non-AD-CSF samples (n=16) and AD CSF samples (n=16) were analyzed using the MILLIPLEX® MAP Human Amyloid Beta and Tau Magnetic Bead Panel (Cat. No. HNABTMAG-68K).

Conclusion

The MILLIPLEX® MAP Human Amyloid Beta and Tau Magnetic Bead Panel (Cat. No. HNABTMAG-68K) is ideal for simultaneous measurement of A β 1-40, A β 1-42, total Tau and phosphorylated Tau (Thr181) for neurodegenerative disease biomarkers in human CSF. This assay is rapid, sensitive, reproducible and requires only 12.5 μ L or less of CSF sample. It is important to note that the presence of either A β 1-40 or A β 1-42 does not interfere with the detection or quantitation of other biomarkers. The ability to determine both the A β 1-40:A β 1-42 ratio as well as the A β 1-42:phosphorylated Tau ratio indicates that this assay is a unique and powerful tool for studying the complexities of the nervous system and the mechanisms of pathology in AD and other neurodegenerative disorders.

References

1. Blennow K, Hampel H, Weiner M, Zetterberg H. Cerebrospinal fluid and plasma biomarkers in Alzheimer disease. *Nat Rev Neurol*. 2010 Mar;6(3):131-44.
2. Herskovits AZ, Locascio JJ, Peskind ER, Li G, Hyman BT. A Luminex assay detects amyloid β oligomers in Alzheimer's disease cerebrospinal fluid. *PLoS One*. 2013 Jul 2;8(7):e67898.
3. Dolev I, Fogel H, Milshtein H, Berdichevsky Y, Lipstein N, Brose N, Gazit N, Slutsky I. Spike bursts increase amyloid- β 40/42 ratio by inducing a presenilin-1 conformational change. *Nat Neurosci*. 2013 May;16(5):587-95.
4. Bombois S, Duhamel A, Salleron J, Deramecourt V, Mackowiak MA, Deken V, Sergeant N, Pasquier F, Buée L, Sablonnière B, Schraen-Maschke S. A new decision tree combining Abeta 1-42 and p-Tau levels in Alzheimer's diagnosis. *Curr Alzheimer Res*. 2013 May 1;10(4):357-64.

Ordering Information

MILLIPLEX® MAP Kits	Cat. No.
Human Amyloid Beta and Tau Panel	HNABTMAG-68K
Human Neuroscience Panel 1	HNS1MAG-95K
Human Neurodegenerative Disease Panel 1	HNDG1MAG-36K
Human Neurodegenerative Disease Panel 2	HNDG2MAG-36K
Human Neurodegenerative Disease Panel 3	HNDG3MAG-36K
Human Neurodegenerative Disease Panel 4	HNDG4MAG-36K
Human Neurological Disorders Panel 3	HND3MAG-39K
Human Neuropeptide Panel	HNP MAG-35K
Human Circadian Stress	HNC SMAG-35K

To Place an Order or Receive Technical Assistance

In the U.S. and Canada, call toll-free
1(800)-645-5476

For other countries across Europe, call
+44 (0) 115 943 0840

For other countries across Europe
and the world, please visit
EMDMillipore.com/offices

For Technical Service, please visit
EMDMillipore.com/techservice

MilliporeSigma
290 Concord Road
Billerica, MA 01821

