

montage™

life science kits



GENOMICS

Montage Plasmid Miniprep₉₆ Kit

Suitable for use with Qiagen® 9600 & 3000 Robotics Systems

User Guide

Notice

The information in this document is subject to change without notice and should not be construed as a commitment by Millipore Corporation or an affiliated corporation. Millipore Corporation or an affiliated corporation assumes no responsibility for any errors that may appear in this document. This manual is believed to be complete and accurate at the time of publication. In no event shall Millipore Corporation or an affiliated corporation be liable for incidental or consequential damages in connection with or arising from the use of this manual.

© 2001 MILLIPORE CORPORATION. PRINTED IN THE UNITED STATES OF AMERICA. ALL RIGHTS RESERVED. THIS BOOK OR PARTS THEREOF MAY NOT BE REPRODUCED IN ANY FORM WITHOUT THE WRITTEN PERMISSION OF THE PUBLISHERS.

P36382 Rev. A, 09/01

Millipore is a registered trademark of Millipore Corporation.

MultiScreen and Montage are trademarks of Millipore Corporation.

Big Dye is a registered trademark of the Perkin-Elmer Corporation or an affiliated company.


MEGABACE is a registered trademark of Amersham Pharmacia Biotech AB.

Qiagen is a registered trademark of Qiagen GMBH Corporation.

Contents

Introduction	1
Automation with Montage Plasmid Miniprep ₉₆ Kit	1
Kit Components	2
Additional Equipment Required	3
Precautions	4
Storage Conditions	4
Procedures for Plasmid DNA Miniprep	5
Bacterial Host Cultures	5
Culture Media and Antibiotics	6
Bacterial Cultures Growth	7
Reagents Preparation	7
Qiagen® Robot (9600 or 3000) Setup	8
Overview of Procedure	9
Standard Protocol	10

Product Performance	13
DNA Yields	14
DNA Sequence Data	16
Troubleshooting	17
References	19
Ordering Information	20
Technical Assistance	22
Standard Warranty	23




Introduction

The Montage Plasmid Miniprep₉₆ Kit provides all of the reagents and materials necessary to purify plasmid DNA using a simple protocol that eliminates lengthy bind/elute methods and centrifugation. Employing unique separation technology, Millipore has developed a line of easy-to-use DNA plasmid miniprep kits that yield DNA suitable for the most sensitive downstream applications. In addition, this new technology has significantly reduced the time required for processing samples. Following bacterial lysis, three short filtration steps are all that is required to prepare 96 clean DNA samples from each plate. The plasmid DNA is retained by Millipore's proprietary size-exclusion membrane while proteins and contaminants are filtered through to waste. Adhere to the protocols provided with this kit to ensure the generation of DNA samples that can be used for applications such as cloning, DNA sequencing, transformation, and PCR.

Automation with Montage Plasmid Miniprep₉₆ Kits

Millipore has enhanced the Montage Plasmid Miniprep₉₆ Kit's compatibility with automated systems by making it possible to collect purified samples from the top side of the plates in the kits and by eliminating lengthy centrifugation steps. This user guide provides information for automation of Montage Plasmid Miniprep₉₆ Kits on Qiagen® robotics systems.

Kit Components



Letter	Part	Function
A	96-well culture block	Growth of host bacteria
B	Montage PLASMID ₉₆ plate	Purification of plasmid DNA
C	Lysate Clearing plate	Clearing of bacterial lysate
D	V-bottom storage plate	Storage of plasmid DNA samples
E	Sealing tape	Sealing of plasmid DNA samples

Letter	Part	Function
F	RNase A	Required additive for Solution 1
G	Solution 1	Cell resuspension
H	Solution 2	Cell lysis
I	Solution 3	Neutralization
J	Solution 4	Wash
K	Solution 5	Resuspension/storage of plasmid DNA

Additional Equipment Required

- Pipettor
- Vacuum manifold (Millipore Cat. No. MAVM 096 OR)
- Centrifuge (for deep well culture block)
- Incubator shaker

Information about obtaining the following items is available through the Millipore web site: www.millipore.com/automation.


- Millipore manifold adapter and tubing assembly
- Software program for plasmid prep application

Precautions

- Montage plates are disposable, single-use-only devices.
- This kit is for research use only. Not for use in clinical applications.
- Avoid contact with Solution 2 to prevent skin irritation.
- Montage plates consist of a polystyrene plate that is sealed to a polyethylene underdrain, forming 96 independent (individually sealed) wells. **Do not separate the underdrain from the polystyrene plate. Separation will result in plate failure and well leakage.**

Storage Conditions

The kit reagents should be stored at 15 °C to 30 °C. However, following the addition of RNase A to Solution 1, this solution must be stored at 2 °C to 8 °C and should be used within six months (not to exceed expiration date on the kit).



Procedures for Plasmid DNA Miniprep

The protocols described below include a modified version of a common alkaline lysis method for isolation of plasmid DNA from bacteria¹. All of the necessary reagents are provided with this kit, including the final storage buffer (Solution 5).

Bacterial Host Cultures

Choosing a bacterial host strain is a critical part of plasmid preparation. Endonuclease I is a 12kDa protein encoded by the *endA* gene of *E. coli*. This protein degrades double-stranded DNA and can reduce the stability of plasmid preparations if contaminating levels of this endonuclease are substantial⁴. Many *E. coli* strains carry a mutation in the *endA* gene that inactivates Endonuclease I. These strains are referred to as EndA negative (EndA⁻) and are preferable hosts when preparing plasmid DNA. The table below lists some common host strains.

Table 1. EndA⁻ and EndA⁺ strains of E. coli .

EndA ⁻	EndA ⁺
DH1	BL21(DE3)
DH5a	CJ236
JM109	HB101
SRB	JM101
XL1-Blue	Q358
XLO	TB1

Culture Media and Antibiotics

Cultures grown for plasmid preparation should always be picked from a single colony taken from a plate containing a selective agent. Antibiotics should be used at every stage of growth³, because many plasmids do not contain the *par* locus that ensures equal segregation of plasmids during cell division. In addition, the stability of many antibiotics decreases during culturing. Therefore, exceeding the prescribed culturing times may result in decreased DNA yields and/or reduce the quality of the purified samples. High copy number plasmids such as pUC, pBluescript, and pGEM² are recommended for use with this kit in order to achieve optimal plasmid yields.

The recommended medium for culturing bacteria prior to purification of plasmid DNA using the Montage Plasmid Miniprep₉₆ Kits is 2x Luria-Bertani (Miller) broth. This culture medium has twice the content of tryptone and yeast extract compared to normal LB (Miller) broth but has the same salt content. The formulation per liter of media is as follows:

Tryptone	20g
Yeast Extract	10g
NaCl	10g

Richer media such as Terrific Broth or Super Broth can be used for culturing to increase cell densities and plasmid yields. However, the kit performance has been optimized using the 2x LB described above and the use of richer media may increase filtration times and/or elevate protein contaminant levels and adversely affect downstream applications.

Bacterial Cultures Growth

Bacterial cultures must be started the day before the intended run of the robot.

Inoculate *E. coli* host into 1 mL aliquots of 2X LB plus antibiotic in sterile 96 deep well blocks (2.2 mL capacity). Cover blocks with the lids provided and secure in incubator. Incubate at 37 °C at 320 rpm for 20–24 hours.

Reagents Preparation

- Add RNase (total contents of tube) to Solution 1, mix thoroughly, and store at 2 °C to 8 °C. (All other solutions should be stored at 15 °C to 30 °C.)
- Bring Solution 2 to room temperature before using in order to dissolve detergent that may have precipitated due to lower temperatures during shipping.
- Screw cap on Solution 2 tightly immediately after use in order to avoid destabilization that may occur on exposure to air.

CAUTION: Avoid skin contact with Solution 2 to prevent irritation.






Qiagen® Robot (9600 or 3000) Deck Setup

The protocol for using the robotics system to purify plasmid DNA starts with the bacterial cell pellet obtained after the culture has been grown and centrifuged, and the supernatant discarded. Set up the robot as follows:

1. Place the PLASMID plate inside the Millipore vacuum manifold.
2. Replace top ring.
3. Place Lysate Clearing plate on the top ring of the manifold.

The robot is now ready for processing the bacterial cultures started the day before.

Overview of Procedure

<p>1</p> <p>Deep well culture block</p>  <ol style="list-style-type: none">1. Culture bacteria.2. Centrifuge.3. Resuspend cells.4. Lyse cells.5. Neutralize.	<p>2</p> <p>Clearing plate</p> <p>PLASMID plate</p>  <ol style="list-style-type: none">6. Place PLASMID plate inside manifold.7. Transfer lysates into Clearing plate.8. Reassemble manifold with Clearing plate on top of manifold. Apply vacuum for 3 minutes.	
<p>3</p> <p>PLASMID plate</p>  <ol style="list-style-type: none">9. Discard Clearing plate. Transfer PLASMID plate to top of manifold. Apply vacuum for 8 minutes.10. Add 200 μL of Solution 4 to each well. Apply vacuum for 3–5 minutes.	<p>4</p> <p>PLASMID plate</p>  <ol style="list-style-type: none">11. Add 50 μL of Solution 5 to each well. Shake for 5 minutes to resuspend plasmids.	<p>5</p> <p>Plate sealing tape</p> <p>V-bottom plate</p>  <ol style="list-style-type: none">12. Transfer samples to V-bottom plate for use or storage.

Standard Protocol

1. Centrifuge the deep well blocks containing the bacterial cultures at $1500 \times g$ for 5–7 minutes. After centrifugation, immediately decant culture supernatant to a container for proper disposal. Invert and tap the plates firmly on several layers of paper towels to remove residual culture supernatant.

NOTE: Failure to remove media will add undesired volume to lysate.

2. Place deep well block on the front left plate shaker position of the Qiagen robot deck. Resuspend pellets by adding 150 μL of Solution 1 to each well then mixing on the plate shaker until cells are totally resuspended.

NOTE: Thorough resuspension of cells is critical for successful lysis.

3. Add 150 μL of Solution 2 to each well and mix immediately and vigorously on the plate shaker.
4. Add 150 μL of Solution 3 to each well and mix immediately and vigorously on the plate shaker. At this point, the bacterial lysate is ready for transfer to the Lysate Clearing plate.

5. Transfer 200 μL of lysate from the deep well block to the Lysate Clearing plate on top of the Millipore manifold. Apply vacuum at 270 mbar for 3 minutes.

NOTE: For processing of entire lysate volume, see “Alternative Protocol” below.

Standard Protocol, continued

6. Discard the Clearing plate and place the PLASMID plate on top of the empty manifold. Apply vacuum at 815 mbar for 8 minutes or until wells are empty. Direct filtrate to waste.

NOTE: Filtration time is sample, temperature, and pressure dependent. The filters will appear shiny even after the wells are empty.

7. Add 200 μL of Solution 4 to each well of the PLASMID plate. Apply vacuum at 815 mbar for 3–5 minutes, or until wells are empty. Direct filtrate to waste.
8. Move the PLASMID plate to front right shaker position. Place V-bottom collection plate on back right shaker position.
9. Recover plasmid by adding 50 μL of Solution 5 to each well of the PLASMID plate. Shake for 5 minutes.
10. Pipette retained plasmid from the wells of the PLASMID plate into the V-bottom plate for storage. Use the sealing tape to seal wells of the V-bottom storage plate. Alternatively, samples may be stored in sealed PLASMID plates for several weeks in the refrigerator.


Alternative Protocol

An alternative protocol may be used to maximize plasmid DNA yields from the same bacterial cultures described above. This protocol makes use of the full lysate generated from the entire 1 mL culture and plasmid DNA yields may reach up to 10 µg per well. However, filtration times may have to be extended for optimal performance, and there may be a decrease in the uniformity of the yields across the plate. See Table 2 in the "DNA Yields" section for typical yields obtained with this alternative protocol.

The following two modifications need to be made to the standard protocol:

- Use 100 µL (instead of 150 µL) of each of Solutions 2, 3, and 4 in steps 2, 3, and 4 of the standard protocol, and
- Load the entire bacterial lysate into the Clearing plate (instead of only 200 µL) in step 5 of the standard protocol.

For all other steps, proceed as described in the standard protocol.



Product Performance

The devices and protocols included in this kit are designed to generate high quality plasmid DNA for automated platforms such as the Qiagen 9600 and 3000 robotics systems.

In the manual protocols, after the centrifugation of cultures, each 96 well block can typically be processed in less than 15 minutes to produce a bacterial lysate. The time required for clarification of the lysate is typically 1–3 minutes, while the filtration of the cleared lysate through the PLASMID plate and the ensuing wash step are generally 5 minutes and 3 minutes, respectively, when following the standard protocol. The total time required for the combined lysis and purification steps is typically less than 25 minutes. The processing time for more than one plate depends on the number of available vacuum manifolds.

For details on automated protocols, please refer to the publications shipped with the instrument. Also, log onto Millipore's web site at www.millipore.com/automation for information and products to help you set up and operate Qiagen robotics systems with the Montage Plasmid Miniprep₉₆ Kit. Contact Millipore if you would like our automation specialists to work with you to optimize your Qiagen automated platform for use with Montage Plasmid Miniprep₉₆ Kit.

DNA Yields

The DNA sample yields are dependent upon several factors, including the type of plasmid employed², the type of media used, and the final cell densities. Typical yields achieved using our protocols are shown in Table 2, below. The agarose gel analysis in Figure 1 demonstrates the purity and reproducibility of the plasmid DNA samples prepared using the Montage Plasmid Miniprep₉₆ kit.


Table 2. Typical yields of plasmid DNA using Montage Plasmid Miniprep₉₆ kits.

Host Strain	Plasmid	Protocol	Yield (µg)
XL1-blue	pUC19	Standard	2.6 ± 0.5
		Alternative	5.1 ± 1.0
	pLH2*	Standard	3.9 ± 0.4
		Alternative	8.2 ± 1.0
JM109	pUC19	Standard	2.7 ± 0.3
		Alternative	5.7 ± 0.7
	pLH2*	Standard	3.9 ± 0.4
		Alternative	8.0 ± 1.1

*The plasmid pLH2 is derived from pUC19. The 2.0 Kb fragment of HindIII digested λ phage was cloned into the HindIII site of pUC19, resulting in a 4.7 Kb plasmid.

DNA Yields, continued

Figure 1. Agarose gel showing plasmid DNA purified using the Montage Plasmid Miniprep₉₆ Kit.




DNA Sequence Data

When using fluorescent dye terminator chemistry for sequencing of plasmid DNA, the purity of the plasmids is a critical determinant of the quality and utility of the DNA sequence data, especially when attempting to miniaturize reactions on automated platforms. Table 3 demonstrates that plasmids prepared using this kit consistently yield high quality sequence data, as judged by PHRED analysis ^{5,6}.

Table 3. Sequence data from plasmids prepared with Montage Plasmid Miniprep₉₆ Kit.

Host Strain	Plasmid	Avg Read Length n=24 (>98% accuracy)	Avg PHRED 20 bases n=24
XL1-blue	PUC19	707	618
	PLH2	658	623
JM109	PUC19	706	602
	PLH2	695	543

*All plasmid DNA samples were prepared with the Standard Protocol. They were sequenced using 1/16th BigDye® sequencing reactions (Applied Biosystems) and were run on a MegaBACE® capillary array sequencer (Molecular Dynamics). 160ng of pLH2 and 100ng of pUC19 were used in the respective reactions. Reactions were cleaned up using Montage SEQ₉₆ plates.




Troubleshooting

This section outlines how to troubleshoot problems you may encounter when using the manual protocols for the Montage Plasmid Miniprep₉₆ Kit.

Problem	Possible Causes	Suggestions
Low plasmid yields	Inadequate resuspension of cell pellets	Ensure that cells are completely resuspended in Solution 1 prior to addition of Solution 2. Failure to do so will result in reduced plasmid yields.
	No antibiotic added to media	Ensure that antibiotics are added at every stage of bacterial culturing.
	Inappropriate culturing times	Adhere to culturing conditions prescribed in the protocols.
	Proliferation of non-transformed cells	Always inoculate cultures from freshly picked colonies grown on a selective plate.
	Low copy number plasmid used	Use high copy number plasmids such as pUC, pBluescript, and pGEM.


Troubleshooting, continued

Problem	Possible Causes	Suggestions
Low plasmid yields	Overdrying PLASMID plate	Stop vacuum when wells appear empty.
Plates separating	No underdrain support grid	Use Millipore manifold or equivalent with support grid.
Wells not filtering uniformly	Airlock due to bubble in the well	Agitate the lysate until bubble surfaces.
Nicking or denaturation of plasmid DNA; poor plasmid quality.	Excessive incubation at the alkaline lysis step (Solution 2).	Do not exceed the incubation time specified in the protocol and avoid excessive shaking or vortexing.



References

1. Birnboim, H.C. and Doly, J. (1979) "A rapid alkaline extraction procedure for screening recombinant plasmid DNA." *Nucl. Acids Res.* 7, 1513-1522.
2. Sambrook, J., Fritsch, E.F., Maniatis, T. eds. (1989). "Molecular Cloning, a laboratory manual," 2nd ed. Cold Spring Harbor Laboratory Press.
3. Yanisch-Perron, D., Vieira, J. and Messing, J., (1985) *Gene* 33, 103.
4. Taylor R.G., Walker D.C., and McInnes R.R. (1993) "*Nucleic Acids Res.*" 21(7):1677-8.
5. Ewing B, Green P. (1998) "Base-calling of automated sequencer traces using phred. II. Error probabilities." *Genome Res.*8 (3):186-94.
6. Ewing B, Hillier L, Wendl MC, Green P. (1998) "Base-calling of automated sequencer traces using phred. I. Accuracy assessment." *Genome Res.*8(3): 175-85.



Ordering Information

This section lists catalogue numbers for the Montage Plasmid Miniprep₉₆ Kit. See “Technical Assistance” for information about contacting Millipore. You can also buy Millipore products on-line at www.millipore.com/purecommerce.

Product	Catalogue No.	Qty/Pk
---------	---------------	--------

Kits

<i>Montage Plasmid Miniprep₉₆ Kit suitable for Qiagen 3000 and 9600 robotics systems:</i> plates, culture blocks, and reagents for 4 × 96 samples	LSKP Q96 04	4/pk
---	-------------	------

<i>Montage Plasmid Miniprep₉₆ Kit suitable for Qiagen 3000 and 9600 robotics systems:</i> plates, culture blocks, and reagents for 24 × 96 samples	LSKP Q96 24	24/pk
--	-------------	-------

Reagents

Solution 1, cell resuspension solution, 1,500 mL	LSKCR1500	1/pk
--	-----------	------

Solution 1, cell resuspension solution, 500 mL	LSKCRS500	1/pk
--	-----------	------

Solution 2, cell lysis solution, 1,500 mL	LSKCL1500	1/pk
---	-----------	------


Solution 2, cell lysis solution, 500 mL	LSKCLS500	1/pk
---	-----------	------

Solution 3, neutralization solution, 1,500 mL	LSKNS1500	1/pk
---	-----------	------

Solution 3, neutralization solution, 500 mL	LSKNS0500	1/pk
---	-----------	------

Ordering Information, continued

Product	Catalogue No.	Qty/Pk
Solution 4, wash solution, 1,500 mL	LSKNF1500	1/pk
Solution 4, wash solution, 500 mL	LSKNF0500	1/pk
Solution 5, storage solution, 500 mL	LSKCTB500	1/pk
RNase A, 30mg. (in 50% glycerol)	LSKPMRN30	1/pk
Accessories		
V-bottom plates	LSKV BP1 00	100/pk
Cell culture blocks with lids, 96 wells, 2.2 mL	LSKC CB0 50	50/pk
Adhesive plate sealing tape	LSKA ST1 00	100/pk
Millipore filtration manifold	MAVM 096 0R	1/pk
Vacuum pressure pump, 115V/60Hz	XX55 000 00	1/pk
Vacuum pressure pump, 220V/50Hz	XX55 220 50	1/pk
Vacuum pressure pump, 100V/50 or 60Hz	XX55 100 00	1/pk
Manifold adapter and tubing assembly	Contact Millipore Technical Service	
Software program for plasmid prep application	Contact Millipore Technical Service	




Technical Assistance

For more information, contact the Millipore office nearest you. In the U.S., call 1-800-MILLIPORE (1-800-645-5476). Outside the U.S., see your Millipore catalogue for the phone number of the office nearest you or go to our web site at www.millipore.com/offices for up-to-date worldwide contact information. You can also visit the tech service page on our web site at <http://www.millipore.com/techservice>.

Millipore Corporation is pleased to provide internet access to Material Safety Data Sheets (MSDS) for its products that contain hazardous materials. To obtain any MSDS documents that may be associated with this product, go to the MSDS page of our website at <http://www.millipore.com/msds.nsf/home>.

This product is a Millipore Corporation product and Millipore is the sole warrantor of the goods. This product is neither endorsed nor warranted by Qiagen GMBH Corporation.



Standard Warranty

Millipore Corporation ("Millipore") warrants its products will meet their applicable published specifications when used in accordance with their applicable instructions for a period of one year from shipment of the products. **MILLIPORE MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** The warranty provided herein and the data, specifications and descriptions of Millipore products appearing in Millipore's published catalogues and product literature may not be altered except by express written agreement signed by an officer of Millipore. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorized and if given, should not be relied upon.

In the event of a breach of the foregoing warranty, Millipore's sole obligation shall be to repair or replace, at its option, the applicable product or part thereof, provided the customer notifies Millipore promptly of any such breach. If after exercising reasonable efforts, Millipore is unable to repair or replace the product or part, then Millipore shall refund to the customer all monies paid for such applicable product or part. **MILLIPORE SHALL NOT BE LIABLE FOR CONSEQUENTIAL, INCIDENTAL, SPECIAL OR ANY OTHER INDIRECT DAMAGES RESULTING FROM ECONOMIC LOSS OR PROPERTY DAMAGE SUSTAINED BY ANY CUSTOMER FROM THE USE OF ITS PRODUCTS.**

MILLIPORE