

Data Sheet

Gamma Stable Opticap® XL Capsule Filters with Millipore Express® SPG Hydrophobic Membrane

Gamma stable sterilizing-grade gas filters with superior flow rates for single-use applications

Gamma Stable Opticap® XL Capsule Filters with Millipore Express® SPG (Sterile Phobic Gamma) hydrophobic membrane are the ultimate choice for optimal gas filtration and venting in single-use applications. The 0.2 µm sterilizing-grade membrane provides sterility assurance and high flow rates. Its ability to maintain flow rates during extended batch durations contributes to higher process confidence.



Typical Applications

Gamma Stable Opticap® XL Capsule Filters with Millipore Express® SPG Hydrophobic membrane are validated for retention using a liquid bacterial challenge test ensuring retention even in humid process conditions. These filters will sterilize air or gas streams in:

- Single-use bioreactors
- Filtration and bag assemblies
- Carboys
- Filling vessels
- Transfer vessels
- Fill-finish applications
- Holding/storage bags
- Mixer bags
- Bottle assemblies

Filter Formats

Gamma stable Opticap® XL capsule filters are available in two sizes:

- Opticap® XL 50 Capsules, EFA: 19.6 cm²
- Opticap® XL 300 Capsules, EFA: 480 cm²

Features

The polyethersulfone membrane based gas filter:

- Can be sterilized by gamma irradiation or autoclaving
- Can be integrity tested using HydroCorr™ Water Flow Integrity Test
- Is validated for liquid bacterial retention
- Maintains superior flow rate over extended process duration

Superior Performance, higher process confidence

Gamma stable Opticap® XL Capsule Filters with Millipore Express® SPG Hydrophobic Membrane are sterilizing-grade devices. The hydrophobic membrane inherently resists wetting—providing superior flow rate over extended operating durations. This makes it ideal for single-use applications such as high density cell culture processes using small scale bioreactors, bag assemblies, carboys, and mixers (Figure 1).

Gamma stable Opticap® XL Capsule Filters with Millipore Express® SPG membrane exhibits superior flow rate (up to 3x) when compared to PVDF membrane based vent filters under typical operating conditions. See Figure 3.



Figure 1. Opticap® XL 50 Capsules with Millipore Express® SPG Hydrophobic Membrane installed on a Mobius® CellReady 3L Bioreactor.




Figure 3. Typical Air Flow Rates. Gamma Stable Opticap® XL Capsule Filters with Millipore Express® SPG Hydrophobic membrane offer superior performance.

These filters are ideal for use in single-use redundant filtration assemblies and standard assemblies with filters, where integrity testing and air blow down of the liquid filter can be performed without breaching sterility (Figure 2). In addition, these devices provide sterilization flexibility and can be sterilized using either gamma irradiation or autoclaving. These devices are qualified for gamma sterilization (up to 45 kGy) and/or 3 autoclave cycles. See the Specification table on the last page for more details.





Figure 2. Single-use redundant filtration assembly with gamma stable Opticap® XL 50 Capsule Filter. Opticap® XL 50 Capsule Filter can be used as shown (in green box above) while integrity testing and drying the assemblies to maintain a sterile barrier.



Regulatory Compliance

Gamma Stable Opticap® XL Capsule Filters with Millipore Express® SPG membrane are designed, developed and manufactured in accordance with a Quality Management System approved by an accredited registering body to ISO® 9001 and ISO® 14001 Quality Systems Standard. Every Gamma Stable Opticap® XL Capsule Filter with Millipore Express® SPG Hydrophobic Membrane is integrity tested during manufacturing to guarantee sterility assurance and is supported with a comprehensive Validation Guide for compliance with regulatory requirements.

A Certificate of Quality, available for each filter lot, certifies that Gamma Stable Opticap® XL Capsule Filters with Millipore Express® SPG Hydrophobic Membrane meet quality assurance lot release criteria. The bubble point of the membrane is correlated to bacterial retention using the ASTM® F838-05 challenge methodology. Each capsule filter is labeled with its catalogue number, lot number, serial number, and indicator for direction of flow.

Opticap® XL capsule filters with Millipore Express® membranes are part of the Mobius® integrated disposable bioprocess solution. No matter what your application step or scale, Mobius® solutions can help you achieve greater process efficiency and productivity with the right combination of single-use

products, application solutions, and expert validation support. From disposable process containers to capsule filters and connectors, to validated, gamma compatible turnkey assemblies, Mobius® solutions provide faster turnaround time and reliable performance right out of the box.

Integritest® Manifold: for Parallel Integrity Testing

As a complement to the Opticap® XL 50 capsule filters, we also offer the Integritest® Manifold that can be used with up to 10 filters to perform fast, reliable integrity testing in parallel (Figure 4). This manifold can be used with Opticap® XL 50 and Aervent® 50 capsule filters to improve your efficiency when performing integrity tests such as the HydroCorr™ Water Flow Integrity Test (resistance to water intrusion) and/or bubble point using an automated integrity tester.

The Integritest® Manifold is designed to offer a high degree of integrity test robustness while improving your testing efficiency. This manifold provides flexibility while maintaining ease-of-use with the following features:

- Easy connectivity (TC compatibility and Luer adaptors)
- Hydrocorr™ test friendly; easy to drain, fast setup
- Easy air blow down of filters
- Leak free; easy to clean and sterilize
- Stainless steel welded construction; compatible with IPA
- Vertical design with minimal footprint and no moving parts

Figure 4.
Integritest® Manifold.



Specifications

	Gamma Stable Opticap® XL 50 Capsule Filters with Millipore Express® SPG Hydrophobic membrane	Gamma Stable Opticap® XL 300 Capsule Filters with Millipore Express® SPG Hydrophobic membrane
Filtration Area, m ²	0.00196 (19.6 cm ²)	0.048 (480 cm ²)
Materials of Construction		
Filter Structural Components	Gamma Stable Polypropylene Hydrophobic Polyethersulfone Gamma Stable Polybutylene terephthalate based blend	
Inlet/Outlet	Stepped Hosebarb Inlet and Outlet for silicone flexible tubing with 6.5 to 9.5 mm internal diameter; female Luer slip interior	HH: Stepped Hosebarb (9/16") FF: Sanitary Flange (3/4")
Maximum Pressure, bar (psid)	Forward: 4.2 bar (60 psid) @ 25 °C; Reverse: 2.1 bar (30 psid) @ 25 °C, intermittent	
Hydrocorr™ Test*	≤0.03 mL/min @ 38 psi	≤0.16 mL/min @ 38 psi
Bacterial Retention	Quantitative retention of 10 ⁷ CFU/cm ² <i>Brevundimonas diminuta</i> (ATCC® 19146) per ASTM® F838-05 (2005) methodology	
Bacterial Endotoxin	Aqueous extraction contains ≤0.25 EU/mL as determined by the Limulus Amebocyte Lysate (LAL) Test.	
Gamma/Autoclave Stability	Gamma compatible to 45 kGy; OR up to 3 autoclave cycles of 60 min @ 126 °C; not in-line steam sterilizable	Gamma compatible to 45 kGy; and up to 3 autoclave cycles of 60 min @ 123 °C; not in-line steam sterilizable
Component Material Toxicity	Component materials meet the criteria of the USP <88> Reactivity Test for Class VI Plastics.	
USP Toxicity	Non-toxic per MEM elution ISO® 10993-5	
Good Manufacturing Practices	These products are manufactured in a facility, which adheres to FDA Good Manufacturing Practices.	
Non-Fiber Releasing	Component materials meet the "non-fiber releasing" criteria as defined in 21 CFR 210.3 (b) (6).	

*If you would like to integrity test these filters using alcohol, please contact your EMD Millipore representative for more information.

Ordering Information

Description	Qty/Pk	Catalogue No.
Gamma stable Opticap® XL 50 Capsule Filters with Millipore Express® SPG Hydrophobic membrane. Inlet and outlet: stepped hose barb with female Luer slip interior	100	KEGBG050HH00
	10	KEGBG050HH10
Gamma stable Opticap® XL 300 Capsule Filters with Millipore Express® SPG 0.2µm membrane, 3pk w/ 9/16" hose barb	3	KEGBG003HH3
Gamma stable Opticap® XL 300 Capsule Filters with Millipore Express® SPG 0.2µm membrane, 3pk w/ 3/4" sanitary flange	3	KEGBG003FF3
Integritest® Manifold: for Parallel Integrity Testing	1	P93795



www.emdmillipore.com/offices

To Place an Order or Receive Technical Assistance

In the U.S. and Canada,
call toll-free 1-800-645-5476

For other countries across Europe and the world,
please visit www.emdmillipore.com/offices

For Technical Service, please visit
www.emdmillipore.com/techservice

EMD Millipore, the M mark, Opticap, Millipore Express, Mobius, Integritest, and Aervent are registered trademarks of Merck KGaA, Darmstadt, Germany. HydroCorr is a trademark of Merck KGaA, Darmstadt, Germany. ISO is a registered trademark of the International Organization for Standardization. ATCC is a registered trademark of the American Type Culture Collection. ASTM is a registered trademark of the American Society for Testing and Materials. Lit. No. DS3214EN00 Rev. 5.0 PS-14-09701 07/2015
© 2015 EMD Millipore Corporation, Billerica, MA, USA All rights reserved.