



SIGMA QUALITY CONTROL TEST PROCEDURE  
**Enzymatic Assay of CHOLESTEROL ESTERASE**  
**(EC 3.1.1.13)**

**PRINCIPLE:**

Abbreviation:

4-AAP = 4-Aminoantipyrine

**CONDITIONS:** T = 37°C, pH = 7.0, A<sub>500nm</sub>, Light path = 1 cm

**METHOD:** Continuous Spectrophotometric Rate Determination

**REAGENTS:**

- A. 400 mM Potassium Phosphate Buffer, pH 7.0 at 37°C  
(Prepare 100 ml in deionized water using Potassium Phosphate, Monobasic, Anhydrous, Sigma Prod. No. P-5379. Adjust to pH 7.0 at 37°C with 1 M KOH.)
- B. 0.9% (w/v) Sodium Chloride Solution  
(Prepare 25 ml in deionized water using Sodium Chloride, Sigma Prod. No. S-9625.)
- C. 8.6 mM Cholesteryl Oleate Solution (Chol-Oleate)  
(Prepare 10 ml by first dissolving the Cholesteryl Oleate, Sigma Prod. No. C-9253, in 1 ml of Polyoxyethylene, 9 Lauryl Ether, Sigma Prod. No. P-9641. Stir gently, with heat, until the solution is clear and colorless. Then add 9 ml of hot Reagent B and continue for 5 minutes. Allow the solution to return to ambient temperature prior to use. The solution clears upon cooling.)

**Enzymatic Assay of CHOLESTEROL ESTERASE  
(EC 3.1.1.13)**

**REAGENTS:** (continued)

- D. 15% (w/v) Taurocholic Acid Solution (Tauro)  
(Prepare 10 ml in deionized water using Taurocholic Acid, Sodium Salt, Sigma Prod. No. T-4009.)
- E. 15% (w/v) Cholic Acid Solution (Chol)  
(Prepare 10 ml in deionized water using Cholic Acid, Sodium Salt, Sigma Prod. No. C-1254.)
- F. 1.76% (w/v) 4-Aminoantipyrine Solution (4-AAP)  
(Prepare 1 ml in deionized water using 4-Aminoantipyrine, Free Base, Sigma Prod. No. A-4382.)
- G. 5% (w/v) Phenol Solution (Phenol)  
(Prepare 10 ml in deionized water using Phenol, Sigma Prod. No. P-4161.)
- H. Cholesterol Oxidase Enzyme Solution (Chol Oxid)  
(Immediately before use, prepare a solution containing 20 - 30 units/ml of Cholesterol Oxidase, Sigma Prod. No. C-1512, in cold Reagent A.)
- I. Peroxidase Enzyme Solution (POD)  
(Immediately before use, prepare a solution containing 40 - 60 units/ml of Peroxidase Type II from Horseradish, Sigma Prod. No. P-8250 in cold deionized water.)
- J. Cholesterol Esterase Enzyme Solution  
(Immediately before use, prepare a solution containing 0.25 - 0.85 unit/ml of Cholesterol Esterase in cold Reagent A.)

**PROCEDURE:**

Pipette (in milliliters) the following reagents into suitable cuvettes<sup>1</sup>:

	<u>Test</u>	<u>Blank</u>
Reagent A (Buffer)	2.10	2.10
Reagent D (Tauro)	0.05	0.05
Reagent E (Chol)	0.05	0.05
Reagent I (POD)	0.10	0.10
Reagent C (Chol-Oleate)	0.50	0.50
Reagent G (Phenol)	0.05	0.05

**Enzymatic Assay of CHOLESTEROL ESTERASE  
(EC 3.1.1.13)**

**PROCEDURES:** (continued)

Mix by inversion and equilibrate to 37°C. Then add:

Reagent F (4-AAP)	0.05	0.05
Reagent H (Chol Oxid)	0.05	0.05

Mix by inversion and obtain the baseline at 500 nm. After approximately 5 minutes add:

	<u>Test</u>	<u>Blank</u>
Reagent J (Cholesterol Esterase)	0.05	-----
Reagent A (Buffer)	-----	0.05

Immediately mix by inversion and record the increase in  $A_{500nm}$  for approximately 10 minutes.<sup>2</sup>  
Obtain the  $\Delta A_{500nm}/\text{minute}$  using the maximum linear rate for both the Test and Blank.

**CALCULATIONS:**

$$\text{Units/ml enzyme} = \frac{(\Delta A_{500nm}/\text{min Test} - \Delta A_{500nm}/\text{min Blank})(3)(df)}{(0.5)(13.78)(0.05)} >$$

3 = Total volume (in milliliters) of assay

df = Dilution factor

0.5 = Conversion factor based on one mole of  $H_2O_2$  produces half a mole of Quinoneimine Dye

0.05 = Volume (in milliliters) of enzyme used

13.78 = Millimolar extinction coefficient of Quinoneimine Dye at 500 nm under the assay conditions

$$\text{Units/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

$$\text{Units/mg protein} = \frac{\text{units/ml enzyme}}{\text{mg protein/ml enzyme}}$$

**UNIT DEFINITION:**

One unit will hydrolyze 1.0  $\mu\text{mole}$  of cholesteryl oleate to cholesterol and oleic acid per minute at pH 7.0 at 37°C in the presence of taurocholate.

**Enzymatic Assay of CHOLESTEROL ESTERASE  
(EC 3.1.1.13)**

**FINAL ASSAY CONCENTRATION:**

In a 3.00 ml reaction mix, the final concentrations are 287 mM potassium phosphate, 0.25% (w/v) taurocholic acid, 0.25% (w/v) cholic acid, 4 - 6 units peroxidase, 1.4 mM cholesteryl oleate, 1.7% (v/v) polyoxyethylene 9 lauryl ether, 0.14% (w/v) sodium chloride, 0.083% (w/v) phenol, 0.03% (w/v) 4-aminoantipyrine, 1 - 1.5 units cholesterol oxidase and 0.013 - 0.043 unit cholesterol esterase.

**REFERENCE:**

Allain, C.C. et al., (1974) *Clinical Chemistry*, **20**, 470-475

**NOTES:**

1. Add the reagents in the order written.
2. The fastest rate is usually between 4-8 minutes after addition of the Cholesterol Esterase.
3. Cholesterol Oxidase Unit Definition: One unit will convert 1.0  $\mu$ mole of cholesterol to 4-cholesten-3-one per minute at pH 7.5 at 25°C.
4. Peroxidase Unit Definition: One unit will form 1.0 mg purpurogallin from pyrogallol in 20 sec at pH 6.0 at 20°C.
5. This assay is based on the cited reference.
6. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

**Sigma warrants that the above procedure information is currently utilized at Sigma and that all Sigma-Aldrich, Inc. products conform to the information in this and other Sigma-Aldrich, Inc. publications. Purchaser must determine the suitability of the information and product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packaging slip.**