

Product Information

Fast β -Glucuronidase, Recombinant

From limpets (*Patella vulgata*), expressed in proprietary host, aqueous solution**SRE0095**



Product Description

CAS Registry Number: 9001-45-0

Enzyme Commission (EC) Number: 3.2.1.31

Synonyms: β -D-Glucuronide glucuronosohydrolase

Glucuronidation (conjugation with glucuronic acid) by the human UDP-glucuronosyltransferase (UGT) family of enzymes plays an important role in the metabolic fate of many drugs and other xenobiotics. This biosynthetic reaction also has a role in the conjugation and excretion of endogenous substrates, such as steroids, bilirubin, and bile acids.¹ UGT activity results in the conjugation of glucuronic acid to substrates containing sulfhydryl, hydroxyl, aromatic amino, or carboxylic acid moieties. The glucuronides formed are more polar (water-soluble) than the parent organic substrates, and are generally excreted through the kidney.

 β -glucuronidase catalyzes the general reaction:General β -Glucuronidase Hydrolysis Reaction β -Glucuronidase Hydrolysis Reaction for Morphine 3- β -glucuronide

β -Glucuronidases (GUS) are routinely used for the enzymatic hydrolysis of glucuronides from urine,^{2,3} plasma,^{4,5} and other fluids⁶ prior to analysis by enzyme immunoassay, mass spectrometry, HPLC, gas chromatography, or other methods. Typically, 1-20 units of glucuronidase are used per μL of plasma, urine, or bile for the enzymatic hydrolysis of glucuronides present in these samples.²⁻⁶ β -Glucuronidase from limpets has been shown to be a superior enzyme for the hydrolysis of the widest range of drug-glucuronides from urine.^{7,8}

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Reagent

This recombinant β -glucuronidase, expressed in a proprietary host, is supplied as a solution in 100 mM potassium acetate, pH 5.2. It is free of detergents, carbohydrates, or any other components that may interfere with sample preparation and analysis. It is highly purified to maximize specific activity, and also to eliminate monoacetylmorphine (MAM) esterase activity (6-monoacetylmorphine \rightarrow morphine). Typical analysis at high enzyme concentration (50 units/ μL) exhibited <1% MAM esterase conversion after 4 hours at 60 °C.

This product is formulated with a higher enzyme concentration for applications that require ultra-fast hydrolysis rates.

Glucuronidase Activity: 300,000-400,000 units/mL

Unit Definition: One Sigma or modified Fishman unit will liberate 1.0 mg of phenolphthalein from phenolphthalein glucuronide per hour at 37 °C at pH 3.8 (30-minute assay).

Storage/Stability


This product should be stored at 2-8 °C.

Product Profile

Hydrolysis of codeine 6-β-glucuronide

Codeine 6-β-glucuronide is known as one of the most recalcitrant substrates in many drug analysis panels. The catalytic efficiency of SRE0095 was determined by monitoring the hydrolysis of codeine 6-β-glucuronide to codeine (Figure 1).

Figure 1. LC-MS Analysis of Enzymatic Hydrolysis of Codeine 6-β-Glucuronide with Cat. No. SRE0095.



Complete hydrolysis by Cat. No. SRE0095 was observed in less than 30 minutes with 50 units/μL at 70 °C. Codeine-glucuronide spike level was at 1000 ng/mL in synthetic urine.


The suitability of SRE0095 was also demonstrated for hydrolysis of several other drug-glucuronide substrates:

- Morphine 6-β-Glucuronide
- Hydromorphone Glucuronide
- Oxycodone Glucuronide
- Lorazepam Glucuronide
- Oxazepam Glucuronide
- Testosterone Glucuronide

Optimal temperature

Range: 60-70 °C (see Figure 2)

Figure 2. Temperature Profile of Cat. No. SRE0095 at pH 3.8.




Activity was measured as μg of phenolphthalein liberated per 10 minutes.

Optimal pH

Range: 3.8-5.0 (see Figure 3)

Figure 3. pH Profile of Cat. No. SRE0095 at 37 °C.



Activity was measured as μg of phenolphthalein liberated per 10 minutes.

References

1. Tephly, T.R. *et al.*, *Adv. Pharmacol.*, **42**, 343-346 (1998).
2. Xu, X. *et al.*, *J. Chromatogr. B. Analyt. Technol. Biomed. Life Sci.*, **780(2)**, 315-330 (2002).
3. Staimer, N. *et al.*, *Anal. Chim. Acta*, **444(1)**, 27-36 (2001).
4. Zhai, P. *et al.*, *Am. J. Physiol. Heart Circ. Physiol.*, **281(3)**, H1223-H1232 (2001).
5. Anderson, A. *et al.*, *Cancer Chemother. Pharmacol.*, **44(5)**, 422-426 (1999).
6. Nobilis, M. *et al.*, *J. Chromatogr. A*, **1031(1-2)**, 229-236 (2004).
7. Combie, J. *et al.*, *Clin. Chem.*, **28(1)**, 83-86 (1982).
8. Combie, J. *et al.*, *Res. Commun. Chem. Pathol. Pharmacol.*, **35(1)**, 27-41 (1982).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://www.sigmaaldrich.com/techservice).

Standard Warranty

The applicable warranty for the products listed in this publication may be found at [SigmaAldrich.com/terms](https://www.sigmaaldrich.com/terms).

Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://www.sigmaaldrich.com/offices).

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

SRE0095dat Rev 05/22 DT, RBG, JMB, GCY, MAM

**Millipore
SIGMA**