Instructions for Use

Leukocyte α -Naphthyl Acetate (Non-Specific Esterase) and Leukocyte Naphthol AS-D Chloroacetate (Specific Esterase) Kit

Procedure No. 90

Intended Use

Sigma-Aldrich Naphthol AS-D Chloroacetate Esterase and α -Naphthyl Acetate Esterase kits are intended as a general laboratory use histological stain. Naphthol AS-D chloroacetate esterase and α -naphthyl acetate esterase reagents are for professional "In Vitro Diagnostic Use" only. This manual, qualitative, histochemical procedure demonstrates naphthol AS-D chloroacetate esterase and α -naphthyl acetate esterase activity in leukocytes of human blood, bone marrow films, or tissue touch preparations. Histological visualization of leukocyte cellular esterases is a unique technique currently commonly used in medicine.

Cellular esterases are ubiquitous and appear to represent a series of different enzymes acting upon select substrates. Under defined reaction conditions, it may be possible to determine hemopoietic cell types, using specific esterase substrates. The described methods provide means to distinguish granulocytes from monocytes.¹⁻⁸

To perform the test, blood, bone marrow films or tissue touch preparations are incubated with either naphthol AS-D chloroacetate or α -naphthyl acetate in the presence of a stable diazonium salt. Enzymatic hydrolysis of ester linkages liberates free naphthol compounds. These couple with the diazonium salt, forming highly colored deposits at the sites of enzyme activity.

Reagents

Dimethyl Formamide (Cat. No. 9010-25ML)

Dimethyl formamide, 100% Danger. Flammable liquid and vapour. May be harmful if swallowed. Harmful in contact with skin. Causes mild skin irritation. Causes serious eye irritation. Toxic if inhaled. May damage the unborn child. Obtain special instructions before use.

Ethylene Glycol Monomethyl Ether (Cat. No. 9011-25ML)

2-MethoxyEthanol, 100% Danger. Flammable liquid and vapour. May be harmful if swallowed. Harmful in contact with skin. Causes mild skin irritation. Causes serious eye irritation. Toxic if inhaled. May damage fertility or the unborn child. Obtain special instructions before use.

Naphthol AS-D Chloroacetate (Cat. No. 905-10CAP)

Naphthol AS-D Chloroacetate, 20 mg/cap.

α -Naphthyl Acetate (Cat. No. 906-10CAP)

Naphthol AS-D Chloroacetate, 20 mg/cap Danger. Causes serious eye damage. Wear protective gloves/eye protection/face protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

TRIZMAL™ 6.3 Buffer Concentrate (Cat. No. 903C-50ML)

TRIZMA® maleate, 200 mmol/L, Chloroform added as preservative

TRIZMAL™ 7.6 Buffer Concentrate (Cat. No. 902C-50ML)

TRIZMA® maleate, 200 mmol/L, Chloroform added as preservative

Mayer's Hematoxylin Solution (Cat. No. MHS1-100ML)

Hematoxylin, certified, CI 75290, 0.1% (w/v), and stabilizers

Acid Hematoxylin Solution (Cat. No. 2852-100ML)

Hematoxylin certified, CI 75290, 1 g/L, and stabilizers, pH 3.3 at 25°C

Fast Blue RR Salt (Cat. No. FBS25-10CAP)

Fast Blue, C.I. 37155. Preweighed capsules. Actual weight per capsule will vary with dye lot purity and has been optimized by assay.

Fast Corinth V Salt (Cat. No. 9015-10CAP)

Fast Corinth V Salt, 18-22 mg/cap. Danger. Harmful if swallowed. Harmful in contact with skin. Harmful if inhaled. May cause cancer. Obtain special instructions before use. Wear protective gloves/protective clothing.

Citrate Concentrate (Cat. No. 3861-20ML)

Citrate buffer, 0.38 mol/L, pH 5.4 when diluted according to procedure

Special Materials Required but Not Provided

- Methanol, Absolute
- Acetone, ACS Reagent
- Sodium Fluoride Solution (Cat. No. 919-25ML) Sodium fluoride, 2 g/dL

Storage and Stability

Dimethyl Formamide, Ethylene Glycol Monomethyl Ether, TRIZMAL™ 6.3 Buffer Concentrate, TRIZMAL™ 7.6 Buffer Concentrate, Mayer's Hematoxylin Solution and Acid Hematoxylin Solution, are stored at room temperature (18–26°C).

Naphthol AS-D Chloroacetate, α -Naphthyl Acetate and Fast Blue RR Salt are stored below 0°C.

Fast Corinth V Salt and Citrate Concentrate are stored refrigerated (2–8°C).

Naphthol AS-D Chloroacetate, α -Naphthyl Acetate, Fast Blue RR Salt and Fast Corinth V Salt are stable until the expiration date shown on the labels.

Citrate Dilute Solution is stable for 1 week if stored tightly capped at room temperature (18–26°C).

TRIZMAL™ 6.3 Buffer Concentrate, TRIZMAL™ 7.6 Buffer Concentrate and Citrate Concentrate are suitable for use in the absence of microbial growth.

Sodium fluoride, 2 g/dL. Store at room temperature (18–26°C). Used if " α -Naphthyl Acetate Esterase with Fluoride Inhibition Procedure" is performed.

Deterioration

Discard Dimethyl Formamide and Ethylene Glycol Monomethyl Ether if colored or turbid.

TRIZMAL™ 6.3 Dilute Buffer Solution and TRIZMAL™ 7.6 Dilute Buffer Solution should be used once then discarded.

Discard Mayer's Hematoxylin and Acid Hematoxylin Solution when the time required for suitable staining exceeds the time recommended in the procedure by more than 5 minutes.

Preparation

Naphthol AS-D Chloroacetate Solution is prepared by dissolving contents of 1 capsule Naphthol AS-D Chloroacetate in 2 mL Dimethyl Formamide. Remove 1 capsule from freezer as needed. Prepare immediately prior to use.

α -Naphthyl Acetate Solution is prepared by dissolving contents of 1 capsule α -Naphthyl Acetate in 2 mL Ethylene Glycol Monomethyl Ether. Remove 1 capsule from freezer as needed. Prepare immediately prior to use.

TRIZMAL™ 6.3 Dilute Buffer Solution is prepared by diluting 1 part TRIZMAL™ 6.3 Buffer Concentrate with 9 parts deionized water. The pH should be 6.3 at 25°C.

TRIZMAL™ 7.6 Dilute Buffer Solution is prepared by diluting 1 part TRIZMAL™ 7.6 Buffer Concentrate with 9 parts deionized water. The pH should be 7.6 at 25°C.

Mayer's Hematoxylin and Acid Hematoxylin Solution should be filtered before use.

Citrate Dilute Solution is prepared by diluting 1 part Citrate Concentrate with 9 parts deionized water. pH 5.4 when diluted.

Citrate-Acetone-Methanol Fixative: To 18 mL Citrate Dilute Solution, add 27 mL ACS grade Acetone and 5 mL Methanol. Store tightly capped at room temperature (18–26°C). Discard after 8 hours.

Precautions

The IVDs included in these kits are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Blood, bone marrow films, tissue-touch preparations, and cytocentrifuge preparations may be used with both α -naphthyl acetate esterase and naphthol AS-D chloroacetate esterase. Either EDTA or heparin will serve as an anticoagulant.⁹ Frozen and paraffin embedded tissues may be used with naphthol AS-D chloroacetate esterase. α -Naphthyl acetate esterase may be used successfully on frozen tissue sections.¹⁰ Blood or bone marrow films may be stored fixed at room temperature (18–26°C) for several weeks or unfixed for several days without appreciable change in activity.^{5,9} Do not ship whole blood for assay at other laboratories. Send fixed or unfixed slides. Slides should be kept cool during transit. Allow films to dry at least 1 hour prior to fixation.

Notes

The described procedures are performed at 37°C. If reagents are not at this temperature, weak or negative reactions may be obtained. It is recommended that temperatures be checked with an accurate thermometer. Controlled temperature water baths are more efficient than warm air incubators and should be used for enzyme cytochemical methods. Heat transfer through glass is more rapid than through plastic, thus, glass Coplin jars should be employed.

Many enzyme systems are sensitive to minute traces of detergent. Washing glassware with dilute bleach followed by rinsing in copious quantities of deionized water will prevent detergent effect upon cellular enzymes.

Results are based on a certain degree of subjective interpretation. Individual laboratories should establish their own normal ranges.

Procedures

Naphthol AS-D Chloroacetate Esterase Procedure

1. Fix slides for 1 minute in Citrate-Acetone-Methanol Fixative at room temperature (18–26°C).
2. Wash thoroughly in deionized water and air dry at least 20 minutes.
3. To 50 mL TRIZMAL™ 6.3 Dilute Buffer Solution, PREWARMED TO 37°C, add with constant stirring, contents of 1 capsule Fast Corinth V Salt.
4. When salt is completely dissolved in buffer, add 2 mL Naphthol AS-D Chloroacetate Solution. The solution will appear quite turbid.
5. Continue mixing for 15–30 seconds, then add to Coplin jar. DO NOT FILTER.
6. Place specimens in staining solution (from Step 5) and incubate at 37°C for 5 minutes. NOTE: PROTECT FROM LIGHT.
7. Remove slides from stain and wash in deionized water for 3 minutes. Discard staining solution.
8. If desired, counterstain in Acid Hematoxylin Solution for 5–10 minutes and wash in tap water.
9. Air dry slides and examine microscopically. If coverslipping is required use only an aqueous mounting media.

α -Naphthyl Acetate Esterase Procedure

1. Fix slides in Citrate-Acetone-Methanol Fixative for 1 minute at room temperature (18–26°C).
2. Wash thoroughly in deionized water and air dry at least 20 minutes.
3. To 50 mL TRIZMAL™ 7.6 Dilute Buffer Solution, PREWARMED TO 37°C, add with constant stirring, contents of 1 capsule Fast Blue RR Salt.
4. When salt is completely dissolved in buffer, add 2 mL α -Naphthyl Acetate Solution. The solution will be yellow and slightly turbid.
5. Continue stirring for 15–20 seconds, then add to Coplin jar. DO NOT FILTER.
6. Place specimens in staining solution (from Step 5) and incubate at 37°C for 30 minutes. NOTE: PROTECT FROM LIGHT.
7. Remove slides from stain and wash for 3 minutes in deionized water. Discard staining solution.
8. If desired, counterstain for 5–10 minutes in Mayer's Hematoxylin Solution, and wash in tap water.
9. Air dry slides and examine microscopically. If coverslipping is required use only an aqueous mounting media.

Double Staining Esterase Procedure

1. Perform α -Naphthyl Acetate Esterase test as described in Procedure. Do not counterstain.
2. Rinse slides 5 minutes in deionized water.
3. Perform Naphthol AS-D Chloroacetate Esterase test as described in procedure Steps 3-9.

α -Naphthyl Acetate Esterase with Fluoride Inhibition Procedure

Although α -naphthyl acetate esterase is found primarily in cells of monocytic lineage when performed as described, it should be recognized that megakaryocytes and erythroid precursors are positive for this enzyme.¹¹ Lymphocytes and some mature granulocytes also show occasional positivity.⁵ To differentiate these cells conclusively from monocytes, sodium fluoride is incorporated with the incubation system. The monocyte enzyme is inactivated in the presence of this compound.¹² The following procedure may be used to perform the fluoride inhibition test.

1. Fix slides in Citrate-Acetone-Methanol Fixative for 1 minute at room temperature (18–26°C).
2. Wash thoroughly in deionized water and air dry at least 20 minutes.
3. Label 2 beakers A and B, and add the following:

	Beaker A	Beaker B
Prewarmed 37°C TRIZMAL™ 7.6 Dilute Buffer	50 mL	50 mL
Add with constant stirring, Fast Blue RR	1 capsule*	1 capsule*
α -Naphthyl Acetate Solution	2 mL	2 mL
Sodium Fluoride Solution	-	2 mL

*Contents of 1 capsule

4. Mix well and pour into Coplin jars labeled A and B.
5. Proceed as described in Steps 6–9 of α -Naphthyl Acetate Esterase Procedure.

Performance Characteristics

Method of Scoring

Scan the film and select a thin area with few erythrocytes. Sites of Naphthol AS-D Chloroacetate Esterase activity will appear as bright red granulation, α -Naphthyl Acetate Esterase as black granulation. Rate from 0 to 4+ on the basis of quantity and intensity of individual dyes within the cytoplasm of the respective cell types. Characteristics of scoring are based somewhat on subjective interpretation. A suggested scoring format is presented in Table 1. Conclusions center on relative presence or absence of staining.

Table 1. Characteristics of Scoring

Cell Rating	Staining Intensity	Interpretation
0+	None	—
1+	Faint to Moderate	±
2+	Moderate to Strong	+
3+	Strong	+
4+	Brilliant	+

Results

Naphthol AS-D Chloroacetate Esterase

This enzyme is usually considered specific for cells of granulocytic lineage. The cells should show red granulation. Activity is weak or absent in monocytes and lymphocytes.

α -Naphthyl Acetate Esterase

Under the assay conditions (pH 7.6), this enzyme is detected primarily in monocytes, macrophages and histiocytes, and is virtually absent in granulocytes. Monocytes should show black granulation. Lymphocytes may occasionally exhibit activity.

α -Naphthyl Acetate Esterase with Fluoride Inhibition

All cells of monocytic lineage will be negative for enzyme activity, with the exception of differentiated histiocytes or specialized macrophages in tissue which may also be resistant to sodium fluoride.¹⁰

Double Staining Esterase

Specimens taken through the double staining procedure will demonstrate the granulocytes with red granulation and monocytes with black granulation.

The expected cellular reactivity of tests for esterase activity is summarized in Table II.

Table II. Cytochemical Reactions of Blood Cells

Cell Type	Naphthol AS-D Chloroacetate Esterase	α -Naphthyl Acetate Esterase
Myeloblasts	±	±
Promyelocytes	+	±
Neutrophils	+	—
Eosinophils	—	—
Basophils	±	—
Monocytes	—	+
Lymphocytes	—	±
Lymphoblasts	—	±
Megakaryocytes	—	+
Erythroblasts	—	±
Plasma Cells	—	±
Mast Cells	+	—
Hairy Cells	—	±
Histiocytes	±	+

The reagent system should be monitored by the use of positive and negative control slides. Positive control slides may be prepared from specific cell lines known to be positive.

Alternately, anti-coagulated blood from normal specimens (preferably with increased monocyte count if using α -naphthyl acetate esterase procedure) may also be used; however, they will provide less intense staining and will have fewer positive cells.

Known negative patient slides may be used as a negative control. If unavailable, staining a specimen in an incubation mixture with the substrate omitted will give the desired results. However, use of the former is highly recommended.

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
9010	Dimethyl Formamide	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
9011	Ethylene Glycol Monomethyl Ether	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
905	Naphthol AS-D Chloroacetate	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
906	α -Naphthyl Acetate	Identity	3 of 3	3 of 3	3 of 3	3 of 3
903C	TRIZMAL™ 6.3 Buffer Concentrate	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
902C	TRIZMAL™ 7.6 Buffer Concentrate	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
MHS1	Mayers Hematoxylin Solution	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
2852	Acid Hematoxylin Solution	Nuclei	3 of 3	3 of 3	3 of 3	3 of 3
FBS25	Fast Blue RR Salt	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
9015	Fast Corinth V Salt	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3
3861	Citrate Concentrate	Sites of Activity	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

390A:

H301 + H311 + H331 Toxic if swallowed, in contact with skin or if inhaled.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H318 Causes serious eye damage.

H341 Suspected of causing genetic defects.

H350 May cause cancer.

H410 Very toxic to aquatic life with long lasting effects.

P273 Avoid release to the environment.

P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.

P301 + P310 IF SWALLOWED: Immediately call a POISON CENTER/ doctor.

P302 + P352 + P312 IF ON SKIN: Wash with plenty of water. Call a POISON CENTER/ doctor if you feel unwell.

P304 + P340 + P311 IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a POISON CENTER/ doctor.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer
	Indicates the Authorised Representative in Switzerland		

References

1. Beckmann H, Neth R, Gaedicke G, Lanbeck G, Schoch G, Wiegers U, Winkler K: Cytology and cytochemistry of the leukemic cell. *Haematol Bluttransfus* 14:26, 1974.
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. *Blood Cells* 1:101, 1975.
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. *Clinics in Haematol* 1:49, 1972.
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. *Am J Clin Pathol* 55:283, 1971.
5. Yam, LT, Li CY, Wolfe HJ, Moy PW: Histochemical study of acute leukemia. *Arch Pathol* 97:129, 1974.
6. Burstone MS: The cytochemical localization of esterase. *J Natl Cancer Inst* 18:167, 1957.
7. Moloney WC, McPherson K, Sliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. *J Histochem Cytochem* 8:200, 1960.
8. Brown BA: *IN Hematology: Principles and Procedures*. Lea and Febiger, Philadelphia, PA, 1984, pp 127–130.
9. Sun T: *Atlas of Cytochemistry and Immunochemistry of Hematologic Neoplasms*. American Society of Clinical Pathologists Press, Chicago, IL, 1985, pp 24, 38.
10. Hayhoe FGJ, Flemans RJ: *IN Color Atlas of Hematological Cytology*. John Wiley & Sons, New York, NY, 1982, pp 34, 111.
11. Li CY, Lam KW, Lam LT: Esterases in human leukocytes. *J Histochem Cytochem* 21:1, 1973.74.

Contact Information

To place an order, please visit our web site at [SigmaAldrich.com](https://www.SigmaAldrich.com). For Technical Service, please visit the tech service page on our web site at [SigmaAldrich.com/techservice](https://www.SigmaAldrich.com/techservice).

Revision History

Rev. 4.0 2014

Rev. 5.0 2016

Rev. 6.0 2023

Transferred to new template with current branding. Specified for professional use in intended use and precautions. Moved aid to diagnosis statement to intended use. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Added Warnings and Hazards. Added CH-REP information.

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271
Darmstadt,
Germany

MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau
Switzerland

The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2023 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

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Referenzen

1. Beckmann H, Neth R, Gaedicke G, Lanbeck G, Schoch G, Wiegers U, Winkler K: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974.
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975.
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972.
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971.
5. Yam, LT, Li CY, Wolfe HJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97:129, 1974.
6. Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957.
7. Moloney WC, McPherson K, Sliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960.
8. Brown BA: IN Hematology: Principles and Procedures. Lea and Febiger, Philadelphia, PA, 1984, pp 127-130.
9. Sun T: Atlas of Cytochemistry and Immunochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, IL, 1985, pp 24, 38.
10. Hayhoe FGJ, Flemans RJ: IN Color Atlas of Hematological Cytology. John Wiley & Sons, New York, NY, 1982, pp 34, 111.
11. Li CY, Lam KW, Lam LT: Esterases in human leukocytes. J Histochem Cytochem 21:1, 1973.74.

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Références

1. Beckmann H, Neth R, Gaedicke G, Lanbeck G, Schoch G, Wieggers U, Winkler K: Cytology and cytochemistry of the leukemic cell. *Haematol Bluttransfus* 14:26, 1974.
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. *Blood Cells* 1:101, 1975.
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. *Clinics in Haematol* 1:49, 1972.
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. *Am J Clin Pathol* 55:283, 1971.
5. Yam, LT, Li CY, Wolfe HJ, Moy PW: Histochemical study of acute leukemia. *Arch Pathol* 97:129, 1974.
6. Burstone MS: The cytochemical localization of esterase. *J Natl Cancer Inst* 18:167, 1957.
7. Moloney WC, McPherson K, Sliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. *J Histochem Cytochem* 8:200, 1960.
8. Brown BA: *IN Hematology: Principles and Procedures*. Lea and Febiger, Philadelphia, PA, 1984, pp 127-130.
9. Sun T: *Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms*. American Society of Clinical Pathologists Press, Chicago, IL, 1985, pp 24, 38.
10. Hayhoe FGJ, Flemans RJ: *IN Color Atlas of Hematological Cytology*. John Wiley & Sons, New York, NY, 1982, pp 34, 111.
11. Li CY, Lam KW, Lam LT: Esterases in human leukocytes. *J Histochem Cytochem* 21:1, 1973.74.

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Bibliografia

1. Beckmann H, Neth R, Gaedicke G, Lanbeck G, Schoch G, Wiegner U, Winkler K: Cytology and cytochemistry of the leukemic cell. *Haematol Bluttransfus* 14:26, 1974.
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. *Blood Cells* 1:101, 1975.
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. *Clinics in Haematol* 1:49, 1972.
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. *Am J Clin Pathol* 55:283, 1971.
5. Yam, LT, Li CY, Wolfe HJ, Moy PW: Histochemical study of acute leukemia. *Arch Pathol* 97:129, 1974.
6. Burstone MS: The cytochemical localization of esterase. *J Natl Cancer Inst* 18:167, 1957.
7. Moloney WC, McPherson K, Sliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. *J Histochem Cytochem* 8:200, 1960.
8. Brown BA: *IN Hematology: Principles and Procedures*. Lea and Febiger, Philadelphia, PA, 1984, pp 127-130.
9. Sun T: *Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms*. American Society of Clinical Pathologists Press, Chicago, IL, 1985, pp 24, 38.
10. Hayhoe FGJ, Flemans RJ: *IN Color Atlas of Hematological Cytology*. John Wiley & Sons, New York, NY, 1982, pp 34, 111.
11. Li CY, Lam KW, Lam LT: Esterases in human leukocytes. *J Histochem Cytochem* 21:1, 1973.74.

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Referencias

1. Beckmann H, Neth R, Gaedicke G, Lanbeck G, Schoch G, Wieggers U, Winkler K: Cytology and cytochemistry of the leukemic cell. *Haematol Bluttransfus* 14:26, 1974.
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. *Blood Cells* 1:101, 1975.
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. *Clinics in Haematol* 1:49, 1972.
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. *Am J Clin Pathol* 55:283, 1971.
5. Yam, LT, Li CY, Wolfe HJ, Moy PW: Histochemical study of acute leukemia. *Arch Pathol* 97: 129, 1974.
6. Burstone MS: The cytochemical localization of esterase. *J Natl Cancer Inst* 18:167, 1957.
7. Moloney WC, McPherson K, Sliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. *J Histochem Cytochem* 8:200, 1960.
8. Brown BA: *IN Hematology: Principles and Procedures*. Lea and Febiger, Philadelphia, PA, 1984, pp 127-130.
9. Sun T: *Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms*. American Society of Clinical Pathologists Press, Chicago, IL, 1985, pp 24, 38.
10. Hayhoe FGJ, Flemans RJ: *IN Color Atlas of Hematological Cytology*. John Wiley & Sons, New York, NY, 1982, pp 34, 111.
11. Li CY, Lam KW, Lam LT: Esterases in human leukocytes. *J Histochem Cytochem* 21:1, 1973.74.

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Referencer

1. Beckmann H, Neth R, Gaedicke G, Lanbeck G, Schoch G, Wiegers U, Winkler K: Cytology and cytochemistry of the leukemic cell. *Haematol Bluttransfus* 14:26, 1974.
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. *Blood Cells* 1:101, 1975.
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. *Clinics in Haematol* 1:49, 1972.
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. *Am J Clin Pathol* 55:283, 1971.
5. Yam, LT, Li CY, Wolfe HJ, Moy PW: Histochemical study of acute leukemia. *Arch Pathol* 97:129, 1974. Burstone MS: The cytochemical localization of esterase. *J Natl Cancer Inst* 18:167, 1957.
6. Moloney WC, McPherson K, Sliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. *J Histochem Cytochem* 8:200, 1960.
7. Brown BA: *IN Hematology: Principles and Procedures*. Lea and Febiger, Philadelphia, PA, 1984, pp 127-130.
8. Sun T: *Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms*. American Society of Clinical Pathologists Press, Chicago, IL, 1985, pp 24, 38.
9. Hayhoe FGJ, Flemans RJ: *IN Color Atlas of Hematological Cytology*. John Wiley & Sons, New York, NY, 1982, pp 34, 111.
10. Li CY, Lam KW, Lam LT: Esterases in human leukocytes. *J Histochem Cytochem* 21:1, 1973.74.

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Referenser

1. Beckmann H, Neth R, Gaedicke G, Lanbeck G, Schoch G, Wiegers U, Winkler K: Cytology and cytochemistry of the leukemic cell. *Haematol Bluttransfus* 14:26, 1974.
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. *Blood Cells* 1:101, 1975.
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. *Clinics in Haematol* 1:49, 1972.
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. *Am J Clin Pathol* 55:283, 1971.
5. Yam, LT, Li CY, Wolfe HJ, Moy PW: Histochemical study of acute leukemia. *Arch Pathol* 97:129, 1974.
6. Burstone MS: The cytochemical localization of esterase. *J Natl Cancer Inst* 18:167, 1957.
7. Moloney WC, McPherson K, Sliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. *J Histochem Cytochem* 8:200, 1960.
8. Brown BA: *IN Hematology: Principles and Procedures*. Lea and Febiger, Philadelphia, PA, 1984, pp 127-130.
9. Sun T: *Atlas of Cytochemistry and Immunochemistry of Hematologic Neoplasms*. American Society of Clinical Pathologists Press, Chicago, IL, 1985, pp 24, 38.
10. Hayhoe FGJ, Flemans RJ: *IN Color Atlas of Hematological Cytology*. John Wiley & Sons, New York, NY, 1982, pp 34, 111.
11. Li CY, Lam KW, Lam LT: Esterases in human leukocytes. *J Histochem Cytochem* 21:1, 1973.74.

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Referências

1. Beckmann H, Neth R, Gaedicke G, Lanbeck G, Schoch G, Wieggers U, Winkler K: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974.
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975.
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972.
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971.
5. Yam, LT, Li CY, Wolfe HJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97:129, 1974.
6. Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957.
7. Moloney WC, McPherson K, Sliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960.
8. Brown BA: IN Hematology: Principles and Procedures. Lea and Febiger, Philadelphia, PA, 1984, pp 127-130.
9. Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, IL, 1985, pp 24, 38.
10. Hayhoe FGJ, Flemans RJ: IN Color Atlas of Hematological Cytology. John Wiley & Sons, New York, NY, 1982, pp 34, 111.
11. Li CY, Lam KW, Lam LT: Esterases in human leukocytes. J Histochem Cytochem 21:1, 1973.74.

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
- Beckmann H, Neth R, Gaedicke G, Lanbeck G, Schoch G, Wiegiers U, Winkler K: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974.
- Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975.
- Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972.
- Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971.
- Yam, LT, Li CY, Wolfe HJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97:129, 1974.
- Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957.
- Moloney WC, McPherson K, Sliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960.
- Brown BA: IN Hematology: Principles and Procedures. Lea and Febiger, Philadelphia, PA, 1984, pp 127-130.
- Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, IL, 1985, pp 24, 38.
- Hayhoe FGJ, Flemans RJ: IN Color Atlas of Hematological Cytology. John Wiley & Sons, New York, NY, 1982, pp 34, 111.
- Li CY, Lam KW, Lam LT: Esterases in human leukocytes. J Histochem Cytochem 21:1, 1973.74.

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- Beckmann H, Neth R, Gaedicke G, Lanbeck G, Schoch G, Wieggers U, Winkler K: Cytology and cytochemistry of the leukemic cell. *Haematol Bluttransfus* 14:26, 1974.
- Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. *Blood Cells* 1:101, 1975.
- Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. *Clinics in Haematol* 1:49, 1972.
- Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. *Am J Clin Pathol* 55:283, 1971.
- Yam, LT, Li CY, Wolfe HJ, Moy PW: Histochemical study of acute leukemia. *Arch Pathol* 97:129, 1974.
- Burstone MS: The cytochemical localization of esterase. *J Natl Cancer Inst* 18:167, 1957.
- Moloney WC, McPherson K, Sliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. *J Histochem Cytochem* 8:200, 1960.
- Brown BA: *IN Hematology: Principles and Procedures*. Lea and Febiger, Philadelphia, PA, 1984, pp 127–130.
- Sun T: *Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms*. American Society of Clinical Pathologists Press, Chicago, IL, 1985, pp 24, 38.
- Hayhoe FGJ, Flemans RJ: *IN Color Atlas of Hematological Cytology*. John Wiley & Sons, New York, NY, 1982, pp 34, 111.
- Li CY, Lam KW, Lam LT: Esterases in human leukocytes. *J Histochem Cytochem* 21:1, 1973.74.

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Reference

1. Beckmann H, Neth R, Gaedicke G, Lanbeck G, Schoch G, Wiegiers U, Winkler K: Cytology and cytochemistry of the leukemic cell. *Haematol Bluttransfus* 14:26, 1974.
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. *Blood Cells* 1:101, 1975.
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. *Clinics in Haematol* 1:49, 1972.
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. *Am J Clin Pathol* 55:283, 1971.
5. Yam, LT, Li CY, Wolfe HJ, Moy PW: Histochemical study of acute leukemia. *Arch Pathol* 97:129, 1974.
6. Burstone MS: The cytochemical localization of esterase. *J Natl Cancer Inst* 18:167, 1957.
7. Moloney WC, McPherson K, Sliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. *J Histochem Cytochem* 8:200, 1960.
8. Brown BA: *IN Hematology: Principles and Procedures*. Lea and Febiger, Philadelphia, PA, 1984, pp 127–130.
9. Sun T: *Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms*. American Society of Clinical Pathologists Press, Chicago, IL, 1985, pp 24, 38.
10. Hayhoe FGJ, Flemans RJ: *IN Color Atlas of Hematological Cytology*. John Wiley & Sons, New York, NY, 1982, pp 34, 111.
11. Li CY, Lam KW, Lam LT: Esterases in human leukocytes. *J Histochem Cytochem* 21:1, 1973.74.

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Referanser

1. Beckmann H, Neth R, Gaedicke G, Lanbeck G, Schoch G, Wiegers U, Winkler K: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974.
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975.
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972.
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971.
5. Yam, LT, Li CY, Wolfe HJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97:129, 1974.
6. Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957.
7. Moloney WC, McPherson K, Sliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960.
8. Brown BA: IN Hematology: Principles and Procedures. Lea and Febiger, Philadelphia, PA, 1984, pp 127–130.
9. Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, IL, 1985, pp 24, 38.
10. Hayhoe FGJ, Flemans RJ: IN Color Atlas of Hematological Cytology. John Wiley & Sons, New York, NY, 1982, pp 34, 111.
11. Li CY, Lam KW, Lam LT: Esterases in human leukocytes. J Histochem Cytochem 21:1, 1973.74.

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Referanslar

1. Beckmann H, Neth R, Gaedicke G, Lanbeck G, Schoch G, Wiegers U, Winkler K: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974.
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975.
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972.
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971.
5. Yam, LT, Li CY, Wolfe HJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97:129, 1974.
6. Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957.
7. Moloney WC, McPherson K, Sliegerman L: Esterase activity in leukocytes demonstrated by the use of naftol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960.
8. Brown BA: IN Hematology: Principles and Procedures. Lea and Febiger, Philadelphia, PA, 1984, s. 127-130.
9. Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, IL, 1985, s. 24, 38.
10. Hayhoe FGJ, Flemans RJ: IN Color Atlas of Hematological Cytology. John Wiley & Sons, New York, NY, 1982, s. 34, 111.
11. Li CY, Lam KW, Lam LT: Esterases in human leukocytes. J Histochem Cytochem 21:1, 1973.74.

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Lähteet

- Beckmann H, Neth R, Gaedicke G, Lanbeck G, Schoch G, Wiegiers U, Winkler K: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974.
- Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975.
- Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972.
- Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971.
- Yam, LT, Li CY, Wolfe HJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97: 129, 1974.
- Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957.
- Moloney WC, McPherson K, Sliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960.
- Brown BA: IN Hematology: Principles and Procedures. Lea and Febiger, Philadelphia, PA, 1984, pp 127-130.
- Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, IL, 1985, pp 24, 38.
- Hayhoe FGJ, Flemans RJ: IN Color Atlas of Hematological Cytology. John Wiley & Sons, New York, NY, 1982, pp 34, 111.
- Li CY, Lam KW, Lam LT: Esterases in human leukocytes. J Histochem Cytochem 21:1, 1973.74.

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Referenties

1. Beckmann H, Neth R, Gaedicke G, Lanbeck G, Schoch G, Wiegers U, Winkler K: Cytology and cytochemistry of the leukemic cell. Haematol Bluttransfus 14:26, 1974.
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. Blood Cells 1:101, 1975.
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. Clinics in Haematol 1:49, 1972.
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. Am J Clin Pathol 55:283, 1971.
5. Yam, LT, Li CY, Wolfe HJ, Moy PW: Histochemical study of acute leukemia. Arch Pathol 97:129, 1974. Burstone MS: The cytochemical localization of esterase. J Natl Cancer Inst 18:167, 1957.
6. Moloney WC, McPherson K, Slierman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. J Histochem Cytochem 8:200, 1960.
7. Brown BA: IN Hematology: Principles and Procedures. Lea and Febiger, Philadelphia, PA, 1984, pp 127-130.
8. Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, IL, 1985, pp 24, 38.
9. Hayhoe FGJ, Flemans RJ: IN Color Atlas of Hematological Cytology. John Wiley & Sons, New York, NY, 1982, pp 34, 111.
10. Li CY, Lam KW, Lam LT: Esterases in human leukocytes. J Histochem Cytochem 21:1, 1973.74.

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Reference

1. Beckmann H, Neth R, Gaedicke G, Lanbeck G, Schoch G, Wiegers U, Winkler K: Cytology and cytochemistry of the leukemic cell. *Haematol Bluttransfus* 14:26, 1974.
2. Bennet JM, Reed CE: Acute leukemia cytochemical profile: Diagnostic and clinical implications. *Blood Cells* 1:101, 1975.
3. Cawley JC, Hayhoe FGJ: Acute leukemia: Cellular morphology, cytochemistry and fine structure. *Clinics in Haematol* 1:49, 1972.
4. Yam LT, Li CY, Crosby WH: Cytochemical identification of monocytes and granulocytes. *Am J Clin Pathol* 55:283, 1971.
5. Yam, LT, Li CY, Wolfe HJ, Moy PW: Histochemical study of acute leukemia. *Arch Pathol* 97:129, 1974.
6. Burstone MS: The cytochemical localization of esterase. *J Natl Cancer Inst* 18:167, 1957.
7. Moloney WC, McPherson K, Sliegerman L: Esterase activity in leukocytes demonstrated by the use of naphthol AS-D chloroacetate substrate. *J Histochem Cytochem* 8:200, 1960.
8. Brown BA: IN Hematology: Principles and Procedures. Lea and Febiger, Philadelphia, PA, 1984, pp 127-130.
9. Sun T: Atlas of Cytochemistry and Immunocytochemistry of Hematologic Neoplasms. American Society of Clinical Pathologists Press, Chicago, IL, 1985, pp 24, 38.
10. Hayhoe FGJ, Flemans RJ: IN Color Atlas of Hematological Cytology. John Wiley & Sons, New York, NY, 1982, pp 34, 111.
11. Li CY, Lam KW, Lam LT: Esterases in human leukocytes. *J Histochem Cytochem* 21:1, 1973.74.

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MDSS GmbH