Instructions for Use

Harris Hematoxylin Solution

Procedure No. HHS

Intended Use

Harris Hematoxylin Solutions are nuclear stains intended for use in Histology and Cytology. Harris Hematoxylin Solutions are for "In Vitro Diagnostic Use". For professional use only. The data obtained from this manual qualitative procedure is used for the determination of chromatin in human specimens. This data can be used as an aid to diagnosis of certain clinical conditions or pathophysiological states as hematoxylin identifies the Nuclei of cells. It should be reviewed in conjunction with other clinical diagnostic tests or information.

Hematoxylin, a common nuclear stain, is isolated from an extract of logwood.¹ The first successful biologic application of hematoxylin was described by Bohmer¹ in 1865. Since then numerous formulations have appeared. Of these, Harris', Gill's, Mayer's and Weigert's have retained popularity.

Before hematoxylin can be used as a nuclear stain, it must be oxidized to hematein and combined with a metallic ion (mordant). Most successful mordants have been salts of aluminum or iron.

Hematoxylin Solutions are regressive stains for use in routine histology and cytology. The positively charged aluminum-hematein complex combines with negatively charged phosphatase of nuclear DNA forming the blue purple color characteristic of hematoxylin stains.

Harris Hematoxylin Solution may also be used in conjunction with Papanicolaou staining procedures for cytology use. See Sigma-Aldrich Procedure No. HT40.

Reagents

Harris Hematoxylin Solution

(Cat. No. HHS: HHS16-500ML; HHS32-1L; HHS80-2.5L; HHS128-4L)

Certified hematoxylin, 7.0 g/L, C.I. 75290, sodium iodate, aluminum ammonium sulfate • 12 H₂O, preservative and stabilizers

Special Materials Required but Not Provided

- Differentiation Solution, Catalog Nos. A3179-1L or A3429-4L
- Eosin Y Solution Counterstains:
Eosin Y Solution, Alcoholic (Cat. No. HT1101: HT110116-500ML; HT110132-1L; HT110180-2.5L; HT1101128-4L)
OR
Eosin Y Solution, Alcoholic, With Phloxine (Cat. No. HT1103: HT110316-500ML; HT110332-1L; HT110380-2.5L; HT1103128-4L)
- Reagent Alcohol (Cat. No. R8382-1GA) or Ethanol, 100%
- Scott's Tap Water Substitute Concentrate, 10x (Cat. No. S5134-6x100ML)
- Xylene or Xylene substitute
- Hydrochloric Acid, Concentrated
- Microscope, Microscope slides, coverslips, and staining dishes

Storage and Stability

Store reagent at room temperature (18-26°C) protected from light. Reagent is stable until expiration date shown on the label.

Deterioration

Discard if staining time becomes excessive or if solution color changes from plum to blue or brown.

Preparation

Filter Harris Hematoxylin Solution before each use.

Scott's Tap Water Substitute is prepared by mixing 1 part of Scott's Tap Water Substitute Concentrate with 9 parts deionized water.

Precautions

These IVDs are intended for in vitro diagnostic use in a clinical laboratory environment. These IVDs are for professional use by qualified personnel only. Sigma-Aldrich IVDs may be operated by laboratory personnel who are trained to handle human specimens that can be infectious, use microscopes and other laboratory equipment and have color perception and visual acuity to distinguish colors and other objects under a microscope.

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations.

Procedure

Specimen Collection

No known test method can offer complete assurance that blood samples or tissue will not transmit infection. Therefore, all blood derivatives or tissue specimens should be considered potentially infectious.

Standard histology texts provide necessary details for specimen collection and storage.^{2,3}

Notes

- Staining times may be varied for individual color preference.
- Other dilute alkaline solutions may be used in place of Scott's Tap Water Substitute.
- A 0.25% acid alcohol solution may be used in place of Differentiation Solution. Prepare by adding 0.25 mL concentrated Hydrochloric Acid to 100 mL of 70% alcohol.
- The times given in the insert are approximate. Personal preferences will vary and the times can be adjusted to suit personal preferences. Stain solutions which are heavily used will lose their staining powers and the staining times should be lengthened or new solutions should be used.⁴
- Some tap water supplies are acidic and unsuitable for use in the "blueing" portion of this procedure. If tap water is acidic, use a dilute alkaline solution.
- Purple or red-brown nuclei are indicative of inadequate "blueing".
- If eosin staining is excessive, nuclear staining may be masked. Proper eosin staining will demonstrate a 3-tone effect. To increase differentiation of eosin, extend time in alcohols or use a first alcohol with a higher water content. The times in the alcohols may be adjusted to obtain the proper degree of Eosin staining.
- Positive control slides should be included in each run.

Procedure

- Prepare a 95% alcohol solution by adding 5 mL deionized water to 95 mL Reagent Alcohol or Ethanol (100%).
- Deparaffinize to water or fix and dehydrate frozen sections.
- Stain in Harris Hematoxylin Solution for 2.0 to 2.5 minutes.
- Rinse slide in running tap water.
- Differentiation Solution for 1-2 dips.
- Rinse slide in running tap water.
- Blue in Scott's Tap Water Substitute for 5-60 seconds.
- Reagent Alcohol, 95% for 30 seconds.
- Eosin Y Solution Counterstain:
Eosin Y Solution, Alcoholic (Cat. No. HT1101) for 30-60 seconds.
OR
Eosin Y Solution, Alcoholic, With Phloxine (Cat. No. HT1103) for 30-60 seconds.
- Dehydrate, clear and mount.

Performance Characteristics

Nuclear chromatin should be blue. Nucleoli should be conspicuous and crisply outlined. Cytoplasm will display various shades of pink to pink-orange depending upon the counterstain used and RBC's will be red.

If observed results vary from expected results, please contact Sigma-Aldrich Technical Service for assistance.

Analytical Performance Characteristics

The analytical performance results for the given tests conducted on all target structures, confirm 100% sensitivity, specificity, and repeatability.

Cat. No	Product Description	Target	Intra-assay Specificity	Intra-assay Sensitivity	Inter-assay Specificity	Inter-assay Sensitivity
HHS	Harris Hematoxylin Solution	Nuclear Chromatin	3 of 3	3 of 3	3 of 3	3 of 3

Warnings and Hazards

Refer to Safety Data Sheet and product labeling for any updated risk, hazard or safety information.

HHS16, HHS32, HHS80, HHS128:

Not a hazardous substance or mixture.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Symbol Definitions

Symbols as defined in EN ISO 15223-1:2021

	Manufacturer		Catalogue Number
	Consult Instructions for Use		Batch Code
	Authorized Representative in the European Community/ European Union		European Union Declaration of Conformity (defined in IVDR 2017/746)
	Use-by Date		In vitro diagnostic medical device
	Temperature Limit		Caution
	Date of Manufacture		Importer

References

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, pp 17
- Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
- Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, 3rd ed., LG Luna, Editor, McGraw Hill, New York, 1968
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p129

Contact Information

To place an order, please visit our web site at SigmaAldrich.com. For Technical Service, please visit the tech service page on our web site at SigmaAldrich.com/techservice.

Revision History

Rev 4.0	2016
Rev 5.0	2022
Rev 6.0	2022

Transferred to new template with current branding. Specified for professional use in intended use and precautions. Moved aid to diagnosis statement to intended use. Revised intended use to align with IVDR guidelines. Updated Material Safety Data Sheet to Safety Data Sheet. Updated contact information. Removed instruction to follow CLSI for specimen collection. Removed EN 980 and changed to EN ISO 15223-1:2021 for symbols. Added adverse event contact information. Added Warnings and Hazards.

Sigma-Aldrich, Inc.,
3050 Spruce Street,
St. Louis, MO 63103 USA
an affiliate of Merck KGaA,
Darmstadt, Germany
+1(314) 771-5765

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Merck KGaA,
64271 Darmstadt,
Germany

The Initial M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

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Références

1. Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, pp 17
2. Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
3. Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, 3rd ed., LG Luna, Editor, McGraw Hill, New York, 1968
4. Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p129

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Riferimenti

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, pp 17.
- Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980.
- Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, 3rd ed., LG Luna, Editor, McGraw Hill, New York, 1968.
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p129.

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Referenser

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, pp 17
- Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
- Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, 3rd ed., LG Luna, Editor, McGraw Hill, New York, 1968
- Theory and Practice of Histological Techniques, omarbetad av Bancroft JD och Gamble, M, Churchill Livingstone, New York, 2002, p129

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Referências

- Conn's Biological Stains, 10th ed., RW Horobin, JA Kiernan, Editors, Taylor & Francis, NY, 2002, pp 17
- Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
- Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, 3rd ed., LG Luna, Editor, McGraw Hill, New York, 1968
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p129

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1. Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, pp 17
2. Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
3. Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, 3rd ed., LG Luna, Editor, McGraw Hill, New York, 1968
4. Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p 129

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- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, pp 17
- Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
- Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, 3rd ed., LG Luna, Editor, McGraw Hill, New York, 1968
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, p 129

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Referanser

- Conn's Biological Stains, 10. utg., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, s. 17
- Theory and Practice of Histotechnology, 2. utg., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
- Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, 3. utg., LG Luna, Editor, McGraw Hill, New York, 1968
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, s. 129

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Referanslar

- Conn's Biological Stains, 10th ed., RW Horobin and JA Kiernan, Editors, Taylor & Francis, NY, 2002, s. 17
- Theory and Practice of Histotechnology, 2nd ed., DC Sheehan, BB Hrapchak, Editors, CV Mosby Co., St. Louis, MO, 1980
- Manual of Histologic Staining Methods of the Armed Forces Institute of Pathology, 3rd ed., LG Luna, Editor, McGraw Hill, New York, 1968
- Theory and Practice of Histological Techniques, Edited by Bancroft JD and Gamble, M, Churchill Livingstone, New York, 2002, s. 129

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MDSS GmbH