

Technical Bulletin

Lectin from *Phaseolus vulgaris* (red kidney bean)

Leucoagglutinin PHA-L, lyophilized powder

L2769

Product Description

Lectins are proteins or glycoproteins of non-immune origin that agglutinate cells and/or precipitate complex carbohydrates. Lectins can bind glycoproteins even in the presence of various detergents.¹ The agglutination activity of these highly specific carbohydrate-binding molecules is usually inhibited by a simple monosaccharide. However, for some lectins, di-, tri-, and even polysaccharides are required.

Lectins are isolated from a wide variety of natural sources, including seeds, plant roots and bark, fungi, bacteria, seaweed and sponges, mollusks, fish eggs, body fluids of invertebrates and lower vertebrates, and from mammalian cell membranes. The precise physiological role of lectins in nature is still unknown. However, lectins have proven very valuable in a wide variety of applications *in vitro*, including:

- blood grouping and erythrocyte polyagglutination studies
- mitogenic stimulation of lymphocytes
- lymphocyte subpopulation studies
- fractionation of cells and other particles
- histochemical studies of normal and pathological conditions

A range of lectins suitable for the above applications is available. Most of our lectins are highly purified by affinity chromatography. However, some are offered as purified or partially purified lectins, suitable for specific applications. Many of the lectins are available conjugated to the following moieties, where conjugation does not alter the specificity of the lectin:

- fluorochromes (for detection by fluorimetry)
- enzymes (for enzyme-linked assays)
- insoluble matrices (for use as affinity media)

Please refer to Table 1 for general information on the most common lectins.

Tri- and tetra-antennary glycopeptides containing outer galactose residues and an α -linked mannose residue substituted at positions C-2 and C-6 are specifically bound by PHA-L.² PHA-L also reacts with Tamm-Horsfall glycoprotein, fetuin, and porcine thyroglobulin (weakly). Agglutination activity is inhibited by galactosyl-(1,4)-N-acetylglucosaminyl- α (1,2)-mannose. However, PHA-L does not react with all glycoproteins that contain this sugar.

Lectins from *Phaseolus vulgaris*, have not been proven to be specific for any particular oligosaccharide.^{3,4} Instead, PHA-L recognizes terminal galactose residues of complex glycans on mammalian glycoproteins such as thyroglobulin.

Five *Phaseolus vulgaris* isolectins with similar amino acid compositions can be separated: L4, L3E1, L2E2, L1E3, and E4.⁵

Several theses⁶⁻¹⁰ and dissertations¹¹⁻¹⁶ cite use of product L2769 in their research protocols.

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

Aggregation is thought to occur in the presence of high concentrations of 2-mercaptoethanol.

Preparation Instructions

This lectin is soluble in phosphate buffered saline, pH 7.2 (1 mg/mL).

References

1. Rueben, L., *et al.*, *Biochemistry*, **16(9)**, 1787-1794 (1977).
2. Cummings, R.D., and Kornfeld, S., *J. Biol. Chem.*, **257(19)**, 11230-11234 (1982).
3. Chrispeels, M.J., and Raikhel, N.V., *Plant Cell*, **3(1)**, 1-9 (1991).
4. Schumacher, U. *et al.*, "Lectins as tools in histochemical techniques: a review of methodological aspects", in *Lectin Reviews* (New York), Vol. 1, pp. 195-201 (1991).
5. Leavitt, R.D. *et al.*, *J. Biol. Chem.*, **252(9)**, 2961-2966 (1977).
6. Sprayberry, Kristen M., "History of stress affects cell-mediated immunity". The Pennsylvania State University, B.S. honors thesis, p. 7 (2018).
7. Thompson, Ameeka Louise, "Characterization of Interferon Regulatory Factor Mucosal and peripheral blood responses to vaccination with the licensed live oral typhoid vaccine Ty21a". University of Liverpool, M.Phil. thesis, p. 35 (2011).
8. Côrte-Real, Beatriz Fernandes, "Effect of amino acid substitutions in alpha-gliadin peptides on their immunogenicity in coeliac disease". Universidade de Lisboa, M.Sc. thesis, p. 20 (2015).
9. Harris, Angela, "Characterization of Interferon Regulatory Factor-7 in Defined Subsets of Human Peripheral Blood Mononuclear Cells and Analysis of the Effect of Knockdown on HIV-1 Infection". University of Manitoba, M.Sc. thesis, p. 47 (2017).
10. Rowley, Allison Annette, "Potential Downstream Immunological Effects of Evolved Disease Tolerance in House Finches". Virginia Polytechnic Institute and State University, M.S. thesis, p. 13 (2020).
11. Alizadeth, Hossein, "An Investigation into Bioactive Proteins and Their Changes During Imbibition, Germination and Development of Red Kidney Bean Seeds (*Phaseolus vulgaris* L.)". University of Canterbury, Ph.D. dissertation, pp. 89, 115 (2011).
12. Jabeen, Almas, "Study of the Suppression of Inflammatory Arthritis at Molecular Level by Natural and Synthetic Inhibitors of TNF- α and IL-1 β ". University of Karachi, Ph.D. dissertation, pp. 42, 59 (2013).
13. Ehmke, Manja, "NLRC5: An NLR member with impact on adaptive and innate immune mechanisms and its protective role in the progression of colitis". Christian-Albrechts-Universität zu Kiel, Dr. rer. nat. dissertation, p. 28 (2015).
14. Rump, Katharina, "Einfluss genetischer Varianten im Aquaporin 5-Promotor-Gen auf Schlüsselmechanismen der Sepsis" ("Influence of genetic variants in the aquaporin 5 promoter gene on key mechanisms of sepsis"). Universität Duisberg-Essen, Dr. rer. nat. dissertation, p. 32 (2015).
15. Wong, Raymond Waiman, "Targeting Viral RNA Processing to Control HIV-1 Infection". University of Toronto, Ph.D. dissertation, p. 191 (2018).
16. Al-Abedi, Raheem Abd Zaid Abed, "Metastatic cell responses to anticancer therapy: The possible role of cancer stem cells". Oxford Brookes University, Ph.D. dissertation, p. 65 (2018).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

Table 1. General lectins guide

Lectin	MW (kDa)	Subunits	Specificity		Mitogenic Activity
			Blood Group	Sugar	
<i>Abrus precatorius</i>					+
• Agglutinin	134	4		Gal	
• Abrin A (toxin)	60	2		Gal	
• Abrin B (toxin)	63.8	2($\alpha\beta$)		Gal	
<i>Agarius bisporus</i>	58.5	-	-	β -gal(1 \rightarrow 3)galNAc	+
<i>Anguilla anguilla</i>	40	2	H	α -L-Fuc	
<i>Arachis hypogaea</i>	120	4	T	β -gal(1 \rightarrow 3)galNAc	
<i>Artocarpus integrifolia</i>	42	4	T	α -gal \rightarrow OMe	+
<i>Bandeiraea simplicifolia</i>					
• BS-I	114	4	A, B	α -gal, α -galNAc	
• BS-I-A ₄	114	4	A	α -galNAc	
• BS-I-B ₄	114	4	B	α -gal	
• BS-II	113	4	acq, B, Tk, T	glcNAc	
<i>Bauhinia purpurea</i>	195	4		β -Gal(1 \rightarrow 4)GalNAc	+
<i>Caragana arborescens</i>	60; 120 ^a	2/4	-	GalNAc	
<i>Cicer arietinum</i>	44	2	-	Fetuin	
<i>Codium fragile</i>	60	4	-	GalNAc	
Concanavalin A	102	4	-	α -Man, α -Glc	+
Succinyl-Concanavalin A	51	2	-	α -Man, α -Glc	+ ^b
<i>Cytisus scoparius</i>	-	-	-	GalNAc, Gal	
<i>Datura stramonium</i>	86	2($\alpha\beta$)	-	(GlcNAc) ₂	
<i>Dolichos biflorus</i>	140	4	A ₁	α -GalNAc	
<i>Erythrina corallodendron</i>	60	2	-	β -Gal(1 \rightarrow 4)GlcNAc	+
<i>Erythrina cristagalli</i>	56.8	2($\alpha\beta$)	-	β -Gal(1 \rightarrow 4)GlcNAc	
<i>Euonymus europaeus</i>	166	4($\alpha\beta$)	B,H	α -Gal(1 \rightarrow 3)Gal	+
<i>Galanthus nivalis</i>	52	4	(h)	Non-reduced α -Man	
<i>Glycine max</i>	110	4	-	GalNAc	+ ^c
<i>Helix aspersa</i>	79	-	A	GalNAc	
<i>Helix pomatia</i>	79	6	A	GalNAc	
<i>Lathyrus odoratus</i>	40-43	4($\alpha\beta$)	-	α -Man	+
<i>Lens culinaris</i>	49	2	-	α -Man	+
<i>Limulus polyphemus</i>	400	18	-	NeuNAc	
• Bacterial agglutinin	-	-	-	GalNAc, GlcNAc	
<i>Lycopersicon esculentum</i>	71	-	-	(GlcNAc) ₃	
<i>Maackia amurensis</i>	130	2($\alpha\beta$)	O	Sialic acid	+
<i>Maclura pomifera</i>	40-43	2($\alpha\beta$)		α -Gal, α -GalNAc	
<i>Momordica charantia</i>	115-129	4($\alpha\beta$)	-	Gal, GalNAc	
<i>Naja mocambique mocambique</i>	-	-	-	-	
<i>Naja naja kaouthia</i>	-	-	-	-	

<i>Narcissus pseudonarcissus</i>	26	2	(h)	α -D-Man	
<i>Persea americana</i>	-	-	-	-	
<i>Phaseolus coccineus</i>	112	4	-	-	
<i>Phaseolus limensis</i>	247 (II)	8	A	GalNAc	+
	124 (III)	4			
<i>Phaseolus vulgaris</i>					
PHA-E	128	4	-	Oligosaccharide	+
PHA-L	128	4	-	Oligosaccharide	+
PHA-P					
PHA-M					
<i>Phytolacca americana</i>	32	-	-	(GlcNAc) ₃	+
<i>Pisum sativum</i>	49	4($\alpha\beta$)	-	α -Man	+
<i>Pseudomonas aeruginosa PA-I</i>	13 – 13.7	-	-	Gal	+ ^c
<i>Psophocarpus tetragonolobus</i>	35	1	-	GalNAc, Gal	
<i>Ptilota plumosa</i>	65; 170	-	B	α -Gal	
<i>Ricinus communis</i>					
• Toxin, RCA60	60	2	-	GalNAc, β -Gal	
• Toxin, RCA120	120	4	-	β -Gal	
<i>Sambucus nigra</i>	140	4($\alpha\beta$)	-	α NeuNAC(2→6)Gal, GalNAc	+ ^c
<i>Solanum tuberosum</i>	50, 100 ^a	1,2	-	(GlcNAc) ₃	
<i>Sophora japonica</i>	133	4	A,B	β -GalNAc	
<i>Tetragonolobus purpureas</i>	120 [A]	4	H	α -L-Fuc	
	58 [BA]	2	H	α -L-Fuc	
	117 [C]	4	H	α -L-Fuc	
<i>Triticum vulgare</i>	36	2	-	(GlcNAc) ₂ , NeuNAC	+
<i>Ulex europaeus</i>					
UEA I	68	-	H	α -L-Fuc	
UEA II	68	-			
<i>Vicia faba</i>	50	4($\alpha\beta$)	-	Man, Glc	+
<i>Vicia sativa</i>	40	4($\alpha\beta$)	-	Glc, Man	+
<i>Vicia villosa</i>	139	4	A ₁ +T _n	GalNAc	
• A4	134	4	A ₁	GalNAc	
• B4	143	4	T _n	GalNAc	
<i>Vigna radiata</i>	160	4	-	α -Gal	
<i>Viscum album</i>	115	4($\alpha\beta$)	-	β -Gal	
<i>Wisteria floribunda</i>	68	2	-	GalNAc	

^a Concentration-dependent molecular weight

^c Mitogenic for neuraminidase-treated lymphocyte

^b Non-agglutinating and mitogenic

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

L2769pis Rev 05/22 IRB,MWM,JRC,NSB,SAG,MAP,GCY

**MILLIPORE
SIGMA**