



## STREPTAVIDIN-PEROXIDASE POLYMER LABELED

Product Number **S9420**

Storage Temperature <0°C

CAS #: N/A

### Product Description

Appearance: Brown powder

Molecular weight: 60 kDa for streptavidin<sup>1</sup>

44 kDa for monomeric peroxidase<sup>2</sup>

Streptavidin (from *Streptomyces avidinii*) is used in biochemical research because of its extremely high affinity for biotin, similar to that of avidin. (See product S4762 and related references.<sup>3</sup> Peroxidase from horseradish (HRP) is widely used in protein conjugation because of its comparative stability and selection of substrates for a variety of applications. Making use of the avidin-biotin or streptavidin-biotin complex is often used over other detection systems due to the ease of biotinylating target proteins.

This product ("strep-poly-HRP") is expected to be useful for highly sensitive detection of biotin conjugates, demonstrating at least a 10-fold increase in sensitivity over the nonpolymerized enzyme conjugate, product S5512. Early tests showed use of S9420 at 5 micrograms/mL detected albumin-biotin (product A6043) at 24 ng/mL, compared to a positive response with S5512 at 390 ng/mL.<sup>4</sup> Use of peroxidase-polymer-labeled streptavidin in ELISA in detecting a biotinylated second antibody detected the allergen at 0.6 ng/mL.<sup>5</sup>

The product is a conjugate of streptavidin (product S4762) and a polymerized horseradish peroxidase (from product P0889). The streptavidin was activated using S-acetylthioglycolic acid N-hydroxysuccinimide ester (SATA, product A9043),<sup>6</sup> then reacted with maleimide-activated poly-peroxidase. The polymerization of HRP produced approximately 20 molecules per polymer, as determined by LALLS and gel filtration chromatography. The conjugate was purified by gel filtration chromatography and lyophilized from a sodium citrate buffer containing trehalose.<sup>4</sup>

### Components

The product is approximately 50% protein (BCA assay), with the balance being primarily trehalose and sodium citrate.<sup>4</sup> Note that native peroxidase is approximately 75% protein, 25% carbohydrate and heme group.<sup>2</sup>

### Preparation Instructions

## Product Information

The strep-poly-HRP dissolves at 0.25 mg/mL in phosphate buffered saline pH 6.0 to give a clear to light tan solution. Stock solutions should be stable at least two weeks at 2-8°C, but it is recommended to freeze aliquots at -20°C. Solutions in 50% glycerol should be stable for at least a year at -20°C.<sup>4</sup>

By analogy to usage of S5512, further dilutions should be made with PBS containing BSA and Tween 80 as blocking agents to reduce non-specific background.<sup>7,8</sup>

### Storage/Stability

The product as shipped should be stored dry at -20°C, retaining at least 95% initial activity for a year.<sup>4</sup>

### Procedure

The usage will be much like that of streptavidin-peroxidase S5512, but with greater sensitivity in detection of biotinylated molecules. Determining a specific titer will depend on the customer's application. However, for product S5512, a first dilution was suggested as 1:100 to 1:1000, followed by incubation at room temperature for an hour, then rinsing thoroughly with buffer before applying HRP substrate.<sup>7,8</sup> Given the 10-fold improvement in sensitivity, we suggest a stock solution of 1 mg/mL, with 5 µg/mL as an effective working concentration, i.e., 1:1000 or 1:10,000 as a first dilution to try.<sup>4</sup>

### Product Profile

Streptavidin activity: 1 to 4 units/mg protein

Unit definition: One unit will bind 1 µg biotin.<sup>9</sup>

Peroxidase activity: 100 to 200 units/mg protein

Unit definition: One unit will form 1 mg purpurogallin in 20 sec from pyrogallol at pH 6.0 at 20°C.

### References

1. Haeuptle, M.-T. et al., *J. Biol. Chem.*, 258, 305 (1983).
2. Welinder, K.J., *Eur. J. Biochem.*, 96, 483 (1979).
3. *Methods in Enzymology*, 184 (1990). Whole volume: "Avidin-Biotin Technology."
4. Dapron, J. and Quinn, T., Sigma production; published as conference poster session.
5. Sander, I. et al., *J. Immunol. Meth.*, 210, 93-101 (1997).
6. Julian, R., *Anal. Biochem.*, 132, 68 (1983).

7. LaRochelle, W.J. and Froehner, S.C., Meth. Enzymol., 184, 433 (1990).

8. Brakel, C.L., Brower, M.S. and Garry, K., Meth. Enzymol., 437 (1990).

9. Green, N.M., Meth. Enzymol., 18A, 418 (1970).

CKV 12/16/98

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.