


# Purification with Millipore MultiScreen<sub>384</sub> SEQ Plates on Beckman-Coulter's Biomek FX for High Quality Direct Sequencing in High-Throughput Genomics

Wei-Hsien Yang<sup>1</sup>, Nathan McDonald<sup>1,3</sup>, Jack Leonard<sup>2</sup>, Christopher Barbagallo<sup>2</sup>, and Hans Fuernkranz<sup>1,3</sup>  
<sup>1</sup>DNA Sciences, Inc., <sup>2</sup>Millipore Corporation, <sup>3</sup>Applied Biosystems



- Throughput: ~3 384-well plates / hour
- Currently using 2 Vacuum Manifolds on one FX Deck
  - Throughput can be significantly increased by using 4 or more manifolds.
- 96-Well Plates get run on MegaBACE 1000

**Benefits of High Quality Sequencing Data**

**Polymorphism Discovery**

- Polymorphism Detection at phred 50 is significantly more accurate and automatable than at phred 30.
- CG polymorphism shown here at base 478 of consensus sequence (yellow tag) was called correctly by software.

**Accuracy of Results**

- Fewer false positives (less automated mis-calls)
- Fewer questionable calls, i.e. less human interference required.
- Data confidence is higher:
  - Phred 30 = 1 mistake in 1,000 bases
  - Phred 40 = 1 mistake in 10,000 bases
  - Phred 50 = 1 mistake in 100,000 bases

Significant cost savings - fewer data analysts are needed

**Problem Detection**

- Easier to spot sequencing peculiarities

- Overall Conclusions**
- The MultiScreen<sub>384</sub>SEQ Clean-up procedure has allowed *DNA Sciences, Inc.* to detect polymorphisms in disease associated genes faster, cheaper, and with greater data confidence.
  - The MultiScreen<sub>384</sub>SEQ set-up on Biomek FX provides a high level of walk-away confidence.
  - Excellent recovery of products from MultiScreen<sub>384</sub>SEQ Plates can allow further miniaturization of sequencing reactions, which leads to more cost savings.
  - To further increase throughput, a Biomek FX with dual pod, 384/96 head design can be implemented.


## Overview

The authors have developed a post-extension clean-up protocol for sequencing using *Millipore's* MultiScreen<sub>384</sub> SEQ plates that is high-throughput, cost effective, amenable to automation and results with high consistency in outstanding sequencing quality scores (phred >= 35 for approx. 90% of samples).

## DNA Sciences, Inc.

a genetics discovery company focused on identifying the genetic basis of disease susceptibility and response to drug treatment.

- Sequencing
  - High-Quality Sequencing for polymorphism discovery, incl.
  - Heterozygote detection (DYEamic ET Dye Terminator on MegaBACE)
- SNP Genotyping
  - High-Throughput multiplexed SNP genotyping (Single Base Extension, SNaPshot on ABI3700)
- STR Genotyping
  - Linkage analysis studies, (LMS HD5 [ABI] on MegaBACE)



*DNA Sciences, Inc.* utilizes a 'Generic Interchangeable Module System (GIMS)' that allows for rapid changes in project type (sequencing or genotyping) without production delays.

*Millipore's* MultiScreen<sub>384</sub>SEQ Plates setup on *Beckman-Coulter's* Biomek FX is an integral part of *DNA Sciences' 'GIMS'* that has produced sequencing quality scores of phred 50.


•The size-exclusion method efficiently removes salts and unincorporated dye terminators through vacuum filtration.

•No sephadex gel.

•No centrifugation.

•MultiScreen<sub>384</sub>SEQ Plates easily adaptable to liquid-handling instruments.

•Consistent high quality results.



## Comparison Study of Ethanol and MultiScreen<sub>384</sub> SEQ Clean-up

**Size of the study:** Forward and Reverse sequences of 18 regions on 96 genomic samples → ~ 3600 samples.

•Prepared 3 sets of the above samples:

**Set 1 =** Half reaction sequences for Ethanol Clean-up

**Set 2 =** Quarter reaction sequences for Ethanol Clean-up

**Set 3 =** Quarter reaction sequences for MultiScreen<sub>384</sub> SEQ Clean-up

Compared the results of Average Phred, > Phred 30 Read length, Pass Rate, and Number of resequenced plates.

Set 1				Set 2				Set 3			
ROID	Avg Phred	Avg Phred 30 RL	Pass Rate	ROID	Avg Phred	Avg Phred 30 RL	Pass Rate	ROID	Avg Phred	Avg Phred 30 RL	Pass Rate
hmar1_4.A.ya	37.92	367.09	97%	hmar1_4.A.ya	40.53	386.33	93%	hmar1_4.A.ya	41.95	413.66	98%
hmar1_6.A.ya	39.35	400.37	98%	hmar1_6.A.ya	41.45	419.22	98%	hmar1_6.A.ya	39.51	395.27	95%
hmar1_8.A.ya	39.47	372.66	93%	hmar1_8.A.ya	42.42	373.77	96%	hmar1_8.A.ya	42.51	373.61	91%
hmar1_9.A.ya	35.87	302.61	93%	hmar1_9.A.ya	39.44	339.58	91%	hmar1_9.A.ya	38.88	354.38	97%
hmar1_10.A.ya	37.21	338.91	95%	hmar1_10.A.ya	37.63	301.23	91%	hmar1_10.A.ya	39.36	379.63	97%
hmar1_2.A.ya	35.83	306.95	88%	hmar1_2.A.ya	39.15	373.14	93%	hmar1_2.A.ya	38.61	380.21	92%
hmar2_2.A.ya	37.26	374.83	99%	hmar2_2.A.ya	39.35	414.91	91%	hmar2_2.A.ya	37.95	385.32	98%
hmar2_3.A.ya	38.51	376.88	98%	hmar2_3.A.ya	40.84	408.05	94%	hmar2_3.A.ya	40.48	430.64	98%
hmar2_3.ya	41.86	410.79	97%	hmar2_3.ya	42.31	408.01	92%	hmar2_3.ya	43.09	445.2	98%
hmar2_1.A.ya	39.38	383.49	95%	hmar2_1.A.ya	35.12	298.71	76%	hmar2_1.A.ya	38.51	352.84	90%
hmar2_2.ya	39.24	332.17	95%	hmar2_2.ya	40.1	338.87	89%	hmar2_2.ya	34.96	270.93	88%
hmar2_1.A.ya	40.2	349.38	95%	hmar2_1.A.ya	38.3	295.85	78%	hmar2_1.A.ya	35.97	277.05	83%
h12a_1.ya	38.89	336.82	98%	h12a_1.ya	37.48	371.95	88%	h12a_1.ya	37.44	365.39	93%
h12a_1.ya	39.78	400.66	98%	h12a_1.ya	41.14	383.16	88%	h12a_1.ya	44.93	452.26	98%
h12a_7.A.ya.r1	13.65	0	0%	h12a_7.A.ya.r1	0	0	0%	h12a_7.A.ya.r1	49.76	411.51	97%
h12a_7.A.ya.r2	0	0	0%	h12a_7.A.ya.r2	0	0	0%	h12a_7.A.ya.r2	no resequencing necessary		
h12a_7.A.ya.r3	41.49	390.22	96%	h12a_7.A.ya.r3	40.01	382.05	97%	h12a_7.A.ya.r3	no resequencing necessary		
h12a_7.A.ya.r1	38.97	383.58	97%	h12a_7.A.ya.r1	40.55	400.58	91%	h12a_7.A.ya.r1	40.03	403.13	95%
h12a_7.A.ya.r2	42.17	436.31	97%	h12a_7.A.ya.r2	34.63	293.07	68%	h12a_7.A.ya.r2	no resequencing necessary		
h12a_7.A.ya.r3	42.92	433.2	97%	h12a_7.A.ya.r3	40.55	409.43	91%	h12a_7.A.ya.r3	no resequencing necessary		
h10_4a.ya	38.64	333.96	93%	h10_4a.ya	41.28	359.4	95%	h10_4a.ya	39.21	360.73	98%
h10_4a.ya	41.31	386.95	91%	h10_4a.ya	43.99	372.77	89%	h10_4a.ya	40.81	386.13	95%
hmag_10.B.ya	35.44	329.11	91%	hmag_10.B.ya	33.29	301.25	92%	hmag_10.B.ya	35.82	335.52	93%
hmag_10.B.ya	39.33	382.28	93%	hmag_10.B.ya	38.64	378.15	97%	hmag_10.B.ya	41.93	422.06	97%
hmag_8.ya	38.52	354.95	96%	hmag_8.ya	39.57	368.2	97%	hmag_8.ya	40.43	381.9	93%
hmag_8.ya	40.85	375.01	97%	hmag_8.ya	40.35	395.19	98%	hmag_8.ya	37.93	357.41	97%
hmag_12.A.ya.r1	41.38	350.53	96%	hmag_12.A.ya.r1	39.56	326.16	89%	hmag_12.A.ya.r1	41.11	368.71	98%
hmag_12.A.ya.r1	41.49	358.03	83%	hmag_12.A.ya.r1	41.38	342.22	88%	hmag_12.A.ya.r1	42.69	393.69	89%
hmag_12.A.ya.r2	38.85	302.66	94%	hmag_12.A.ya.r2	37.92	303.81	97%	hmag_12.A.ya.r2	39.44	314.03	91%
hmag_12.A.ya.r3	42.53	332.03	98%	hmag_12.A.ya.r3	43.03	333.07	92%	hmag_12.A.ya.r3	42.83	357.21	92%
hmag_12.A.ya.r1	41.1	404.63	99%	hmag_12.A.ya.r1	39.49	408.77	94%	hmag_12.A.ya.r1	no resequencing necessary		
hmag_12.A.ya.r2	38.82	372.05	98%	hmag_12.A.ya.r2	38.69	387.47	98%	hmag_12.A.ya.r2	35.67	335.67	95%
hmag_2.ya	36.96	239.65	95%	hmag_2.ya	32.34	184.32	88%	hmag_2.ya	31.25	181.96	83%
hmag_2.ya	38.94	249.41	94%	hmag_2.ya	32.5	186.19	88%	hmag_2.ya	32.04	252.27	93%
hmf_pro_42.ya	41.45	370.51	91%	hmf_pro_42.ya	42.97	413.47	97%	hmf_pro_42.ya	39.82	364.73	98%
hmf_pro_42.ya	29.84	229.71	86%	hmf_pro_42.ya	32.5	186.19	88%	hmf_pro_42.ya	32.04	252.27	93%
hmf_pro_42.ya	35.72	223.55	92%	hmf_pro_42.ya	37.88	283.87	97%	hmf_pro_42.ya	40.47	322.09	97%
hmf_pro_42.ya	39.18	299.04	94%	hmf_pro_42.ya	38.74	286.11	96%	hmf_pro_42.ya	35.27	271.95	97%
hmf_pro_42.ya.r1	34.87	359.52	98%	hmf_pro_42.ya.r1	35.97	353.19	98%	hmf_pro_42.ya.r1	35.40	338.16	95%
hmf_pro_42.ya.r2	35.41	359.24	98%	hmf_pro_42.ya.r2	36.45	357.45	98%	hmf_pro_42.ya.r2	no resequencing necessary		
hmf_pro_42.ya.r1	29.79	174.48	36%	hmf_pro_42.ya.r1	32.57	193.73	87%	hmf_pro_42.ya.r1	32.73	282.93	92%
hmf_pro_42.ya.r2	29.41	117.87	80%	hmf_pro_42.ya.r2	27.02	62.57	8%	hmf_pro_42.ya.r2	no resequencing necessary		
hmf_pro_42.ya.r3	33.95	382.83	98%	hmf_pro_42.ya.r3	36.48	346.12	96%	hmf_pro_42.ya.r3	no resequencing necessary		
<b>Averages:</b>	<b>37</b>	<b>323</b>	<b>87%</b>	<b>Averages:</b>	<b>37</b>	<b>321</b>	<b>83%</b>	<b>Averages:</b>	<b>39</b>	<b>353</b>	<b>94%</b>

## Summary of Results

	1/2 rxn ethanol	1/4 rxn ethanol	1/4 rxn MultiScreen <sub>384</sub> SEQ
<b>Average Phred</b>	37	37	39
<b>Average &gt; Phred 30 Readlength</b>	323	321	353
<b>% Pass Rate*</b>	87%	83%	94%
<b># of Resequenced Plates**</b>	7	8	0

\*Pass rate = number of samples with Phred score > 30 and at least 100 bp of Phred 30 Read length / total number of samples

\*\*Sample plates resequenced due to < 80% pass rate or contamination of samples.

## Comparison Study Conclusions


MultiScreen<sub>384</sub>SEQ Plates Clean-up procedure resulted in

- Higher Phred score
- Longer good quality sequence
- Higher sample pass rate
- No cross contamination of samples
- No resequences needed

Ethanol Clean-up procedure, although produces satisfactory results,

- Lacks consistency in performance, especially among different users.
- Has higher occurrences of cross-contamination → more resequences needed.

## Contamination Issue with Ethanol Clean-up Procedure



Cross contamination occurs after inversion of ethanol and extension products mixture in sealed plate during ethanol clean-up procedure.

This is demonstrated in the picture below where the "pink wells" are the results of crossing over of red food coloring + ethanol mixture from source wells into the adjacent wells during the inversion process.

Phred 40	Phred 50	Ave. Phred	Genotype
334	238	42.7	CG
328	230	41.4	CG
345	250	43.4	CG
111	46	35.0	CG
342	234	43.2	CC
292	157	40.2	CC
249	238	42.7	GG
377	301	43.4	GG

Sample of MultiScreen<sub>384</sub>SEQ Phred 50 sequences for Polymorphism Detection

