



Product Information

Z-Phe-Phe Fluoromethyl Ketone

Product Number **C 9109**

Storage Temperature $-20\text{ }^{\circ}\text{C}$

Product Description

Molecular formula: $\text{C}_{27}\text{H}_{27}\text{FN}_2\text{O}_4$

Molecular weight: 462.5

Z-Phe-Phe Fluoromethyl Ketone (Z-FF-FMK) is an irreversible inhibitor of cathepsin B and cathepsin L. (Z = benzyloxycarbonyl)

Among the intracellular proteinases, thiol proteinases such as cathepsin B and cathepsin L, active at slightly acidic pH, play a role in lysosomal protein catabolism. Both have identical substrate specificity but 4 M urea at pH 5.0 will inactivate cathepsin B, while cathepsin L retains its activity under these conditions.¹

Cathepsin B is a lysosomal cysteine protease involved in cellular protein turnover.² Cathepsin B has a high abundance and exhibits both endopeptidase and peptidyl dipeptidase activity.^{3,4} This enzyme has been implicated in several pathological conditions including arthritis and tumor metastasis.

Cathepsin L is the most powerful of the lysosomal proteinases and has a higher specific activity than cathepsin B in the degradation of physiological protein substrates. Cathepsin L, but not cathepsin B, can also generate kinins from high and low molecular weight kininogens in vitro.⁵

FMK is a trapping group responsible for irreversible inhibition but is non-cytotoxic. Inhibition occurs when the the FMK group covalently bonds to the $-\text{SH}$ of an adjacent cysteine residue on the target protein.

Preparation Instructions

Prepare stock 20 mM solutions in dry ($\geq 99.9\%$) DMSO to maintain product stability. Also soluble in DMF.

Effective final concentration is estimated to be $50\text{ }\mu\text{M}$.

Storage/Stability

Store product at $-20\text{ }^{\circ}\text{C}$. It is reported to be stable at room temperature for one year in a desiccator.

Store stock solutions at $-20\text{ }^{\circ}\text{C}$ for 6-8 months. Allowing the material to warm to room temperature before use to ensure stability.

References

1. Kamboj, R. C., et al. A selective colorimetric assay for cathepsin L using Z-Phe-Arg-4-methoxy-beta-naphthylamide. *Biochimie*, **75**, 873-878 (1993).
2. Kirschke H, Barrett AJ. 1987. Lysosomal cysteine proteases. In: Glaumann H, Ballard FJ, eds. *Lysosomes: Their role in protein breakdown*. London: Academic Press. pp 193-238.
3. Aronson, N. N., Jr., and Barrett, A. J., The specificity of cathepsin B. Hydrolysis of glucagon at the C-terminus by a peptidyl dipeptidase mechanism. *Biochem. J.*, **171**, 759-765 (1978).
4. Fosang, A. J., et al., The interglobular domain of cartilage aggrecan is cleaved by PUMP, gelatinases, and cathepsin B. *J. Biol. Chem.*, **267**, 19470-19474 (1992).
5. Desmazes, C., et al., Cathepsin L, but not cathepsin B, is a potential kininogenase. *Biol. Chem.*, **382**, 811-815 (2001).

DLV/jwm 11/01

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.