

Supel™ Tox SPE Cartridges

Fast and Simple Cleanup for Mycotoxin Analysis



Supel™ Tox SPE Cartridges

Specialty Products for Mycotoxin Analysis

Features and Benefits

- Remove interferences associated with mycotoxin analysis
- Better reproducibility than the current industry standard, immunoaffinity columns
- Time associated with sample preparation is up to ten times less than that of immunoaffinity columns
- Basic and quick methodology requiring little additional method development
- Improved shelf life over immunoaffinity columns due to thermally stable format. No refrigeration required for shipping and storage.

The need for a quick, simplistic sample cleanup approach prior to chromatographic mycotoxin analysis has brought about SPE cartridges that significantly decrease sample prep time, increase reproducibility, and are more user friendly as compared to the industry standard immunoaffinity columns. In addition, the Supel Tox SPE approach requires less equipment and fewer consumables, providing an additional cost savings.

Unlike the multiple step “bind and elute” strategy implemented when using immunoaffinity columns, the Supel Tox AflaZea, DON and Tricho SPE cartridges employ an “interference removal” strategy



which saves time by eliminating wash steps prior to elution of aflatoxin and zearalenone, deoxynivalenol, and tricothecenes (type A and B), respectively. Cartridges removing interferences associated with analysis of fumonisins (B₁ and B₂) and ochratoxin A are also available as a part of our Supel Tox product offering. Refer to the table below for further details.

Sample Purification Strategy for the Various Supel Tox SPE Cartridges


SPE Cartridge	Analyte(s)	Matrix	Purification Strategy
Supel Tox AflaZea	Aflatoxin and Zearalenone	Grains, feeds, TMR samples, peanuts, peanut products and aqueous solutions	Interference removal
Supel Tox DON	Deoxynivalenol (DON)	Wheat, flour and corn	Interference removal
Supel Tox Tricho	Tricothecenes (Type A and B)	Grains and complex matrices	Interference removal
Supel Tox TrichoBind	Tricothecenes (Type A and B)	Grains, feeds and other complex matrices	Multifunctional bind and elute
Supel Tox FumoniBind	Fumonisins (B ₁ and B ₂)	Grains and cereals	Multifunctional bind and elute
Supel Tox OchraBind	Ochratoxin A	Grain and feed samples	Multifunctional bind and elute

Comparison of “Bind and Elute” Principle to “Interference Removal” Principle


In the “bind and elute” strategy, the analyte of interest becomes bound to the SPE sorbent. The unwanted matrix is removed, and the analyte is ultimately eluted. In “interference removal,” the unwanted matrix components remain bound to the sorbent, while the purified

analyte is eluted. This strategy reduces the number of purification steps, allowing for time savings and increase in sample throughput. The Supel™ Tox AflaZea, DON and Tricho SPE cartridges employ the “interference removal” strategy.

“Bind and Elute” Principle



“Interference Removal” Principle




Comparison of Supel™ Tox SPE Cartridges to Immunoaffinity Columns for Matrix Removal Prior to Mycotoxin Analysis

- Ten times faster – ten times sample throughput
- Vast improvement in reproducibility

The following is an experiment comparing the use of the Supel Tox AflaZea SPE cartridge to a leading immunoaffinity column for the analysis of aflatoxins B₁, B₂, G₁ and G₂ in peanut paste. The results illustrate Supel Tox AflaZea SPE cartridges are superior to immunoaffinity columns in terms of process simplicity, time required for sample preparation, and control of variation.

Cleanup of Aflatoxins in Peanut Paste: Immunoaffinity Versus Supel™ Tox (n = 3)




Comparison of Sample Preparation Methods Using Supel Tox AflaZea SPE vs Immunoaffinity for Aflatoxins in Peanut Paste

	Immunoaffinity	Supel Tox AflaZea SPE Cartridge
Sample Prep Time (post-extraction to pre-analysis)	<ul style="list-style-type: none"> • 60 minutes • 8 samples/day (if processing 1 at a time) 	<ul style="list-style-type: none"> • 6 minutes • 80 samples/day (if processing 1 at a time)
Ease of Use	<ul style="list-style-type: none"> • Large volumes of liquid • Controlled drop rates • Numerous complicated steps • Additional buffer salts required • Must be refrigerated, brought to room temp before use 	<ul style="list-style-type: none"> • Small volumes of liquid • Vacuum filtration used • Steps few and not complicated • No additional reagents required • Column does not require special storage conditions
Procedure (post-extraction to analysis)	<p>Stage 1 (15 minutes)</p> <ol style="list-style-type: none"> 1. Configure manifold for waste collection 2. Add 1 mL sample to 17 mL of phosphate buffered saline, vortex 3. Uncap/mount/drain cartridges by gravity 4. Apply reservoirs, load sample onto cartridges <p>Stage 2 (15 minutes)</p> <ol style="list-style-type: none"> 5. Rinse interferences 6. Reconfigure manifold for sample collection 7. Elute/collect sample <p>Stage 3 (30 minutes)</p> <ol style="list-style-type: none"> 8. Evaporate sample to dryness 9. Reconstitute sample and vortex 10. Transfer 0.2 mL sample to vial 11. Dilute sample and vortex <p>Analysis</p>	<p>Purify and Transfer (6 minutes)</p> <ol style="list-style-type: none"> 1. Configure manifold for sample collection 2. Mount cartridges 3. Load 2 mL sample 4. Elute and collect under vacuum 5. Transfer 0.2 mL sample to vial 6. Dilute sample and vortex <p>Analysis</p>

Analysis of Aflatoxin B₁, B₂, G₁ and G₂ in Corn Meal, Wheat Meal and Raw Peanut Paste Matrices


The Supel Tox AflaZea SPE cartridge was used to demonstrate the cleanup of corn meal, wheat meal and raw peanut paste matrices prior to the analysis of Aflatoxin B₁, B₂, G₁ and G₂. In the procedure, corn meal, wheat meal and raw peanut paste extracts were spiked with Aflatoxin Standard Mix (Cat. No. 34036) and run through the Supel Tox AflaZea SPE cartridges. The eluate was collected and diluted for HPLC analysis. Results illustrate Supel Tox AflaZea SPE cartridges provide exceptional cleanup of corn meal, wheat meal and raw peanut paste matrices with high analyte recoveries and low %RSD for aflatoxin analysis.

Analyte Recovery and %RSD of Aflatoxins from Corn, Wheat and Peanut Paste (n = 3)



Aflatoxins in Corn Meal, Wheat Meal and Raw Peanut Paste After Supel Tox Cleanup

column: Ascentis® Express C18, 10 cm x 2.1 mm I.D., 2.7 μm (53823-U)
 mobile phase: (A) water; (B) acetonitrile; (C) methanol; (700:120:120,A:B:C) with 0.780 g potassium bromide and 230 μL of nitric acid
 flow rate: 0.400 mL/min
 column temp.: 35 °C
 detector: Fluorescence detector, Excitation: 360 nm, Emission: 440 nm
 injection: 40 μL




Analysis of Corn Meal and Wheat Meal Matrices for Deoxynivalenol

Deoxynivalenol (DON), also referred to as vomitoxin, is a concern primarily in the livestock industry, as it is known to cause a reduction in food intake in farm animals. DON is a type B trichothecene which may occur in a variety of grains.

Cleanup of corn meal and wheat meal matrices for the analysis of DON was illustrated with the use of Supel™ Tox DON SPE cartridges. The grain extracts were spiked with deoxynivalenol standard (Cat. No. 34124) and applied to the cartridges. The eluate was collected, evaporated, and reconstituted for HPLC analysis. Blank unspiked matrix samples were also run through cartridges and analyzed. For the recovery calculations, the peak area of DON observed in the blanks was subtracted from that of the spiked samples.



The corresponding bar chart and chromatograms illustrate that in addition to effectively removing the unwanted corn and wheat meal matrices, Supel™ Tox DON SPE cartridges also provide high analyte recoveries and reproducibility for the analysis of DON.

Analyte Recovery and %RSD of Deoxynivalenol (DON) from Corn and Wheat (n = 3)



Deoxynivalenol in Corn and Wheat Meal After Supel Tox Cleanup

column: Ascentis® Express C18, 10 cm x 2.1 mm I.D., 2.7 µm (53823-U)
mobile phase: (A) water; (B) acetonitrile; (C) methanol; (92:4:4, A:B:C)
flow rate: 0.400 mL/min
column temp.: 35 °C
detector: UV 220 nm; 360 nm
injection: 40 µL



Supel™ Tox Product Offering

Supel Tox SPE Cartridges

Description	Qty.	Cat. No.
Supel Tox AflaZea SPE Cartridge, 6 mL	30	55314-U
Supel Tox DON SPE Cartridge, 6 mL	30	55316-U
Supel Tox Tricho SPE Cartridge, 6 mL	30	55308-U
Supel Tox TrichoBind SPE Cartridge, LRC	25	55307-U
Supel Tox FumoniBind SPE Cartridge, LRC	25	55315-U
Supel Tox OchraBind SPE Cartridge, LRC	25	55318-U

Featured and Related Products

Description	Qty.	Cat. No.
Ascentis Express HPLC Columns (2.7 µm particles)		
C18, 10 cm x 2.1 mm I.D.	1	53823-U
C18, 15 cm x 2.1 mm I.D.	1	53825-U
RP-Amide, 10 cm x 2.1 mm I.D.	1	53913-U
RP-Amide, 15 cm x 2.1 mm I.D.	1	53914-U
F5, 10 cm x 2.1 mm I.D.	1	53569-U
F5, 15 cm x 2.1 mm I.D.	1	53571-U
Ascentis Express Guard Columns (2.7 µm particles)		
Holder	1	53500-U
C18 Cartridge, 5 mm x 2.1 mm I.D.	3	53501-U
RP-Amide Cartridge, 5 mm x 2.1 mm I.D.	3	53514-U
F5 Cartridge, 5 mm x 2.1 mm I.D.	3	53594-U
Laboratory Equipment		
Coors™ Büchner Funnel, capacity 186 mL	1	Z247324
Individual Filter Adapters, size 3	1	Z254274
Visiprep™ Solid Phase Extraction Manifold, 12-port Model	1	57030-U
Visiprep Solid Phase Extraction Manifold, 24-port Model	1	57250-U
Visiprep DL Solid Phase Extraction Manifold, 12-port Model	1	57044
Visiprep DL Solid Phase Extraction Manifold, 24-port Model	1	57265
Disposable Valve Liners, PTFE	100	57059
Visi™-1 Single SPE Tube Processor	1	57080-U
SPE Tube Adapters for 1, 3, 6 mL tubes	12	57020-U
Mycotoxin Reference Materials (Neats)		
3-Acetyldeoxynivalenol	5 mg	32927
15-Acetyldeoxynivalenol	5 mg	32928
Aflatoxin B ₁	5 mg	32754*
Aflatoxin B ₂	5 mg	32755
Aflatoxin G ₁	5 mg	32756
Aflatoxin G ₂	5 mg	32757*
Deoxynivalenol	5 mg	32943
Fumonisin B ₁	5 mg	32936
Fusarenon X	5 mg	33438
Ochratoxin A	5 mg	32937
T-2 Toxin	5 mg	33947*
Zearalenone	5 mg	32939
Mycotoxin Single Component Standard Solutions		
3-Acetyldeoxynivalenol, 100 µg/mL in acetonitrile	2 mL	34132
15-Acetyldeoxynivalenol, 100 µg/mL in acetonitrile	2 mL	34133
Aflatoxin B ₁ , 2 µg/mL in acetonitrile	2 mL	34029
Aflatoxin B ₂ , 0.5 µg/mL in acetonitrile	2 mL	34034
Aflatoxin G ₁ , 2 µg/mL in acetonitrile	2 mL	34032
Aflatoxin G ₂ , 0.5 µg/mL in acetonitrile	2 mL	34033
Deoxynivalenol, 100 µg/mL in acetonitrile	2 mL	34124
Fumonisin B ₁ , 50 µg/mL in acetonitrile/water	2 mL	34139
Fumonisin B ₂ , 50 µg/mL in acetonitrile/water	2 mL	34142
Fumonisin B ₃ , 50 µg/mL in acetonitrile/water	1 mL	32606
Fusarenon X, 100 µg/mL in acetonitrile	2 mL	34130

*Not available in the U.S.A.

Description	Qty.	Cat. No.
Nivalenol, 100 µg/mL in acetonitrile	2 mL	34131
Ochratoxin A, 10 µg/mL in acetonitrile	2 mL	34037
T-2 Toxin, 100 µg/mL in acetonitrile	2 mL	34071*
Zearalenone, 50 µg/mL in acetonitrile	1 mL	46916-U
Zearalenone, 100 µg/mL in acetonitrile	2 mL	34126
Mycotoxin Mixture Standard Solutions for Multi-analyte Detection		
Aflatoxin Mix 4:	2 mL	34036
B ₁ , 2 µg/mL in acetonitrile	1 mL	
G ₁ , 2 µg/mL in acetonitrile		
B ₂ , 0.5 µg/mL in acetonitrile		
G ₂ , 0.5 µg/mL in acetonitrile		
Aflatoxin Mix 4:	2 mL	33415
B ₁ , 20 µg/mL in acetonitrile		
G ₁ , 20 µg/mL in acetonitrile		
B ₂ , 20 µg/mL in acetonitrile		
G ₂ , 20 µg/mL in acetonitrile		
Japanese Aflatoxin Mix:	5 x 1 mL	40139-U
B ₁ , 25 µg/mL in acetonitrile		
G ₁ , 25 µg/mL in acetonitrile		
B ₂ , 25 µg/mL in acetonitrile		
G ₂ , 25 µg/mL in acetonitrile		
Fumonisin Mix:	2 mL	34143
FB ₁ , 50 µg/mL in acetonitrile/water		
FB ₂ , 50 µg/mL in acetonitrile/water		
B-Trichothecene Mix:	2 mL	34134
DON, 100 µg/mL in acetonitrile		
NIV, 100 µg/mL in acetonitrile		
3-AcDON, 100 µg/mL in acetonitrile		
15-AcDON, 100 µg/mL in acetonitrile		
Certified Matrix Reference Materials for Mycotoxin Analysis (CRMs)		
Compound Feed (Aflatoxin blank)	50 g	BCR375
Aflatoxin B ₁ , <1 µg/kg		
Compound Feedingstuff (Aflatoxins, very low level)	2 x 75 g	ERMBC375
Aflatoxin B ₁ , 2.6 µg/kg		
Aflatoxin B ₂ , 0.2 µg/kg		
Aflatoxin G ₁ , 0.4 µg/kg		
Aflatoxin G ₂ , <0.2 µg/kg		
Compound Feedingstuff (Aflatoxins, high level)	2 x 75 g	ERMBC376
Aflatoxin B ₁ , 12.9 µg/kg		
Aflatoxin B ₂ , 0.68 µg/kg		
Aflatoxin G ₁ , 0.5.2 µg/kg		
Maize (Zearalenone, very low level)	60 g	ERMBC716
Zearalenone, <5 µg/kg		
Maize (Zearalenone, low level)	60 g	ERMBC717
Zearalenone, 83 µg/kg		
Maize Flour (Deoxynivalenol, blank)	150 g	BCR377
Deoxynivalenol, <50 µg/kg		
Peanut Butter (Aflatoxins, low level)	1 ea.	BCR401R
Aflatoxin B ₁ , <0.2 µg/kg		
Aflatoxin B ₂ , <0.2 µg/kg		
Aflatoxin G ₁ , <0.2 µg/kg		
Aflatoxin G ₂ , <0.2 µg/kg		
Wheat (Ochratoxin A, blank)	55 g	BCR471
Ochratoxin A, <0.6 µg/kg		
Mycotoxin Solutions (CRMs)		
Aflatoxin B ₁ , 3.79 µg/g in acetonitrile	4 mL	ERMAC057
Aflatoxin B ₂ , 3.80 µg/g in acetonitrile	4 mL	ERMAC058
Aflatoxin G ₁ , 3.78 µg/g in acetonitrile	4 mL	ERMAC059
Aflatoxin G ₂ , 3.80 µg/g in acetonitrile	4 mL	ERMAC060
4-Deoxynivalenol, 25.1 µg/g in acetonitrile	1 ea.	IRMM315
Nivalenol, 24.0 µg/g in acetonitrile	1 ea.	IRMM316
Zearalenone, 9.95 µg/g in acetonitrile	4 mL	ERMAC699

Sigma-Aldrich® Worldwide Offices

Argentina

Free Tel: 0810 888 7446
Tel: (+54) 11 4556 1472
Fax: (+54) 11 4552 1698

Australia

Free Tel: 1800 800 097
Free Fax: 1800 800 096
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Austria

Tel: (+43) 1 605 81 10
Fax: (+43) 1 605 81 20

Belgium

Tel: (+32) 3 899 13 01
Fax: (+32) 3 899 13 11

Brazil

Free Tel: 0800 701 7425
Tel: (+55) 11 3732 3100
Fax: (+55) 11 5522 9895

Canada

Free Tel: 1800 565 1400
Free Fax: 1800 265 3858
Tel: (+1) 905 829 9500
Fax: (+1) 905 829 9292

Chile

Tel: (+56) 2 495 7395
Fax: (+56) 2 495 7396

People's Republic of China

Free Tel: 800 819 3336
Tel: (+86) 21 6141 5566
Fax: (+86) 21 6141 5567

Czech Republic

Tel: (+420) 246 003 200
Fax: (+420) 246 003 291

Denmark

Tel: (+45) 43 56 59 00
Fax: (+45) 43 56 59 05

Finland

Tel: (+358) 9 350 9250
Fax: (+358) 9 350 92555

France

Free Tel: 0800 211 408
Free Fax: 0800 031 052
Tel: (+33) 474 82 28 88
Fax: (+33) 474 95 68 08

Germany

Free Tel: 0800 51 55 000
Free Fax: 0800 64 90 000
Tel: (+49) 89 6513 0
Fax: (+49) 89 6513 1169

Hungary

Tel: (+36) 1 235 9055
Fax: (+36) 1 235 9068

India

Telephone

Bangalore: (+91) 80 6621 9400
New Delhi: (+91) 11 4358 8000
Mumbai: (+91) 22 4087 2364
Pune: (+91) 20 4146 4700
Hyderabad: (+91) 40 3067 7450
Kolkata: (+91) 33 4013 8000

Fax

Bangalore: (+91) 80 6621 9550
New Delhi: (+91) 11 4358 8001
Mumbai: (+91) 22 2579 7589
Pune: (+91) 20 4146 4777
Hyderabad: (+91) 40 3067 7451
Kolkata: (+91) 33 4013 8016

Ireland

Free Tel: 1800 200 888
Free Fax: 1800 600 222
Tel: +353 (0) 402 20370
Fax: +353 (0) 402 20375

Israel

Free Tel: 1 800 70 2222
Tel: (+972) 8 948 4222
Fax: (+972) 8 948 4200

Italy

Free Tel: 800 827 018
Tel: (+39) 02 3341 7310
Fax: (+39) 02 3801 0737

Japan

Tel: (+81) 3 5796 7300
Fax: (+81) 3 5796 7315

Korea

Free Tel: (+82) 80 023 7111
Free Fax: (+82) 80 023 8111
Tel: (+82) 31 329 9000
Fax: (+82) 31 329 9090

Luxembourg

Tel: (+32) 3 899 1301
Fax: (+32) 3 899 1311

Malaysia

Tel: (+60) 3 5635 3321
Fax: (+60) 3 5635 4116

Mexico

Free Tel: 01 800 007 5300
Free Fax: 01 800 712 9920
Tel: (+52) 722 276 1600
Fax: (+52) 722 276 1601

The Netherlands

Tel: (+31) 78 620 5411
Fax: (+31) 78 620 5421

New Zealand

Free Tel: 0800 936 666
Free Fax: 0800 937 777
Tel: (+61) 2 9841 0555
Fax: (+61) 2 9841 0500

Norway

Tel: (+47) 23 17 60 00
Fax: (+47) 23 17 60 10

Poland

Tel: (+48) 61 829 01 00
Fax: (+48) 61 829 01 20

Portugal

Free Tel: 800 202 180
Free Fax: 800 202 178
Tel: (+351) 21 924 2555
Fax: (+351) 21 924 2610

Russia

Tel: (+7) 495 621 5828
Fax: (+7) 495 621 6037

Singapore

Tel: (+65) 6779 1200
Fax: (+65) 6779 1822

Slovakia

Tel: (+421) 255 571 562
Fax: (+421) 255 571 564

South Africa

Free Tel: 0800 1100 75
Free Fax: 0800 1100 79
Tel: (+27) 11 979 1188
Fax: (+27) 11 979 1119

Spain

Free Tel: 900 101 376
Free Fax: 900 102 028
Tel: (+34) 91 661 99 77
Fax: (+34) 91 661 96 42

Sweden

Tel: (+46) 8 742 4200
Fax: (+46) 8 742 4243

Switzerland

Free Tel: 0800 80 00 80
Free Fax: 0800 80 00 81
Tel: (+41) 81 755 2511
Fax: (+41) 81 756 5449

Thailand

Tel: (+66) 2 126 8141
Fax: (+66) 2 126 8080

United Kingdom

Free Tel: 0800 717 181
Free Fax: 0800 378 785
Tel: (+44) 1747 833 000
Fax: (+44) 1747 833 313

United States

Toll-Free: 800 325 3010
Toll-Free Fax: 800 325 5052
Tel: (+1) 314 771 5765
Fax: (+1) 314 771 5757

Vietnam

Tel: (+84) 8 3516 2810
Fax: (+84) 8 6258 4238

Internet

sigma-aldrich.com



Enabling Science to
Improve the Quality of Life

Order/Customer Service (800) 325-3010 • Fax (800) 325-5052
Technical Service (800) 325-5832 • sigma-aldrich.com/techservice
Development/Custom Manufacturing Inquiries **SAFC**® (800) 244-1173
Safety-related Information sigma-aldrich.com/safetycenter

World Headquarters
3050 Spruce St.
St. Louis, MO 63103
(314) 771-5765
sigma-aldrich.com

©2013 Sigma-Aldrich Co. LLC. All rights reserved. SAFC, SIGMA-ALDRICH and SUPELCO are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. Ascentis is a registered trademark of Sigma-Aldrich Co. LLC. Solutions within and Supel are trademarks of Sigma-Aldrich Co. LLC. Supelco brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at www.sigmaaldrich.com and/or on the reverse side of the invoice or packing slip.

PFK
11644 / T413026
1023