

Product Information

Angiotensin II Receptor Type I, Human

Product Number **A 8602**

Product Description

Angiotensin II Receptor Type I, human is a frozen aliquot of membranes from KAN-TS cells.

Reagent

Angiotensin II Receptor Type I, human membranes are suspended in 50 mM TRIS-HCl containing 10% glycerol, 100 mM NaCl, 1 mM MgCl₂, 0.1% BSA, 0.1 mM bacitracin (pH 7.2 at 25°C).

Procedure

Thaw vials quickly; dilute with incubation buffer; homogenize; store on ice until addition to assay tubes

Incubation buffer:

50 mM TRIS-HCl, containing 100 mM NaCl, 1 mM MgCl₂, 0.1 mM Bacitracin, 0.1 % BSA (pH 7.2 at 25°C).

Binding Protocol

Membranes:

Dilute in incubation buffer (0.2 mL of membrane + 20 mL incubation buffer).

Assay mixture:

200 µl diluted membranes
 25 µl radioligand
 25 µl buffer or unlabeled ligand

Radioligand:

Sar¹, Tyr⁴-[¹²⁵I], Ile⁸-Angiotensin II at a final concentration of 0.06 nM for competition studies.

Unlabeled ligand:

Angiotensin II at a concentration of 1.0 µM

Incubation time:


180 minutes at 25 °C

Separation:

Over GF/B filter (presoaked in 0.1 % polyethylamine prepared in incubation buffer) then washed 5x with 1ml per tube of ice cold 50 mM NaCl, pH 7.4 at 4 °C.

Results

Typical results using standard binding assay above. Results may vary from lot to lot.



Storage/Stability

Store tightly sealed at -70 °C. The receptor can retain its original specific activity for several months when stored at -70 °C in its original packing solution. Repeated freeze-thaw of this product is not recommended.

Precautions

While no human toxicity data is available for this substance, it should be handled with care. Precautions should be taken to avoid contact by all routes of exposure.

References

- Bennett, J.P. and Snyder, S.H., *Jrnl. Biol. Chem.* **251**, 7423-7430 (1976).
- Wiest, S.A., et al., *Jrnl. Cardio. Pharmacol.*, **17**, 177-184 (1991).

LRP/JK 11/2003

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.