

- ▶ Solubility results in <4 hours
- ▶ Automation compatible
- ▶ High drug recovery
- ▶ QC released for filtration flow time, low extractables and extended incubations

MultiScreen[®] Solubility Filter Plate

96-well filter plate designed and optimized for aqueous solubility assays

Determining compound solubility in water has become an essential early measurement in the drug discovery process. Poor water solubility can cause problems in many different *in vitro* testing techniques, leading to unreliable results and/or reproducibility problems. Insoluble precipitates have been shown to cause false positives in bioassays potentially wasting valuable resources. Water solubility also influences absorption and thus can be used to help predict the ADME properties of a molecule.

Increased Screening Throughput and Efficiency

The MultiScreen Solubility filter plate is a high throughput, 96-well system used to classify or quantify the aqueous solubility of compound libraries stored in solvent as a concentrated solution.

The filtration-based protocol is fast and efficient. The protocol has been validated to screen hundreds of samples per day and saves time and sample in comparison with the shake-flask method. Typical MultiScreen Solubility filter plate analysis time is less than 4 hours and uses <100 µg of sample per analysis. This is an improvement as compared to the shake-flask method that can take several days and require milligrams of sample for each analysis.

MultiScreen Solubility filter plates are automation compatible. The plates are in compliance with SBS guidelines for easy handling by robotics. They also include a space for barcoding.

Reliable, Reproducible Results

MultiScreen Solubility filter plates are specifically designed and optimized for quantitative solubility assays. The plates demonstrate reproducible results in well-to-well, plate-to-plate and lot-to-lot comparisons.

Additionally, MultiScreen Solubility filter plate results show good agreement with shake-flask solubility values. Both methods determine solubility by a direct measure of compound in solution. This correlation with the solubility assay standard gives MultiScreen Solubility filter plate users an advantage over other methods, including nephelometry which measures insoluble precipitates rather than compound in solution.


Validated for High Drug Recovery

The MultiScreen Solubility filter plate incorporates low-binding membrane and low-binding plate materials to yield the high drug recovery needed for solubility results. The plate is validated versus a panel of 9 drugs for >80% drug recovery at 10 μM concentration in 5% DMSO/PBS. Example results at 1 μM drug are shown in Figure 2.

Filtration-based Solubility Results in Less than 4 Hours

The MultiScreen Solubility filter plate is tested to ensure reliable, discrete filtrate transfer with no cross-talk. The plate can be incubated for 6 hours or longer with shaking and have no drip out. As soon as vacuum is applied, the wells empty in less than 1 minute.

Protocol



1. Add compound dissolved in organic solvent to aqueous buffer.
2. Shake for 90 minutes to allow insoluble compound to precipitate.
3. Apply vacuum to filter solution into collection plate. Precipitates remain on membrane. Analyze filtrate in collection plate to quantitate amount of compound still in solution.

Performance

Results Correlate to Shake-flask




Figure 1. Fifteen (15) compounds were evaluated by an ASTM Shake-flask method and a protocol using the MultiScreen Solubility filter plate. Shake-flask solubility was determined for solid compounds added to PBS. MultiScreen Solubility filter plate results were determined for 10 μL of 10 mM DMSO stocks added to 190 μL PBS for a final DMSO concentration of 5%. The maximum concentration on the Y axis is given as 500 μM due to the MultiScreen Solubility filter plate protocol that has an upper limit for the amount of compound added.

High Drug Recovery




Figure 2. Soluble drugs were dissolved in 5% DMSO/PBS at 1 μ M, incubated in the MultiScreen Solubility filter plate and filtered into a receiver plate. The results are reported as percent drug recovery as determined by radiometric analysis.

Reproducible Results




Figure 3. Four compounds were analyzed in 45 wells per plate, 2 plates each over a 4 day period. The average and standard deviation from a given plate on a specific day is plotted in μ M units.

Ordering Information

Each MultiScreen Solubility filter plate includes a lid and a 96-well filter plate.

Description	Qty/Pk	Catalogue No.
MultiScreen Solubility filter plate	10	MSS LB PC 10

Required Equipment

Millipore Vacuum Manifold	1	MAVM 096 OR
---------------------------	---	-------------

Related Literature

Application note AN1730EN00: Quantitative method to determine drug aqueous solubility: optimization and correlation to standard methods

Application note AN1731EN00: Performance and correlation of a 96-well high throughput screening method to determine aqueous drug solubility

Protocol note PC2445EN00: Determination of aqueous compound solubility using a 96-well filter plate to remove precipitated solids prior to UV/Vis spectroscopic analysis

To Place an Order or Receive Technical Assistance

For additional information call your nearest Millipore office:

In the U.S. and Canada, call toll-free **1-800-MILLIPORE (1-800-645-5476)**

In the U.S., Canada and Puerto Rico, fax orders to **1-800-MILLIFX (1-800-645-5439)**

Internet: www.millipore.com

Technical Service:

www.millipore.com/techservice

pure e-commerce

Now you can buy Millipore products online @

www.millipore.com/purecommerce

Millipore and MultiScreen are registered trademarks of Millipore Corporation.

Lit. No. PF1315EN00 Rev.- 8/03 03-233

Printed in U.S.A.

© 2003 Millipore Corporation, Billerica, MA. All rights reserved.

Millipore Worldwide

AUSTRALIA

Tel. 1 800 222 111
or (02) 9888 8999
Fax (02) 9878 0788

AUSTRIA

Tel. 0820 87 44 64
Fax 0820 87 44 44

BALTIC COUNTRIES

Tel. +358 203 05 645
Fax +358 203 05 644

BELGIUM

Tel. 070 225 645
Fax 070 225 644

BRAZIL

Tel. (011) 5548-7011
Fax (011) 5548-7923

CANADA

Tel. (800) 645-5476
Fax (800) 645-5439

CHINA, PEOPLE'S REPUBLIC OF

Beijing:
Tel. (86-10) 8519 1250
(86-10) 6518 1058
Fax (86-10) 8519 1255
Guangzhou:
Tel. (86-20) 8752 0187
(86-20) 8752 0173
Fax (86-20) 8752 0172
Hong Kong:
Tel. (852) 2803 9111
Fax (852) 2513 0313

Shanghai:
Tel. (86-21) 5306 9100
Fax (86-21) 5306 0838

CZECH REPUBLIC

Tel. 02 2051 3841
Fax 02 2051 4298

DENMARK

Tel. 7010 5645
Fax 7010 5644

EASTERN EUROPE, C.I.S., AFRICA, MIDDLE EAST AND GULF

Tel. +33 3 88 38 9536
Fax +33 3 88 38 9539

FINLAND

Tel. 0203 05 645
Fax 0203 05 644

FRANCE

Tel. 0825 045 645
Fax 0825 045 644

GERMANY

Tel. 01 805 045 645
Fax 01 805 045 644

HUNGARY

Tel. 01 209 04 33
Fax 01 209 02 95

INDIA

Tel. 080 839 46 57
Fax 080 839 63 45

IRELAND

Tel. +44 1923 816375
Fax +44 1923 818297

ITALY

Tel. 848 845 645
Fax 848 845 644

JAPAN

Tel. (03) 5442-9711
Fax (03) 5442-9736

KOREA

Tel. (822) 551-0230
Fax (822) 551-0228

LUXEMBOURG

Tel. +32 2726 88 40
Fax +32 2726 98 84

MALAYSIA

Tel. 03-7957-1322
Fax 03-7957-1711

MEXICO

Tel. (55) 5576 9688
Fax (55) 5576 8706

THE NETHERLANDS

Tel. 0900 7 645 645
Fax 0900 7 645 644

NORWAY

Tel. 810 62 645
Fax 810 62 644

POLAND

Tel. 22 669 12 25
Fax 22 663 70 33

PORTUGAL

Tel. +34 91 728 39 60
Fax +34 91 729 29 09

PUERTO RICO

Tel. (787) 273-8495
Fax (787) 747-6553

SINGAPORE

Tel. 6842-1822
Fax 6842-4988

SPAIN

Tel. 901 516 645
Fax 902 011 644

SWEDEN

Tel. 0771 200 645
Fax 0771 200 644

SWITZERLAND

Tel. 0848 645 645
Fax 0848 645 644

TAIWAN

Tel. 886-2-2792-9333
Fax 886-2-2792-6555

U.K.

Tel. 0870 900 46 45
Fax 0870 900 46 44

U.S.A.

Tel. (800) 645-5476
Fax (800) 645-5439

IN ALL OTHER COUNTRIES

Millipore International
Tel. +1 781 533 8622
Fax +1 781 533 8630

MILLIPORE