

# Development of Multiplex Immunoassay Panels for Simultaneous Quantification of Bone Metabolism Markers using Luminex® xMAP® Technology



Wen-Rong Lie, Hank Hwang, Jin Wang, Rick Ryan, Jehangir Mistry Millipore Corporation, 15 Research Park Drive, St. Charles, MO 63304


upstate • CHEMICON • Linco  
The expertise of Upstate®, Chemicon® and Linco® is now a part of Millipore

## Abstract

Biochemical markers of bone metabolism play an important role in the assessment of bone diseases such as osteoporosis, arthritis, chronic inflammatory disorders, and bone metastasis with cancers. Millipore has developed various multiplexed immunoassay panels for the simultaneous measurement of multiple bone metabolism markers in mouse, rat, and human, using the Luminex® xMAP® platform. These bead-based sandwich assays are rapid, sensitive, and reproducible and require  $\leq 25 \mu\text{l}$  of serum/plasma sample or tissue culture supernatants per well in 96-well plates. The mouse bone panels allow simultaneous quantification of the following biomarkers in any combinations in the same sample well: Mouse Panel 1 - OPG, Osteocalcin, ACTH, Insulin, Leptin, TNF $\alpha$ , IL-6, IL-1 $\beta$  and Mouse Panel 2 - RANKL, Osteocalcin, ACTH, Insulin, Leptin, TNF $\alpha$ , IL-6, IL-1 $\beta$ . The dynamic ranges for the mouse analytes are 10-40,000 pg/mL (OPG, Osteocalcin, RANKL, Leptin), 24-100,000 pg/mL (Insulin), and 2-10,000 pg/mL (ACTH, TNF $\alpha$ , IL-6, IL-1 $\beta$ ). The rat bone panels allow simultaneous measurement of the following biomarkers: Rat Panel 1 - OPG, ACTH, Insulin, Leptin, Rat Panel 2 - RANKL, ACTH, Insulin, Leptin, and Rat Panel 3 - Osteocalcin, PTH, Osteopontin. The dynamic ranges for the rat analytes are 10-40,000 pg/mL (OPG, PTH, Leptin), 5-20,000 pg/mL (Osteocalcin), 4-15,000 pg/mL (RANKL), 24-100,000 pg/mL (insulin), and 2-10,000 pg/mL (ACTH). The human bone metabolism panels allow for simultaneous measurement of the following biomarkers in any combinations: Human Panel 1A serum samples - OPG, Osteocalcin, Osteopontin, PTH, Leptin, ACTH, Insulin and Human Panel 1B (tissue culture samples) - OPG, Osteocalcin, Osteopontin, PTH, Leptin, Adiponectin, Insulin, ACTH, TNF $\alpha$ , IL-6, IL-1 $\beta$ . We have also developed a sensitive immunoassay for measuring the serum levels of human RANKL. The dynamic range for this Human RANKL Single-plex assay is 5-20,000 pg/mL. The above bone metabolism panels exhibited acceptable analytical performance characteristics in terms of sensitivity, intra- and inter-assay precision, linearity of dilution, spike recovery, and antibody pair specificity. In conclusion, various multiplexed assay panels were developed and validated for the measurement of multiple bone metabolism markers in mouse, rat, and human samples. The availability of these multiplexed Bone Metabolism Panels and the RANKL, OPG, and Osteocalcin Single-plex assays may provide a powerful tool in studying biological functions of these biomarkers as well as the pathological roles of these molecules.

### Luminex® xMAP® Technology

## Bone Panel Assay Method



**Beads:** Capture antibodies were covalently coupled to the carboxylate-modified microsphere beads using EDC and Sulfo-NHS according to the manufacturer's instructions.

**Assay Procedures:** The multiplex assay was performed in a 96-well Millipore MultiScreen filter plate. The protocol is as follows: wet the plate with 200  $\mu\text{l}$  assay buffer  $\rightarrow$  25  $\mu\text{l}$  assay buffer, 25  $\mu\text{l}$  standards/samples, 25  $\mu\text{l}$  matrix/buffer, 25  $\mu\text{l}$  beads  $\rightarrow$  incubate at 4°C, overnight  $\rightarrow$  wash plate 3 times  $\rightarrow$  add 50  $\mu\text{l}$  biotinylated secondary Abs, RT, 1 hr  $\rightarrow$  add 50  $\mu\text{l}$  Streptavidin-Phycoerythrin (SAPE), RT, 30 min  $\rightarrow$  wash plate 3 times  $\rightarrow$  read on Luminex instruments

### Mouse Bone Metabolism Panels

7 Mouse Bone Panels: 4 Multiplex and 3 Single-plex

Mouse Panel 1A <sup>†</sup>	Mouse Panel 1B <sup>*</sup>	Mouse Panel 2A <sup>†</sup>	Mouse Panel 2B <sup>*</sup>	Mouse Single-plex
OPG Insulin Leptin ACTH IL-6 TNF $\alpha$ IL-1 $\beta$	OPG Insulin Leptin ACTH IL-6 TNF $\alpha$ OC IL-1 $\beta$	RANKL Insulin Leptin ACTH IL-6 TNF $\alpha$ E-1 $\beta$	RANKL Insulin Leptin ACTH IL-6 TNF $\alpha$ OC E-1 $\beta$	Mouse OC <sup>‡</sup> OPG <sup>†</sup> RANKL <sup>†</sup>

<sup>†</sup> Serum/plasma <sup>\*</sup> Tissue culture samples or  $>1:20$  diluted serum/plasma <sup>‡</sup> Special release

### Mouse Bone Panels: Assay Characteristics

**ASSAY CONDITIONS:** Overnight assay  
25  $\mu\text{l}$  or less sample volume

**STANDARD CURVE RANGE:**  
10 - 40,000 pg/mL - Leptin, OPG, Osteocalcin, RANKL  
24 - 100,000 pg/mL - Insulin  
2 - 10,000 pg/mL - ACTH, IL-6, TNF $\alpha$

Analyte	SENSITIVITY (Minimum Detectable Concentration, MinDC, pg/mL):				ACCURACY (%):			
	1A	1B	2A	2B	1A	1B	2A	2B
RANKL	—	—	3.3	2.7	—	—	119	78
OPG	2.3	2.4	—	—	—	95	90	—
Insulin	18.6	17.4	16.2	14.4	—	92	101	92
Leptin	3.0	3.2	4.5	4.4	—	87	93	93
ACTH	1.8	1.4	1.7	1.3	—	117	92	115
IL-6	0.6	0.6	0.6	0.6	—	86	87	89
TNF $\alpha$	1.0	0.8	0.8	0.6	—	122	92	133
Osteocalcin	—	8.5	—	4.7	—	—	99	97

**LINEARITY OF DILUTION:**  
1A 108-137%  
1B 87-123%  
2A 86-129%  
2B 98-125%

**PRECISION:**  
Intra-assay <4%  
Inter-assay <11%

**CROSS REACTIVITY:** The antibody pairs in each panel do not show cross reactivity to the other analytes within the panel.

### Rat Bone Metabolism Panels

6 Rat Bone Panels: 3 Multiplex and 3 Single-plex

Rat Panel 1 <sup>†</sup>	Rat Panel 2 <sup>†</sup>	Rat Panel 3 <sup>*†</sup>	Rat Single-plex
OPG Insulin Leptin ACTH	RANKL Insulin Leptin ACTH	OC OPN PTH	Rat OC <sup>*</sup> OPG <sup>†</sup> RANKL <sup>†</sup>

<sup>†</sup> Serum/plasma or tissue culture <sup>\*</sup> Tissue culture or  $>1:50$  diluted serum/plasma <sup>‡</sup> Special release

### Rat Bone Panels: Assay Characteristics

**ASSAY CONDITIONS:** Overnight assay  
25  $\mu\text{l}$  or less sample volume

**STANDARD CURVE RANGE:**  
2.4 - 10,000 pg/mL - ACTH  
3.7 - 15,000 pg/mL - RANKL  
9.8 - 40,000 pg/mL - OPG, Leptin  
24 - 100,000 pg/mL - Insulin

Analyte	SENSITIVITY (Minimum Detectable Concentration, MinDC, pg/mL):		ACCURACY (%):
	Panel 1	Panel 2	
OPG	2.3	—	—
Insulin	18.6	16.2	—
Leptin	3.0	4.5	—
ACTH	1.8	1.7	—
RANKL	—	1.1	—

Analyte	ACCURACY (%):	
	Panel 1	Panel 2
OPG	78.6	—
Insulin	75.8	76.1
Leptin	78.0	87.7
ACTH	99.2	101.8
RANKL	—	118.4

**LINEARITY OF DILUTION:**  
Panel 1 82-117%  
Panel 2 94-120%

**PRECISION:**  
Intra-assay <4%  
Inter-assay <9%

**CROSS REACTIVITY:** The antibody pairs in each panel do not show cross reactivity to the other analytes within the panel.

### Human Bone Metabolism Panels

3 Human Bone Panels: 2 Multiplex and 1 Single-plex

Human Panel 1A <sup>†</sup>	Human Panel 1B <sup>*</sup>	Human Single-plex
OPG OC OPN PTH Leptin ACTH Insulin	OPG OC OPN PTH E-1 $\beta$ TNF $\alpha$ IL-6 Leptin ACTH Adiponectin Insulin	RANKL <sup>‡</sup>

<sup>†</sup> 1:4 diluted serum/plasma <sup>\*</sup> Tissue culture samples <sup>‡</sup> 2 diluted serum/plasma

### Human Bone Metabolism Panel 1A (in Serum Matrix)

Analyte	Sensitivity (pg/ml)	Accuracy (% Recovery)	Linearity (% Recovery)	Intra-Assay (% CV)	Inter-Assay (% CV)
OPG	1.42	96.4	105.0	6.5	6.3
OC	39.3	98.4	85.6	5.6	10.7
Leptin	119.7	136.3	138.0	4.8	6.2
OPN	67.8	79.7	126.5	5.4	5.8
PTH	0.30	98.1	117.5	7.0	9.6
ACTH	2.1	89.2	120.3	4.9	5.6
Insulin	46.3	95.2	104.1	3.5	6.3

**Sensitivity:** Minimum Detectable Concentration (MinDC)  
**Accuracy:** Percentage of measured analytes in 4 independent samples spiked with 3 known concentration of analytes  
**Dilution Linearity:** Average percent recovery of 4 independent samples diluted 1:2, 1:4, and 1:8 in the kit provided matrix  
**Precision:** Intra-assay precision is generated from the mean of the %CV from 10 reportable results across 2 different concentrations of analytes in a single assay. Inter-assay precision is generated from the mean of the %CV from 2 reportable results across 2 different concentrations of analytes across 6 different assays.

**Cross Reactivity:** The antibody pair in each panel do not show cross reactivity to the other analytes within the panel.

### Human Bone Metabolism Panel 1B (in Assay Buffer)

Analyte	Sensitivity (pg/ml)	Accuracy (% Recovery)	Linearity (% Recovery)	Intra-Assay (% CV)	Inter-Assay (% CV)
IL-1 $\beta$	0.13	93.3	109.4	3.5	7.3
IL-6	0.09	94.5	94.5	3.9	6.1
OPG	0.90	102.5	104.0	3.2	4.5
OC	39.8	101.2	91.4	3.7	6.7
Leptin	18.8	94.0	118.4	3.7	4.7
OPN	32.8	103.7	100.9	4.0	5.6
PTH	0.32	98.1	104.3	5.2	9.8
TNF $\alpha$	0.06	101.2	101.7	4.8	5.9
ACTH	1.3	96.6	96.2	3.4	7.8
Adiponectin	22.2	96.9	97.4	3.3	10.2
Insulin	38.5	100.0	96.8	2.7	5.5

### Human RANKL Single-plex (in Serum Matrix)

Analyte	Sensitivity (pg/ml)	Accuracy (% Recovery)	Linearity (% Recovery)	Intra-Assay (% CV)	Inter-Assay (% CV)
RANKL	4.8	93.2	121.7	<5.0	<6.0

### Bone Panel Applications

#### Bone Biomarkers in Drug Development

Drug Screening & Lead Identification	Clinical I, II, III, IV
In vitro Cell-based Assay	In vivo Animal Disease Models & PK/PD
Osteoblasts Osteoclasts Bone Culture Synovial Fibroblasts Stem Cell Culture	Inflammatory Arthritis Models LPS Shock Models Osteoporosis Models Transgenics & Knockouts* Cancer Models Bone Safety
	Clinical Trials
	Rheumatoid Arthritis Osteoarthritis Osteoporosis Bone Loss & Bone Healing Paget's Disease Cancer & Bone Metastases

\*M. Shih et al. SA184, ASBMR (2006)

### Summary

- We have developed mouse, rat, and human specific multiplex immunoassay panels for simultaneous measurement of bone metabolism markers. Analytic parameters evaluated show the robustness of the assay.
- The assay is rapid, sensitive, and reproducible and only 25  $\mu\text{l}$  of sample is required.
- The availability of Millipore multiplex bone metabolism panels provides a useful tool in studying biological functions of these bone biomarkers and the pathological roles of these molecules.

The M logo is a trademark of Millipore Corporation.  
Upstate, Chemicon, Linco and LINCOPlex are registered trademarks of Millipore Corporation.  
Luminex and xMAP are registered trademarks of Luminex Corporation.  
©2007 Millipore Corporation, St. Charles, MO 63304 U.S.A. 10/07 All rights reserved.