



## Product Information

### RNA Marker template set

Catalog Number **R4142**  
Storage Temperature  $-20\text{ }^{\circ}\text{C}$

#### Product Description

The RNA Marker template set provides a mixture of 7 linearized DNA templates. Each template contains the promoter for T7 RNA polymerase and upon transcription with T7 RNA polymerase results in 7 transcripts of the following lengths: 100, 200, 300, 400, 600, 800, and 1,000 bases. The concentration of each template has been adjusted so all 7 bands are approximately equal in intensity. The transcripts can be labeled with radioisotopes, biotin, or any non-radioactive tag compatible with T7 RNA polymerase. The markers are useful for RNase protection assays or as low molecular mass RNA markers for electrophoresis.

The product is supplied as a solution in 10 mM Tris, pH 7.5, with 1 mM EDTA.

#### Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

#### Storage/Stability

Store the product at  $-20\text{ }^{\circ}\text{C}$ .

The RNA transcripts can be stored at  $-20\text{ }^{\circ}\text{C}$  for several weeks and reheating the markers prior to use is not necessary.

#### Procedure

##### Non-labeled Markers

For non-labeled markers, add 0.5  $\mu\text{g}$  of the RNA Marker template set to a 20  $\mu\text{l}$  reaction containing 0.5 mM of each of the ribonucleotide triphosphates with 10 units of T7 RNA polymerase (Catalog Number R0884).

##### Labeled Markers

1. Prepare a 20  $\mu\text{l}$  reaction mix as follows:
  - 2  $\mu\text{l}$  10 $\times$  Transcription Buffer (400 mM Tris-HCl, pH 8.0, with 80 mM  $\text{MgCl}_2$ , 500 mM NaCl, and 20 mM spermidine)
  - 1  $\mu\text{l}$  200 mM dithiothreitol
  - 1  $\mu\text{l}$  10 mM ATP
  - 1  $\mu\text{l}$  10 mM GTP
  - 1  $\mu\text{l}$  10 mM UTP
  - 1  $\mu\text{l}$  10 mM CTP
  - 1  $\mu\text{l}$  RNA Marker template set (diluted if necessary to 0.5  $\mu\text{g}/\mu\text{l}$ )
  - 1-3  $\mu\text{l}$   $\alpha\text{-}^{32}\text{P}$ -UTP or  $\alpha\text{-}^{32}\text{P}$ -CTP, 800 Ci/mmol (10 mCi/ml in aqueous solution)
  - 1  $\mu\text{l}$  T7 RNA polymerase (10 units/ $\mu\text{l}$ )

Bring final volume to 20  $\mu\text{l}$  with RNase-free water (Catalog Number W4502)

Note: Use of an RNase inhibitor may be helpful in transcription reaction.

2. Incubate at  $37\text{ }^{\circ}\text{C}$  for 1 hour.
3. Add 1  $\mu\text{l}$  (2 units/ $\mu\text{l}$ ) of RNase-free DNase I (Catalog Number D7291) to degrade the DNA template, mix well, and incubate at  $37\text{ }^{\circ}\text{C}$  for 15 minutes.
4. Add an equal volume of gel loading buffer (80% formamide, 0.1% xylene cyanole, 0.1% bromophenol blue, and 2 mM EDTA).
5. Heat for 3 minutes at  $95\text{ }^{\circ}\text{C}$  to inactivate the enzyme and denature the transcripts.
6. Separate transcripts by electrophoresis on a 5% polyacrylamide/8 M urea gel.

Notes:

The 100 nucleotide band runs between the xylene cyanole and the bromophenol blue. The remaining bands migrate slower than the xylene cyanole.

Approximate exposure times for radiolabeled markers:

10 minutes for 10-20  $\mu$ l load

1 hour for 1-3  $\mu$ l load

12-16 hours for 1-3  $\mu$ l load of 1:10 dilution

Using an intensifying screen:

3 hours for 2-5  $\mu$ l load of 1:10 dilution

12-16 hours for 1-5  $\mu$ l load of 1:50 dilution

Increase volume of marker loaded on the gel proportionally as the  $^{32}$ P-label decays.

**References**

1. Sambrook, J., *et al.*, Molecular Cloning: A Laboratory Manual, Cold Spring Harbor Laboratory Press (Cold Spring Harbor, NY: 1989) p. 10.27-10.37.
2. Krieg, P.A., and Melton, D.A., Nucleic Acids Res., **12**, 7057-7070 (1984).
3. Melton, D.A., Proc. Natl. Acad. Sci. USA, **82**, 144-148.

JB,KH,MAM 01/07-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.