



# **Alkaline Phosphatase Detection Kit**

For 100 Tests

**Catalog No. SCR004**

**FOR RESEARCH USE ONLY  
Not for use in diagnostic procedures**

USA & Canada

Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233  
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered Worldwide  
[www.Millipore.com](http://www.Millipore.com) • [custserv@Millipore.com](mailto:custserv@Millipore.com) • [techserv@Millipore.com](mailto:techserv@Millipore.com)

This page left blank intentionally.

---

## Application

Embryonic stem (ES) cells are totipotent cells derived from the inner cell mass (ICM) of preimplantation mammalian embryos and are capable of unlimited, undifferentiated proliferation *in vitro* [1,2,3]. Undifferentiated murine ES cells can be maintained *in vitro* for extended periods in media containing the cytokine, leukemia-inhibitory factor (LIF) or MILLIPORE®'s proprietary ES cell culture reagent, ESGRO® [4,5]. The undifferentiated state of ES cells can be characterized by high level of expression of Alkaline Phosphatase (AP) [6], the expression of surface markers including SSEA and TRA antigens and the transcription factor Oct-4.

MILLIPORE®'s Alkaline Phosphatase Detection Kit (Catalog number SCR004) is a specific and sensitive tool for the phenotypic assessment of ES cell differentiation by the determination of AP activity.

Available separately from MILLIPORE® are the monoclonal antibodies TRA-2-49 (Catalog number MAB4349) and TRA-2-54 (Catalog number MAB4354), which permit the detection of Liver / Bone / Kidney isozyme of AP as expressed by ES cells. When used in flow cytometry, these reagents permit a quantitative assessment of AP expression by ES cells [7].

*For Research Use Only; Not for use in diagnostic procedures*

---

## Related Products

The following related products are available from MILLIPORE® as separate items:

1. ES Cell Characterization Kit (Catalog # SCR001)
2. ES Marker Sample Kit (Catalog # SCR002)
3. ES Cell 3D Collagen Culture Kit (Catalog # SCR003)
4. SSEA-1 Monoclonal Antibody, purified 100µg (Catalog # MAB4301)
5. SSEA-3 Monoclonal Antibody, purified 100µg (Catalog # MAB4303)
6. SSEA-4 Monoclonal Antibody, purified 100µg (Catalog # MAB4304)
7. TRA-1-60 Monoclonal Antibody, purified 100µg (Catalog# MAB4360)
8. TRA-1-81 Monoclonal Antibody, purified 100µg (Catalog# MAB4381)
9. TRA-2-49 Monoclonal Antibody, purified 100µg (Catalog # MAB4349)
10. TRA-2-54 Monoclonal Antibody, purified 100µg (Catalog # MAB4354)
11. Murine LIF, 5µg (Catalog # LIF2005), 10µg (Catalog # LIF2010)
12. ESGRO®, 1 x 10<sup>6</sup> units (Catalog # ESG1106), 1 x 10<sup>7</sup> units (Catalog # ESG1107)

---

## Storage

The Alkaline Phosphatase Detection Kit consists of two components used for determining AP activity. When stored at 2° to 8°C, the kit components are stable up to the expiration date. Do not freeze or expose to elevated temperatures. Discard any remaining reagents after the expiration date.

---

## Kit Components

1. Fast Red Violet solution (0.8g/L stock) (Part No. 90239). Two 15mL bottles.
2. Naphthol AS-BI phosphate solution (4mg/mL) in AMPD buffer (2mol/L), pH 9.5 (Part No. 90234). One 15mL bottle.

---

## Materials Required But Not Supplied

1. Fixative (e.g. 4% Paraformaldehyde)
2. 1 x Rinse Buffer (e.g. TBST: 20mM Tris-HCl, pH 7.4, 0.15M NaCl, 0.05% Tween-20)
3. Hematoxylin (optional)
4. Microscope

---

## Preparation of Reagents

1. **Naphthol/Fast Red Violet Solution:** Mix Fast Red Violet (FRV) with Naphthol AS-BI phosphate solution and water in a 2:1:1 ratio (FRV:Naphthol:water).

---

## Staining Protocol

### Alkaline Phosphatase Staining Procedure

1. Culture ES cells for five days prior to analyzing AP activity, at low to medium density (*NOTE: This time-period is critical if activity levels of AP needs to be observed. According to our findings, five days of culturing are optimal for good AP stain visualization*).
2. On day five, aspirate media and fix ES or EC cells with a fixative (e.g. 4% Paraformaldehyde in PBS) for 1-2 minutes.


*Note: Do not overfix. Fixing cells longer than 2 minutes will result in the inactivation of alkaline phosphatase.*

3. Aspirate fixative and rinse with 1 X Rinse Buffer. DO NOT allow wells to dry.
4. Prepare reagents for Alkaline Phosphatase staining as described in “Preparation of Reagents” section.

5. Add enough stain solution to cover each well (e.g. 0.5mL for a well of a 24-well plate). Incubate in dark at room temperature for 15 minutes.
6. Aspirate staining solution and rinse wells with 1 X Rinse Buffer. Cover cells with 1 X PBS to prevent drying and then count the number of colonies expressing AP (red stem cell colonies), versus the number of differentiated colonies (colorless).
7. AP staining criteria: Greater than 90% of colonies should remain undifferentiated and express alkaline phosphatase in the well containing  $10^3$  Units of LIF. P value shall be  $\geq 0.05$ .

---

## Staining with Alkaline Phosphatase Detection Kit



**Alkaline Phosphatase staining of ES cells.** High magnification revealed **(A)** Undifferentiated ES cells (mouse MBL.5 cell line) – cultured for five days in media containing MILLIPORE®'s LIF/ESGRO®. A concentration of  $10^3$  Units/mL is used for inhibition of differentiation. **(B)** Differentiated ES cells – cultured at low-medium density for three days in media without any LIF/ESGRO®. **(C)** Differentiated ES cells – cultured at low-medium density for six days in media without any LIF/ESGRO®.

---

## References

1. Evans M, Kaufman M. *Nature* **292**, 154 (1981).
2. Martin G. *Proc Natl Acad Sci USA* **78**, 7634 (1981).
3. Morrison SJ, Shah NM, Anderson DJ. *Cell* **88**, 287-298 (1997).
4. Smith AG, Heath JK et al. *Nature* **336**, 688-690 (1988).
5. Williams RL, Hilton DJ et al. *Nature* **336**, 684-687 (1988).
6. Pease S, Braghetta P et al. *Develop Biol* **141**, 344 (1990).
7. Draper J, et al. *J.Anat.* **200**, 249 - 258 (2002)

---

## Warranty

**EMD Millipore Corporation** (“EMD Millipore”) warrants its products will meet their applicable published specifications when used in accordance with their applicable instructions for a period of one year from shipment of the products. **EMD MILLIPORE MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** The warranty provided herein and the data, specifications and descriptions of EMD Millipore products appearing in EMD Millipore’s published catalogues and product literature may not be altered except by express written agreement signed by an officer of EMD Millipore. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorized and if given, should not be relied upon.

In the event of a breach of the foregoing warranty, EMD Millipore Corporation’s sole obligation shall be to repair or replace, at its option, the applicable product or part thereof, provided the customer notifies EMD Millipore Corporation promptly of any such breach. If after exercising reasonable efforts, EMD Millipore Corporation is unable to repair or replace the product or part, then EDM Millipore shall refund to the Company all monies paid for such applicable Product. **EMD MILLIPORE CORPORATION SHALL NOT BE LIABLE FOR CONSEQUENTIAL, INCIDENTAL, SPECIAL OR ANY OTHER DAMAGES RESULTING FROM ECONOMIC LOSS OR PROPERTY DAMAGE SUSTAINED BY ANY COMPANY CUSTOMER FROM THE USE OF ITS PRODUCTS.**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

(c) 2009 - 2016: Merck KGaA, Darmstadt. All rights reserved. No part of these works may be reproduced in any form without permission in writing

Cat No. SCR004

August 2016  
Version 5: 41447



