

# sealed with an



The new solvent bottle  
with advanced sealing  
technology

Meet the most innovative packaging system for solvents: **the new M-Bottle**. Our ingenious bottle and cap introduces an enhanced design featuring advanced sealing technology that eliminates leaks, assuring you that your solvent handling is safer and easier than ever before.

## SUPERIOR SAFETY


- Advanced sealing technology effectively eliminates leaks
- Additional headspace ensures greater safety with solvents that have higher vapor pressure at ambient temperatures
- Thicker glass reduces risk of breakage
- Also coming soon: PE-coated M-Bottle with added environmental benefits

## EASY HANDLING

- New cap with larger ribs allows easier removal and closure
- Redesigned pouring ring facilitates handling and dispensing

## GREATER CONVENIENCE

- Markings on bottle shoulder indicate a nominal volume of 4L
- Designed to fit in standard storage spaces and bottle carriers



MilliporeSigma is the U.S. and Canada Life Science business of Merck KGaA, Darmstadt, Germany.

**Sigma-Aldrich**<sup>®</sup>

Lab & Production Materials

**Supelco**<sup>®</sup>

Analytical Products

# The new M-Bottle for solvents



**New neck design**  
captures tamper-  
evident ring

**Innovative sealing technology**  
effectively prevents leaks

**Additional finger  
clearance**  
easier, safer pouring

**Additional headspace**  
greater safety for  
solvents with higher  
vapor pressure at  
ambient temperatures

**Optimized shape**  
fits in all standard  
bottle carriers and  
storage spaces



## To place an order or receive technical assistance

Order/Customer Service: [SigmaAldrich.com/order](https://SigmaAldrich.com/order)  
Technical Service: [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice)  
Safety-related Information: [SigmaAldrich.com/safetycenter](https://SigmaAldrich.com/safetycenter)

### [SigmaAldrich.com/m-bottle](https://SigmaAldrich.com/m-bottle)

We have built a unique collection of life science brands  
with unrivalled experience in supporting you scientific advancements.

**Millipore® Sigma-Aldrich® Supelco® Milli-Q® SAFC® BioReliance®**

MilliporeSigma  
400 Summit Drive  
Burlington, MA 01803

