



QUALITY CONTROL RANGES
MILLIPLEX[®] MAP Human MMP Magnetic Bead Panel 1
Kit Catalog # HMMP1MAG-55K

Control Catalog # HMMP1MG-6055-1
Lot# HMMPM-107 and HMMPM-207

<u>Hormone</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
HMMP-3	Control I	757 - 1572	pg/mL
	Control II	7962 - 16537	pg/mL
HMMP-12	Control I	631 - 1311	pg/mL
	Control II	5280 - 10967	pg/mL
HMMP-13	Control I	358 - 743	pg/mL
	Control II	3262 - 6776	pg/mL

Note: The Quality Control Ranges are generated with same day assay format using serum matrix provided in the kit. Quality Control values in culture media are not tested.