



**QUALITY CONTROL RANGES**  
**MILLIPLEX<sup>®</sup> MAP Human Apolipoprotein Magnetic Bead Panel Kit**  
**Catalog # AP0MAG-62K**

**Standard Catalog # HAP0-8062**  
**Lot# 862M-9K**

<b>Apolipoproteins</b>	<b>Standard Value</b>		<b>Apolipoproteins</b>	<b>Standard Value</b>
Apo AI	51460 ng/mL		Apo AII	4344 ng/mL
Apo CIII	3466 ng/mL		Apo B	24920 ng/mL
Apo E	1848 ng/mL		Apo CII	1755 ng/mL

**Control Catalog # HAP0-6062**  
**Lot# AP0M-109 and AP0M-209**

<u>Apolipoproteins</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
Apo-AI	Control I	AP0M-109	225 - 467	ng/mL
	Control II	AP0M-209	1108 - 2302	ng/mL
Apo-CIII	Control I	AP0M-109	15 - 31	ng/mL
	Control II	AP0M-209	71 - 148	ng/mL
Apo-E	Control I	AP0M-109	5.3 - 11	ng/mL
	Control II	AP0M-209	31 - 64	ng/mL
Apo-AII	Control I	AP0M-109	10 - 21	ng/mL
	Control II	AP0M-209	64 - 134	ng/mL
Apo-B	Control I	AP0M-109	69 - 143	ng/mL
	Control II	AP0M-209	371 - 771	ng/mL
Apo-CII	Control I	AP0M-109	5.8 - 12	ng/mL
	Control II	AP0M-209	32 - 67	ng/mL

HAP0-6062: AP0M-109.209 (23 June 2020)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • [emdmillipore.com](http://emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.