


## Product Information



### 60403 Potassium ionophore I

(Valinomycin)

Selectophore®, function tested

## Electrochemical Transduction Ion-Selective Electrodes

### Application 1 and Sensor Type<sup>1-5</sup>

Assay of K<sup>+</sup> activity in diluted urine, whole blood, plasma, serum and aqueous solutions with solvent polymeric membrane electrodes based on Potassium ionophore I.

### Recommended Membrane Composition

- 1.00 wt% Potassium ionophore I ([60403](#))
- 65.50 wt% Bis(1-butylpentyl) decane-1,10-diyl diglutarate (ETH 469) ([30585](#))<sup>a)</sup>
- 0.50 wt% Potassium tetrakis(4-chlorophenyl)borate ([60591](#))
- 33.00 wt% Poly(vinyl chloride) high molecular weight ([81392](#))

<sup>a)</sup> The use of bis(1-butylpentyl) adipate (BBPA) or bis(2-ethylhexyl)sebacate (DOS) leads to membrane electrodes of similar performance.

### Recommended Cell Assembly

Reference || sample solution || liquid membrane | 0.004 M KCl | AgCl, Ag

### Electrode Characteristics and Function

Selectivity coefficients  $\log K_{K,M}^{Pot}$ .

	Required <sup>b)</sup>	Found <sup>c)</sup>
$\log K_{K,H}^{Pot}$	<-2.8	-3.4
$\log K_{K,Na}^{Pot}$	<-3.6	-4.1
$\log K_{K,Mg}^{Pot}$	<-2.8	-5.7
$\log K_{K,Ca}^{Pot}$	<-2.9	-5.2
Stability: Drift [mV h <sup>-1</sup> ]		0.01
Standard deviation [mV]	<-0.46	0.03
Reproducibility [mV]		0.16
Lifetime: $\log P_{TLC}^{d)}$ ionophore	>8.4	8.6
$\log P_{TLC}^{d)}$ plasticizer	>12.8	10.8

<sup>b)</sup> for measurements in blood (1% interference, worst case)<sup>6,7</sup>

<sup>c)</sup> membrane without potassium tetrakis(4-chlorophenyl)borate

<sup>d)</sup> lipophilicity, determined by thin-layer chromatography<sup>8</sup>



### Application 2 and Sensor Type<sup>4,9,10</sup>

Assay of K<sup>+</sup> activity in undiluted urine, whole blood, plasma, serum, and aqueous solutions with solvent polymeric membrane electrodes based on Potassium ionophore I.

#### Recommended Membrane Composition

2.50 wt%	Potassium ionophore I ( <a href="#">60403</a> )
83.00 wt%	Siloprene K 1000 ( <a href="#">85417</a> )
14.50 wt%	Siloprene Crosslinking Agent K11 ( <a href="#">85418</a> )

#### Recommended Cell Assembly

Reference || sample solution || liquid membrane | 0.1 M KCl | AgCl, Ag

#### Electrode Characteristics and Function

Selectivity coefficients  $\log K_{K,M}^{Pot}$  as obtained by the separate solution method.

	Required <sup>a)</sup>	Found
$\log K_{K,H}^{Pot}$	<-2.8	-4.4
$\log K_{K,Li}^{Pot}$	<-1.7 <sup>b)</sup>	-4.3
$\log K_{K,Na}^{Pot}$	<-3.6	-4.0
$\log K_{K,Mg}^{Pot}$	<-2.8	-4.3
$\log K_{K,Ca}^{Pot}$	<-2.9	-4.2
Lifetime: $\log P_{TLC}^{c)}$ ionophore	>8.4	8.6

<sup>a)</sup> for measurements in blood (1% interference, worst case)<sup>6,7</sup>

<sup>b)</sup> therapeutical Li<sup>+</sup> concentrations

<sup>c)</sup> lipophilicity, determined by thin-layer chromatography<sup>8</sup>

### Application 3 and Sensor Type<sup>11-14</sup>

Assay of K<sup>+</sup> activity for cardiovascular application with biocompatible solvent polymeric membrane electrodes and related microfabricated sensor arrays based on Potassium ionophore I.

#### Recommended Membrane Composition

1.00 wt%	Potassium ionophore I ( <a href="#">60403</a> )
0.50 wt%	Potassium tetrakis(4-chlorophenyl)borate ( <a href="#">60591</a> )
49.50 wt%	Bis(2-ethylhexyl)sebacate ( <a href="#">84818</a> )
49.00 wt%	Poly(vinyl chloride) carboxylated (1.8% carboxyl content) ( <a href="#">81395</a> )

#### Electrode Characteristics and Function

Selectivity coefficients  $\log K_{K,M}^{Pot}$  as obtained by the separate solution method (0.1 M solutions of the chlorides).

$\log K_{K,Li}^{Pot}$	-3.98	$\log K_{K,Ca}^{Pot}$	-4.30
$\log K_{K,Na}^{Pot}$	-3.56	$\log K_{K,NH_4}^{Pot}$	-1.81

Slope of linear regression:

59-60 mV/dec

Detection limit (KCl ion background of 140 mM Na<sup>+</sup>):

4·10<sup>-6</sup> M K<sup>+</sup>

Membrane resistance (10<sup>-2</sup> M KCl at pH 7.0 with Tris buffer):

3.17 MΩ

### Application 4 and Sensor Type<sup>15,16</sup>

Assay of K<sup>+</sup> activity with double-matrix membrane ion-selective electrodes based on Potassium ionophore I

#### Recommended Membrane Composition

2.00 wt%	Potassium ionophore I ( <a href="#">60403</a> )
0.50 wt%	Potassium tetrakis(4-chlorophenyl)borate ( <a href="#">60591</a> )
64.70 wt%	Bis(2-ethylhexyl)sebacate ( <a href="#">84818</a> )
32.80 wt%	Poly(vinyl chloride) high molecular weight ( <a href="#">81392</a> )



The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada



The vibrant M and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources. © 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved.

## Recommended Cell Assembly

Ag, AgCl | 3 M KCl || 0.3 M NH<sub>4</sub>NO<sub>3</sub> || sample solution || liquid membrane | 0.1 M KCl | AgCl, Ag

## Electrode Characteristics and Function

Selectivity coefficients  $\log K_{K,M}^{Pot}$  as obtained by the separate solution method (0.1 M solutions of the chlorides).

$\log K_{K,Li}^{Pot}$	-4.3
$\log K_{K,Na}^{Pot}$	-4.2

Slope of linear regression: 58 mV (1.8·10<sup>-5</sup> to 10<sup>-1</sup> M K<sup>+</sup>)  
Detection limit: 3·10<sup>-6</sup> M K<sup>+</sup>  
Lifetime: >1 month

## Ion-selective Field Effect Transistors

### Application 1 and Sensor Type<sup>17,18,19</sup>

Assay of K<sup>+</sup> activity with Urushi matrix ion-selective field effect transistors of good durability based on Potassium ionophore I.

### Recommended Membrane Composition

- 0.50 wt% Potassium ionophore I ([60403](#))
- 49.25 wt% Bis(2-ethylhexyl) phthalate ([80030](#))
- 0.25 wt% Potassium tetrakis((4-chlorophenyl)borate) ([60591](#))
- 50.00 wt% Urushi (polymer from lacquer tree)

### Electroanalytical Characteristics

Selectivity coefficients  $\log K_{K,M}^{Pot}$  as obtained by the mixed solution method.

$\log K_{K,Mg}^{Pot}$	-4.6	$\log K_{K,Ca}^{Pot}$	-4.5
$\log K_{K,Na}^{Pot}$	-3.3	$\log K_{K,NH_4}^{Pot}$	-1.7

Slope of linear regression: 53 mV (10<sup>-4</sup> to 1 M K<sup>+</sup>)

### Application 2 and Sensor Type<sup>18,19,20</sup>

Assay of K<sup>+</sup> activity with silicon rubber matrix ion-selective field effect transistors based on Potassium ionophore I.

### Recommended Membrane Composition

- 3.0 wt% Potassium ionophore I ([60403](#))
- 88.0 wt% Siloprene K 1000 ([85417](#))
- 9.0 wt% Siloprene crosslinking agent K 11 ([85418](#))

### Electroanalytical Characteristics

Selectivity coefficients  $\log K_{K,M}^{Pot} < -3.7$  as obtained by the fixed interference method (0.1 M of the chlorides).

Lifetime: ~2 months



## Ion-selective Conductometric Microsensors

### Application<sup>21</sup>

Assay of K<sup>+</sup> activity with ion-selective conductometric microsensors (ISCOM). Detection is accomplished by measurement of the bulk conductance of the solvent polymeric membrane based on Potassium ionophore I.

### Recommended Membrane Composition

- 5.0 wt% Potassium ionophore I ([60403](#))
- 30.0 wt% Poly(vinyl chloride) high molecular weight ([81392](#))
- 65.0 wt% 2-Nitrophenyl octyl ether ([73732](#))

### Electroanalytical Characteristics

Selectivity coefficients  $\log K_{K,M}^{Pot}$  as obtained by the mixed solution method.

$\log K_{K,Li}^{Pot}$	-3.0	$\log K_{K,Cs}^{Pot}$	-0.5
$\log K_{K,Na}^{Pot}$	-2.6	$\log K_{K,Rb}^{Pot}$	0.6
$\log K_{K,NH_4}^{Pot}$	-1.3		

Membrane specific conductivity:	0.5 to 5 $\mu$ S/cm
Frequency range:	0.5 to 20 kHz
Membrane thickness:	~2 mm
Detection limit:	>10 <sup>-6</sup> M K <sup>+</sup>
Response time:	~1 s

## Optical Transduction

### Application 1 and Sensor Type<sup>22,23,24</sup>

Assay of K<sup>+</sup> activity in aqueous pH buffered solutions and in diluted blood plasma with solvent polymeric optode membranes based on Chromoionophore I (ETH 5294) and Potassium ionophore I.

### Recommended Membrane Composition

- 0.48 wt% Chromoionophore I ([27086](#))
- 1.00 wt% Potassium ionophore I ([60403](#))
- 0.44 wt% Potassium tetrakis(4-chlorophenyl)borate ([60591](#))
- 66.05 wt% Bis(1-ethylhexyl)sebacate ([84818](#))
- 32.03 wt% Poly(vinyl chloride) high molecular weight ([81392](#))

### Recommended pH Buffer

0.16 M sodium acetate, adjusted with acetic acid to pH 5.1 for recording the calibration curve to pH 5.5 for diluting blood plasma samples.<sup>25</sup>

### Optode Characteristics and Function

Selectivity coefficients  $\log K_{K,M}^{Opt}$  as obtained by the fixed interference method in pH buffered solutions.

$\log K_{K,Na}^{Opt}$	-3.5	$\log K_{K,Ca}^{Opt}$	-3.7
$\log K_{K,Mg}^{Opt}$	-4.0	$\log K_{K,Li}^{Opt}$	-3.7

### Application 2 and Sensor Type<sup>26</sup>

Assay of K<sup>+</sup> activity in aqueous pH buffered solutions with solvent polymeric fluorescent optode membranes based on Chromoionophore I (ETH 5294) and Potassium ionophore I. LEDs or diode lasers may be used as light sources.

### Recommended Membrane Composition

- 2.98 wt% Chromoionophore I ([27086](#))
- 13.43 wt% Potassium ionophore I ([60403](#))
- 2.98 wt% Potassium tetrakis(4-chlorophenyl)borate ([60591](#))
- 44.78 wt% Bis(2-ethylhexyl)sebacate ([84818](#))
- 17.91 wt% WM-3 plasticizer
- 17.92 wt% Poly(vinyl chloride) high molecular weight ([81392](#))



The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada



## Recommended pH Buffer

0.1 M TRIS at pH 7.38

## Optode Characteristics

Membrane range:	5 $10^{-6}$ to $10^{-1}$ K <sup>+</sup>
Membrane thickness:	~2 $\mu$ m
Detection limit:	~5 $10^{-6}$ M K <sup>+</sup>
Response time:	~1 to 3 min

## Application 3 and Sensor Type<sup>27</sup>

Assay of ammonia gas in aqueous solutions with solvent polymeric optode membranes based on Chromoionophore III (ETH 5350) and Potassium ionophore I. The optode membrane is covered with a gas-permeable membrane.

## Recommended Membrane Composition

3.50 wt%	Potassium ionophore I ( <a href="#">60403</a> )
1.52 wt%	Chromoionophore III ( <a href="#">27088</a> )
1.52 wt%	Potassium tetrakis(4-chlorophenylphenyl)borate ( <a href="#">60591</a> )
62.31 wt%	Bis(2-ethylhexyl)sebacate ( <a href="#">84818</a> )
31.15 wt%	Poly(vinyl chloride) high molecular weight ( <a href="#">81392</a> )

## Optode Characteristics

Selectivity coefficients  $\log K_{NH_3,G}^{Opt}$  as obtained by the separate solution method.

$\log K_{NH_3,MeNH_2}^{Opt}$	-0.4	$\log K_{NH_3,Me_2NH}^{Opt}$	-2.0
$\log K_{NH_3,EtNH_2}^{Opt}$	-1.7	$\log K_{NH_3,Me_3N}^{Opt}$	-2.9
$\log K_{NH_3,PrNH_2}^{Opt}$	-1.9		

Reproducibility:  $10^{-4}$  NH<sub>3</sub> ~3.2%;  $10^{-3}$  M NH<sub>3</sub> ~1.6%

<sup>1</sup> D. Ammann, P. Anker, E. Metzger, U. Oesch, W. Simon, in: Ion Measurements in Physiology and Medicine, Eds. M. Kessler, D.K. Harrison, J. Höper, Springer-Verlag, Berlin, Heidelberg 120 (1985)

<sup>2</sup> Relationship between activity and concentration measurements of plasma potassium. D. M. Band, J. Kratochvil, P. A. Poole Wilson, T. Treasure, Analyst 103, 246 (1978).

<sup>3</sup> Flexible valinomycin electrodes for on-line determination of intravascular and myocardial K<sup>+</sup>. J. L. Hill, L. S. Gettes, M.R. Lynch, N.C. Hebert, Am. J. Physiol. 235, H455 (1978).

<sup>4</sup> Determination of K<sup>+</sup> in diluted and undiluted urine with ion-selective electrodes. H.-B. Jenny, C. Riess, D. Ammann, B. Magyar, R. Asper, W. Simon, Mikrochim. Acta 74, 309. (1980).

<sup>5</sup> On-line continuous potentiometric measurement of potassium concentration in whole blood during open-heart surgery. H. F. Osswald, R. Asper, W. Dimai, W. Simon, Clin. Chem. 25, 39 (1979).

<sup>6</sup> Ion selective electrodes in clinical chemistry. A. Lewenstam, Anal. Proc. 28, 106 (1991)

<sup>7</sup> U. Oesch, P. Anker, D. Ammann, W. Simon, in: Ion-Selective Electrodes, ed. E. Pungor, I. Buzás, Akadémiai Kiadó, Budapest (1985).

<sup>8</sup> Lifetime of neutral-carrier-based liquid membranes in aqueous samples and blood and the lipophilicity of membrane components. O. Dinten, U. E. Spichiger, N. Chaniotakis, P. Gehrig, B. Rusterholz, W. E. Morf, W. Simon, Anal. Chem. 63, 596 (1991).

<sup>9</sup> Determination of [K<sup>+</sup>] in blood serum with a valinomycin-based silicone rubber membrane of universal applicability to body fluids. P. Anker, H.-B. Jenny, U. Wuthier, R. Asper, D. Ammann, W. Simon, Clin. Chem. 29, 1447 (1983).

<sup>10</sup> D. Ammann, P. Anker, H.-B. Jenny, W. Simon, in: Ion-Selective Electrodes, Eds. E. Pungor, I. Buzás, Akadémiai Kiadó, Budapest (1981).

<sup>11</sup> Ion-selective microchemical sensors with reduced preconditioning time: Membrane biostability studies and applications in blood analysis. V. V. Cosofret, M. Erdösy, E. Lindner, T.A. Johnson. R. P. Buck, W. J. Kao, M. R. Neuman, J. M. Anderson, Anal. Lett. 27, 3039 (1994).

<sup>12</sup> Microfabricated sensor arrays sensitive to pH and K<sup>+</sup> for ionic distribution Measurements in the Beating Heart. V. V. Cosofret, M. Erdösy, T. A. Johnson, R. P. Buck, R. B. Ash, M. R. Neuman, Anal. Chem. 67, 1647 (1995).

<sup>13</sup> Electroanalytical and biocompatibility studies on carboxylated poly(vinyl chloride) membranes for microfabricated array sensors. V. V. Cosofret, M. Erdösy, R. P. Buck, W. J. Kao, J. M. Anderson, E. Lindner, M. R. Neuman, Analyst 119, 2283 (1994).



- <sup>14</sup> Electroanalytical and Surface Characterization of Encapsulated Implantable Membrane Planar Microsensors. V. V. Cosofret, M. Erdösy, T. A. Johnson, B. A. Bellinger, R. P. Buck, R. B. Ash, M. R. Neuman, *Anal. Chim. Acta* 314, 1 (1995).
- <sup>15</sup> Microfibre matrix-supported ion-selective PVC membranes. M. Knoll, K. Cammann, C. Dumschat, M. Borchardt, G. Högg, *Sens. Actuators B20*, 1 (1994).
- <sup>16</sup> Disposable ion-selective electrodes. M. Borchardt, C. Dumschat, K. Cammann, M. Knoll, *Sens. Actuators B24-25*, 721 (1995).
- <sup>17</sup> Urushi matrix sodium, potassium, calcium and chloride-selective field-effect transistors. S. I. Wakida, M. Yamane, K. Higashi, K. Hiro, Y. Ujihara, *Sens. Actuators B1*, 412 (1990).
- <sup>18</sup> J. Janata, R. Huber, *Ion-Selective Electrode Rev.* 1, 31 (1979).
- <sup>19</sup> The operation of an ISFET as an electronic device. P. Bergveld, *Sensors and Actuators* 1, 17 (1981).
- <sup>20</sup> New membrane materials for potassium-selective ion-sensitive field-effect transistors. P. Van der Wal, M. Skowronska-Ptasinska, A. van den Berg, P. Bergveld, E. J. R. Sudhölter, D. N. R. Reinhoudt, *Anal. Chim. Acta* 231, 41 (1990).
- <sup>21</sup> Ion-selective conductometric microsensors based on the phenomenon of specific salt extraction. A. A. Shul'ga, B. Ahlers, K. Cammann, *J. Electroanal. Chem.*, 395, 305 (1995).
- <sup>22</sup> K. Seiler, *Ion-selective Optode Membranes*, monograph, describing theory, preparation and application of ion-selective optode membranes as well as recent developments in this field. With 237 references. published by Fluka Chemie GmbH, Buchs, Switzerland (1993); K. Seiler, *Ionenselektive Optodenmembranen*, dt. Monographie, herausgegeben von Fluka Chemie GmbH, Buchs, Switzerland (1993)
- <sup>23</sup> Optical quantification of sodium, potassium and calcium ions in diluted human plasma based on ion-selective liquid membranes. U.E. Spichiger, K. Seiler, K. Wang, G. Suter, W.E. Morf, W. Simon, *International Congress on Optical Sciences and Engineering*, The Hague, The Netherlands, 11-15 March 1991. *SPIE Proceedings* Nr. 1510, pp. 118-130 (1991).
- <sup>24</sup> Optodes in clinical chemistry: potential and limitations. U.E. Spichiger, D. Freiner, E. Bakker, T. Rosatzin, W. Simon, *Sens. Actuators B11*, 263 (1993).
- <sup>25</sup> D.D. Perrin, B. Dempsey, *Buffers for pH and Metal Ion Control*. Chapman & Hall, London, New York (1983).
- <sup>26</sup> An Enantio-Selective Optode for the  $\beta$ -Blocker Propranolol. H. He, G. Uray, O. Wolfbeis, *Proc. SPIE-Int. Soc. Opt. Eng.* 1368, 165 (1990).
- <sup>27</sup> Ammonia-gas-selective optical sensors based on neutral ionophores. S. Ozawa, P. C. Hauser, S. S. S. Tan, W. E. Morf, W. Simon, *Anal. Chem.* 63, 640 (1991).

