

**Product No. I-4770**  
**Lot 044H0297**

**Anti-Human Interleukin-9**  
**Developed in Goat**  
IgG Fraction of Antiserum

Anti-Human Interleukin-9 is developed in goat using human, recombinant interleukin-9 (rhIL-9), expressed in the insect cell line *Sf* 21, as the immunogen. The product is purified by Protein G affinity chromatography. Goat Anti-Human IL-9 is provided lyophilized from phosphate buffered saline, pH 7.4, to which no preservatives have been added.

#### **Description**

Interleukin-9 is a T cell-derived<sup>1</sup> glycoprotein having multiple functions on cells of the lymphoid, myeloid, and mast cell lineages.<sup>2</sup> Anti-Human IL-9 neutralizes the bioactivity of rhIL-9 but not that of recombinant mouse IL-9. In indirect ELISA and immunoblotting, this antibody displays less than 20% cross-reactivity with recombinant, mouse IL-9. The product shows no cross-reactivity with other tested cytokines using indirect ELISA.

#### **Performance**

Anti-Human IL-9 was tested for its ability to neutralize the bioactivity of rhIL-9 in a cell proliferation assay using the factor-dependent human cell line MO7e.<sup>3</sup> The ND<sub>50</sub> of the antibody is defined as the concentration of antibody resulting in a one-half maximal inhibition of bioactivity of rhIL-9, which is present at five times its own EC<sub>50</sub> (the concentration of rhIL-9 producing a one-half maximal bioactivity without antibody). In this bioassay, rhIL-9 was pre-incubated with various dilutions of the antibody for 1 hour at 37°C in a 96-well microtiter plate. Then, MO7e cells were added to each well to give a final concentration of 1 x 10<sup>5</sup> cells/ml in 0.1 ml containing 2.5 ng/ml rhIL-9. This was incubated for 72 hours at 37°C in a 5% CO<sub>2</sub> humidified incubator and then pulsed for 4 hours with <sup>3</sup>H-thymidine. Cells were harvested onto glass filters and the <sup>3</sup>H-thymidine incorporation into DNA was measured.

#### **Product Information**

Mass/vial:	1 mg
Immunogen:	Human, recombinant IL-9 (rhIL-9)
Host animal:	Goat
Formulation:	Lyophilized from PBS without additives
Endotoxin:	<10 ng/vial by LAL method
Bioactivity:	ND <sub>50</sub> = 24.8 µg/ml
Indirect ELISA:	0.5 µg/ml detects 0.3 ng/well of rhIL-9.
Indirect Immunoblotting:	1 µg/ml antibody detects rhIL-9 at 5 ng/lane under reducing conditions and at 20 ng/lane under non-reducing conditions.
Sterility:	0.2 µm-filtered, aseptic fill

#### **Reconstitution and Use**

To one vial of lyophilized powder, add 1 ml of sterile-filtered PBS to produce a 1 mg/ml stock solution of Anti-Human IL-9. If aseptic technique is used, no further filtration should be needed for use in cell culture environments.

#### **Storage**

Prior to reconstitution, store at -20°C for 6 months. Reconstituted product may be stored at 0-5°C for up to one month. For prolonged storage, freeze in working aliquots at -20°C. Avoid repeated freezing and thawing.

#### **References**

1. Yang, Y., et al., *Blood*, **74**, 1880 (1989).
2. Renauld, J., et al., *Int. Rev. Exp. Path.*, **34A**, 99 (1993).
3. Avanzi, G., et al., *Br. J. Haematol.*, **69**, 359 (1988).

Prepared for Sigma.

Sigma warrants that its products conform to the information contained in this and other Sigma publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Issued 06/94.