



MOUSE ANTI-ECHOVIRUS BLEND MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB9670
LOT NUMBER:	
QUANTITY:	100 μ L
SPECIFICITY:	Reacts with the Echovirus 4, 6, 9, 11, 30 and 34.
ISOTYPE:	IgG
APPLICATIONS:	IFA at 1:100. Dilute with buffer pH 7.5-8.0 containing protein or other stabilizing medium. Final working dilutions must be determined by end user.
CROSS REACTIVITIES:	No cross-reactivity to CMV, HSV 1 & 2, Para 1, 2, & 3, Poliovirus 1 & 3 or Mumps.
FORMAT:	Ascites, with 0.1% sodium azide.
STORAGE/HANDLING:	Maintain at -20°C in convenient, undiluted aliquots for up to 12 months. Avoid repeated freeze/thaw cycles.
REFERENCE:	<i>Current Microbiology</i> (1994) 28 :49-52

For research use only; not for use as a diagnostic.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ l or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*