



QUALITY CONTROL RANGES

HUMAN FIBROBLAST GROWTH FACTOR-23

ELISA

CATALOG # EZHFGF23-32K

STANDARD, LYOPHILIZED

CATALOG # E8032-K

STANDARD	HUMAN FGF-23 VALUE
LOT NUMBER: FGF32-10K	6869 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED

CATALOG # E6032-K

Reconstitute each vial with **0.5 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	FGF-110	66 – 138	pg/mL
QC 2	FGF-210	512 – 1064	pg/mL

E6032-K: FGF-110.210 (11 JUNE 2020)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.