



QUALITY CONTROL RANGES
MILLIPLEX[®] MAP Human Circulating Cancer Biomarker Magnetic Bead Panel 2
Kit Catalog# HCCBP2MAG-58K
Control Catalog # HCC-6058-2

<u>Hormone</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
ECM1	Control 1	HCC2 – 106A	0.08 - 0.17	ng/mL
	Control 2	HCC2 – 206A	3.7 - 7.6	ng/mL
Vitamin D BP	Control 1	HCC2 – 106A	2.8 - 5.7	ng/mL
	Control 2	HCC2 – 206A	112 - 232	ng/mL
Antithrombin III	Control 1	HCC2 – 106A	3.3 - 6.8	ng/mL
	Control 2	HCC2 – 206A	121 - 250	ng/mL
CFH	Control 1	HCC2 – 106A	5.5 - 11	ng/mL
	Control 2	HCC2 – 206A	207 - 429	ng/mL

HCC-6058-2: HCC2-106A.206A (05 November 2021)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.