

Product Information Sheet

Hanks' Balanced Salts

Modified, powder, suitable for cell culture

H1387

Product Description

Although there have been many modifications to the original formulas in efforts to produce fully defined media, salt solutions still play an important role in tissue culture. A salt solution's basic function, to maintain the pH and osmotic balance in the medium and to provide the cells with water and essential inorganic ions, is as valuable today as when it was first developed a century ago.

Components	g/L
Calcium Chloride (anhydrous)	0.1396
Magnesium Sulfate (anhydrous)	0.09767
Potassium Chloride	0.4
Potassium Phosphate Monobasic (anhydrous)	0.06
Sodium Chloride	8.0
Sodium Phosphate Dibasic (anhydrous)	0.04788
D-Glucose	1.0

Precautions and Disclaimer

For research or further manufacturing uses only.
Not intended for direct use in humans or animals.

Preparation Instructions

Powdered salts are hygroscopic and should be protected from moisture. The entire contents of each package should be used immediately after opening. Preparing a concentrated salt solution is not recommended as precipitates may form. Supplements can be added prior to filtration or introduced aseptically to sterile salt solution.

1. Measure out 90% of final required volume of water. Water temperature should be 15-20 °C.
2. While gently stirring the water, add the powdered medium. Stir until dissolved. Do NOT heat.
3. Rinse original package with a small amount of water to remove all traces of powder. Add to solution in step 2.
4. To the solution in step 3, add 0.35 g sodium bicarbonate or 4.7 mL of sodium bicarbonate solution [7.5% w/v] for each liter of final volume of medium being prepared. Stir until dissolved.
5. While stirring, adjust the pH of the medium to 0.1-0.3 pH units below the desired pH since it may rise during filtration. The use of 1 N HCl or 1 N NaOH is recommended.
6. Add additional water to bring the solution to final volume.
7. Sterilize immediately by filtration using a membrane with a porosity of 0.22 microns.
8. Aseptically dispense medium into sterile container.

Storage and Stability

Store the dry powdered salts at 2-8 °C under dry conditions and liquid medium at 2-8 °C in the dark. Deterioration of the powdered medium may be recognized by any or all of the following: color change, granulation/clumping, insolubility. Deterioration of the liquid medium may be recognized by any or all of the following: pH change, precipitate or particulates, cloudy appearance and color change. The nature of supplements added may affect storage conditions and shelf life of the medium. Product label bears expiration date.

Materials Required (Not Provided)

- Water for tissue culture use [Cat. No. W3500]
- Sodium Bicarbonate [Cat. No. S5761] or Sodium Bicarbonate Solution, 7.5% [Cat. No. S8761]
- 1 N Hydrochloric Acid [Cat. No. H9892]
- 1 N Sodium Hydroxide [Cat. No. S2770]
- Medium additives as required

References

1. Hanks, J. (1976) Hanks' Balanced Salt Solution and pH Control. Tissue Culture Association Manual. 3, 3.

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

H1387pis Rev 09/22

**MILLIPORE
SIGMA**