



## QUALITY CONTROL RANGES

### HUMAN FIBROBLAST GROWTH FACTOR-23

#### ELISA

#### CATALOG # EZHFGF23-32K

#### STANDARD, LYOPHILIZED

#### CATALOG # E8032-K

STANDARD	HUMAN FGF-23 VALUE
LOT NUMBER: FGF32-12K	6869 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5 mL** distilled or deionized water.

#### QUALITY CONTROLS 1 AND 2, LYOPHILIZED

#### CATALOG # E6032-K

Reconstitute each vial with **0.5 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	FGF-112	57 – 119	pg/mL
QC 2	FGF-212	469 – 975	pg/mL

E6032-K: FGF-112.212 (30 May 2024)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) [www.emdmillipore.com](http://www.emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA