



- **EZLabpure™** Aspirator Bottles are manufactured using high purity virgin USP Class VI Medical grade polypropylene
- Very strong and durable for extended user over time. They are equipped with inner-cap for a tight seal. Ideal for containing liquid powders and granulates. Comes with high quality & durable stopcock or tap. These are used for Packing, Shipping & Storage purpose
- All Can be used at a temperature below 0°C
- All plastic consumables are autoclavable at 121°C for 20 minutes, 15 psi/1atm. Do not use a dry cycle as this may cause deformation or weakening of the plastic. Allow temperature within the autoclave to return to at least 50°C before removing product

Caution

- Do not store strong oxidizing agents
- Do not place any plastic labware in a flame
- Do not mix any chemicals that may result in a thermic reaction, which can cause product failure
- Consult local fire codes prior to storage of flammable liquids in carboys/bottles/aspirators
- Consult OSHA guidelines prior to handling or lifting filled carboys/bottles/aspirators

Foxx Life Sciences offers a premium line of high quality, polypropylene aspirator bottles. **EZLabpure™** Aspirator Bottles come with strong easy-grip handle and Spigot. Our versatile bottles are leak-proof and aim to take full advantage of storage efficiency at minimal effort. Overall functionality is available for many applications, in five liter, ten liter, or twenty liter sizes, with spigots and a high density poly ethylene (HDPE) cap.

CE	Yes
Autoclavability	Yes, From 0°C to 121°C
Translucence	Yes
Resistant	To be effected by radiations
Sterilization	Gamma
Leak Proof	Yes

- 1 HDPE (High density poly ethylene) aspirator cap
- 2 Handle design to make handling, lifting, and pouring bottle easier
- 3 Important instructions and warnings for the **EZLabpure™** Aspirator Bottle
- 4 High quality and durable stopcock with on and off engraving





Product Number	Description	Capacity (L)	Cap Size (mm)	Autoclavable
155-J897-FLS	5L PP Aspirator Bottle with White Cap and Spigot	5	68	Yes
155-K897-FLS	10L PP Aspirator Bottle with White Cap and Spigot	10	63.5	Yes
155-L897-FLS	20L PP Aspirator Bottle with White Cap and Spigot	20	63.5	Yes



EZLabpure™ Round Carboys

Foxx Life Sciences offers a line of round polypropylene (PP) laboratory carboys with the patented VersaCap® technology and a wide mouth option. These durable, rugged and shatterproof round carboys are a great alternative to glass carboys. They can be used for mixing, storing, and dispensing liquids in volumes from ten to twenty liters. The round carboys are available with and without a spigot. EZLabpure™ Round Carboys are compatible with Foxx's patented VersaCap® that is leak-proof, has replaceable platinum-cured silicone gaskets, and provides the ultimate flexibility. The open VersaCap® allows a variety of configurations by changing the adapter insert, enabling users to quickly change the top connections of each carboy or bottle to match their application needs.