



**QUALITY CONTROL RANGES**  
**MILLIPLEX<sup>®</sup> MAP Human Apolipoprotein Magnetic Bead Panel Kit**  
**Catalog # AP0MAG-62K**

**Standard Catalog # HAP0-8062**  
**Lot# 862M-8K**

<b>Apolipoproteins</b>	<b>Standard Value</b>		<b>Apolipoproteins</b>	<b>Standard Value</b>
Apo AI	54560 ng/mL		Apo AII	5772 ng/mL
Apo CIII	4332 ng/mL		Apo B	28840 ng/mL
Apo E	2156 ng/mL		Apo CII	1990 ng/mL

**Control Catalog # HAP0-6062**  
**Lot# AP0M-108 and AP0M-208**

<u>Apolipoproteins</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
Apo-AI	Control I	AP0M-108	241 - 517	ng/mL
	Control II	AP0M-208	1238 - 2570	ng/mL
Apo-CIII	Control I	AP0M-108	19 - 39	ng/mL
	Control II	AP0M-208	83 - 172	ng/mL
Apo-E	Control I	AP0M-108	4.6 - 14	ng/mL
	Control II	AP0M-208	35 - 75	ng/mL
Apo-AII	Control I	AP0M-108	15 - 30	ng/mL
	Control II	AP0M-208	90 - 187	ng/mL
Apo-B	Control I	AP0M-108	72 - 149	ng/mL
	Control II	AP0M-208	385 - 800	ng/mL
Apo-CII	Control I	AP0M-108	7.4 - 15	ng/mL
	Control II	AP0M-208	37 - 77	ng/mL

HAP0-6062: AP0M-108.208 (11 March 2020)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • [emdmillipore.com](http://emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.