



RABBIT ANTI-PROHORMONE CONVERTASE PC2 POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB15610
LOT NUMBER:	
QUANTITY:	100 μ L
SPECIFICITY:	Reacts with Prohormone Convertase PC2
IMMUNOGEN:	Amino acids 112-159 of human preproconvertase PC2.
APPLICATIONS:	Immunohistochemistry/immunocytochemistry: 1:50-1:250 Western blot: 1:50 RIA: 1:8,000 Optimal working dilutions must be determined by end user.
SPECIES REACTIVITIES:	Human, rat and mouse.
FORMAT:	Rabbit antiserum
PRESENTATION:	Liquid, no preservatives.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 12 months after date of receipt. Avoid repeated freeze/thaw cycles.
REFERENCES:	<i>J. Biol. Chem.</i> (1995) 270 :10136-10146. <i>Mol. Endocrinol.</i> (1996) 10 :331-341.

Important Note: *During shipment, small volumes of antibody will occasionally become entrapped in the seal of the product vial. For antibodies with volumes of 200 μ l or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial on a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.

©2002 - 2009: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209
www.millipore.com