



QUALITY CONTROL RANGES
MILLIPLEX® MAP Canine Pituitary Magnetic Panel
Catalog # CANPIT-96K
Control Catalog # CANP-6096 Lot# CANP-103 and
CANP-7096 Lot# CANP-203

<u>Hormone</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
ACTH	CANP-6096	CANP-103	31 - 64	pg/mL
	CANP-7096	CANP-203	252 - 524	pg/mL
FSH	CANP-6096	CANP-103	30 – 62	ng/mL
	CANP-7096	CANP-203	254 - 528	ng/mL
GH	CANP-6096	CANP-103	0.15 – 0.31	ng/mL
	CANP-7096	CANP-203	1.3 – 2.7	ng/mL
TSH	CANP-6096	CANP-103	47 – 98	pg/mL
	CANP-7096	CANP-203	397 – 825	pg/mL
LH	CANP-6096	CANP-103	1.1 – 2.3	ng/mL
	CANP-7096	CANP-203	9.5 - 20	ng/mL
Prolactin	CANP-6096	CANP-103	15 – 32	ng/mL
	CANP-7096	CANP-203	128 - 266	ng/mL

Note: The Quality Control Ranges are generated with overnight assay format using serum matrix provided in the kit.

Quality Control values in culture media are not tested.

CANP-6096 Lot# CANP-103 and CANP-7096 Lot# CANP-203 (22 Feb 2022)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA