

**RAT ANTI-NEUROFILAMENT H
MONOCLONAL ANTIBODY**

CATALOG NUMBER: MAB5448

LOT NUMBER:

QUANTITY: 1 mL

SPECIFICITY: Neurofilament H (NF-H), carboxy terminal (phosphate dependent (+)).

ISOTYPE: IgG

CLONE: TA51

APPLICATIONS: Western blot
Immunohistochemistry
Optimal working dilutions must be determined by end user.

SPECIES REACTIVITY: Human and mouse. Other species have not yet been tested.

FORMAT: Culture supernatant from hollow fiber system.

PRESENTATION: Liquid. Contains no preservative.

STORAGE/HANDLING: Maintain at -20°C in undiluted aliquots for up to 12 months after date of receipt. Avoid repeated freeze/thaw cycles.

REFERENCES: Trojanowski, J.K., et al., *American J. Pathology* (1989) **135**:747-58.
Tu, P., et al., *PNAS.USA* (1996) **93**:3155-3160.
Paulson, H.L., et al., *Neuron* (1997) **19**:333-344.
Myrayama, H., et al., *American J. Pathology* (1999) **155**:877-885.
De Girolamo, L.A., et al., *J. Neurochemistry* (2000) **75**:133-140.

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

For research use only; not for use as a diagnostic.

© 2002-2004: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (909) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com