

Product Change Notification

Ultracel® 100K Regenerated Cellulose Membrane Disc Products

The manufacture of Ultracel® PL-100 (100,000 MWCO) membrane discs has transferred from our U.S. manufacturing facility in Bedford, MA. to our state-of-the-art production line located in Jaffrey, NH, in order to increase manufacturing efficiency.

As part of this transfer, the polyolefin backing for the Ultracel® PL-100 membrane has changed to a similar polyolefin backing manufactured by a different supplier. This backing is being used for other ultrafiltration membranes manufactured in Jaffrey, NH and will provide consistency and efficiency across the ultrafiltration membrane product line.

The manufacture and sales of Ultracel® YM-100 (100,000 MWCO) membrane discs have been discontinued and these products have been replaced by Ultracel® PL-100 products with equivalent filtration properties. Replacement catalogue numbers are listed below.

Discontinued Cat. No.	Description	Replacement Cat. No.
14412	100K MWCO Ultrafiltration Discs, 25 mm, 10/pk	PLHK02510
14422AM	100K MWCO Ultrafiltration Discs, 44.5 mm, 10/pk	PLHK04310
14432	100K MWCO Ultrafiltration Discs, 63.5 mm, 10/pk	PLHK06210
14442AM	100K MWCO Ultrafiltration Discs, 76 mm, 10/pk	PLHK07610
14451	100K MWCO Ultrafiltration Discs, 90 mm, 5/pk	PLHK09005
14461	100K MWCO Ultrafiltration Discs, 150 mm, 5/pk	PLHK15005

All changes have been fully validated and there are no changes to product claims. However, since we cannot evaluate specific customer applications, we recommend that you check your application to confirm that no modifications are required.

If you have questions regarding this change, please contact the office nearest you. In the U.S., call 1-800-221-1975. Outside the U.S., go to our web site at www.millipore.com/offices for up-to-date worldwide contact information. You can also visit the tech service page on our web site at www.millipore.com/techservice.

