



QUALITY CONTROL RANGES

MILLIPLEX[®] MAP Human Complement Magnetic Bead Panel 1 Kit Catalog# HCMP1MAG-19K Control Catalog # HCMP1-6019-1

<u>Analyte</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
Complement C2	Control 1	HCM1-106	11 – 24	ng/mL
	Control 2	HCM1-206	43 – 90	ng/mL
Complement C4b	Control 1	HCM1-106	9.8 – 20	ng/mL
	Control 2	HCM1-206	47 – 97	ng/mL
Complement C5	Control 1	HCM1-106	21 – 43	ng/mL
	Control 2	HCM1-206	83 – 172	ng/mL
Complement C5a	Control 1	HCM1-106	29 – 60	pg/mL
	Control 2	HCM1-206	116 – 240	pg/mL
Complement Factor D	Control 1	HCM1-106	0.54 – 1.1	ng/mL
	Control 2	HCM1-206	2.1 – 4.3	ng/mL
MBL	Control 1	HCM1-106	1.1 – 2.2	ng/mL
	Control 2	HCM1-206	4.1 – 8.5	ng/mL
Complement Factor I	Control 1	HCM1-106	6.1 – 13	ng/mL
	Control 2	HCM1-206	23 – 49	ng/mL

HCMP1-6019-1: HCM1-106.206 (11 August 2022)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA