



QUALITY CONTROL RANGES

MILLIPLEX[®] MAP Human Neurodegenerative Disease Magnetic Bead Panel 1 Catalog # HNDG1MAG-36K

Control Catalog # HNDG1-6036-1 Control Lot # HNDG1-107A.207A

<u>Hormone</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
APOA-1	Control I	HNDG1-107A	12 - 25	ng/mL
	Control II	HNDG1-207A	66 - 137	ng/mL
APOC-3	Control I	HNDG1-107A	0.18 - 0.38	ng/mL
	Control II	HNDG1-207A	0.92 - 1.9	ng/mL
APOE	Control I	HNDG1-107A	0.75 - 1.7	ng/mL
	Control II	HNDG1-207A	3.7 - 7.7	ng/mL
HCFH	Control I	HNDG1-107A	5.0 - 10	ng/mL
	Control II	HNDG1-207A	27 - 56	ng/mL
HCC3	Control I	HNDG1-107A	0.78 - 1.6	ng/mL
	Control II	HNDG1-207A	4.0 - 8.4	ng/mL
HA2MG	Control I	HNDG1-107A	7.8 - 16	ng/mL
	Control II	HNDG1-207A	41 - 85	ng/mL

HNDG1-6036-1: HNDG1-107A.207A (28 July 2021)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.