

M Air T™ Sampling Head

■ A sampling system especially designed for barrier environments



Precisely sample airborne microorganisms in barrier environments and aseptic manufacturing areas with the M Air T Sampling Head. The unit includes the unique M Air T microperforated sieve design and double-wrapped, pre-sterilized M Air T agar cassettes.

Easy-Loading Agar Cassettes

M Air T agar cassettes are easier to handle than agar strips — especially when wearing gloves, which means there's much less risk of inadvertent contamination.

Double-Wrapped, Pre-sterilized Media

M Air T cassettes are available in double-wrap packages to minimize the risk of contamination within the barrier environment. In addition, the packages, which contain four, gamma-irradiated media cassettes, are qualified and validated for impermeability to sanitizing agents.

Sampling Above the "Turbulence Zone"

Inside a barrier environment with vertical laminar flow, some turbulence can exist at and just above the work surface. The M Air T Sampling Head is supported on a stand above this "turbulence zone."

Autoclave-Compatible

The M Air T Sampling Head is manufactured from stainless steel with sanitary Tri-Clover® connections.

Compatible with a Variety of Pumps

The M Air T Sampling Head is compatible with the pumps commonly used for sampling barrier environments. Please contact the nearest Millipore office for additional information.

MILLIPORE

Ordering Information

Description	Qty/Pk	Catalogue No.
M Air T Sampling Head	1	ATBH EAD 01
M Air T Agar Cassettes, double wrapped, gamma irradiated	24	ATSM TTB 24
M Air T Calibration Kit	1	ATAC ALO 01
M Air T Validation Protocol	1	ATAS OVG 01

To Place an Order or Receive Technical Assistance

For additional information call your nearest Millipore office:

In the U.S. and Canada, call toll-free **1-800-MILLIPORE (1-800-645-5476)**

In the U.S., Canada and Puerto Rico, fax orders to **1-800-MILLIFX (1-800-645-5439)**

On the Internet

<http://www.millipore.com>

E-mail: tech_service@millipore.com



Now you can buy Millipore products online @

www.millipore.com/purecommerce

MILLIPORE

Millipore and M Air T are trademarks of Millipore Corporation or an affiliated company.

Tri-Clover is a trademark of Ladish Company, Inc.

Lit. No. PF1533EN00 6/00 00-156

Printed in U.S.A.

Copyright © 2000 Millipore Corporation or an affiliated company. All rights reserved.

Millipore Worldwide

AUSTRALIA

Tel. 1 800 222 111
or (02) 9888 8999
Fax (02) 9878 0788

AUSTRIA, CENTRAL EUROPE, C.I.S., AFRICA, MIDDLE EAST AND GULF
Tel. (43) 1 877-8926
Fax (43) 1 877-1654

BALTIC COUNTRIES
Tel. +358 9 804 5110
Fax +358 9 256 5660

BELGIUM AND LUXEMBOURG
Tel. (02) 726 88 40
Fax (02) 726 98 84

BRAZIL
Tel. (011) 548-7011
Fax (011) 548-7923

CANADA
Tel. 1-800-645-5476
Fax 1-800-645-5439

CHINA, PEOPLE'S REPUBLIC OF
Beijing:
Tel. (8610) 8519-1250
(8610) 8519-1058
Fax (8610) 8519-1255
Guangzhou:
Tel. (8620) 8755-4021
Fax (8620) 8755-4350
Hong Kong:
Tel. (852) 2803-9111
Fax (852) 2513-0313
Shanghai:
Tel. (8621) 5306-9100
Fax (8621) 5306-0838

CZECH REPUBLIC

Tel. 02-2051 3841
02-2051 3842
Fax 02-2051 4298

DENMARK
Tel. 70 10 00 23
Fax 70 10 13 14

FINLAND
Tel. (09) 804 5110
Fax (09) 256 5660

FRANCE
Tel. (01) 30.12.70.00
Fax (01) 30.12.71.80

GERMANY
Tel. (06196) 494-0
Fax (06196) 43901

HUNGARY
Tel. 1-205 9784
Fax 1-205 9792

INDIA
Tel. (91) 80-839 46 57
Fax (91) 80-839 63 45

ITALY
Vimodrone (Milano):
Tel. (02) 25.07.81
Fax (02) 26.50.324
Roma:
Tel. (06) 52.03.600
Fax (06) 52.95.735

JAPAN
Tel. (03) 5442-9711
Fax (03) 5442-9736 Analytical
9737 BioProcess
9734 Lab Water

KOREA

Tel. (822) 551-0990
Fax (822) 551-0228

MALAYSIA
Tel. 603-757 1322
Fax 603-757 1711

MEXICO
Tel. (525) 576-9688
Fax (525) 576-8706

THE NETHERLANDS
Tel. 076-5022000
Fax 076-5022436

NORWAY
Tel. 22 67 82 53
Fax 22 66 04 60

POLAND
Tel. 22-669 12 25
22-663 70 31
Fax 22-663 70 33

PUERTO RICO
Tel. (787) 273-8495
Fax (787) 747-6553

SINGAPORE
Tel. (65) 842 1822
Fax (65) 842 4988

SPAIN AND PORTUGAL
Madrid:
Tel. 917 283 960
Fax 917 292 909
Barcelona:
Tel. 934 525 530
Fax 934 516 048

SWEDEN

Tel. 08-628 6960
Fax 08-628 6457

SWITZERLAND
Volketswil:
Tel. (01) 908-30-60
Fax (01) 908-30-80
Lausanne:
Tel. (021) 641 2550
Fax (021) 641 2551

TAIWAN
Taipei:
Tel. (886-2) 700-1742
Fax (886-2) 755-3267
Hsin Chu City:
Tel. (886-3) 571-0178
Fax (886-3) 572-9520

U.K. AND IRELAND
Tel. (01923) 816375
Fax (01923) 818297

U.S.A.
Tel. (781) 533-6000
Fax (781) 533-3110

IN ALL OTHER COUNTRIES
Millipore Intertech (U.S.A.)
Tel. +1 (781) 533-8622
Fax +1 (781) 533-8630