



QUALITY CONTROL RANGES
MILLIPLEx[®] MAP Human Kidney Injury Magnetic Bead Panel 2
Kit Catalog # HKI2MAG-99K
Control Catalog # HKI2-6099-2
Lot # HKI2-107 and HKI2-207

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
EGF	Control 1	0.07 - 0.14	ng/mL
	Control 2	0.65 - 1.4	ng/mL
Lipocalin-2/NGAL	Control 1	0.13 - 0.27	ng/mL
	Control 2	1.3 - 2.7	ng/mL
Albumin	Control 1	23 - 48	ng/mL
	Control 2	258 - 535	ng/mL
Clusterin	Control 1	31 - 65	ng/mL
	Control 2	315 - 654	ng/mL
Cystatin C	Control 1	0.06 - 0.13	ng/mL
	Control 2	1.0 - 2.2	ng/mL
OPN	Control 1	4.3 - 8.9	ng/mL
	Control 2	42 - 87	ng/mL
α 1-Microglobulin	Control 1	6.2 - 13	ng/mL
	Control 2	64 - 132	ng/mL

**Note: The Quality Control Ranges are generated using assay buffer provided in the kit.
Quality Control Ranges in culture media are not tested.**

HKI2-6099-2: HKI2-107.207 (21 July 2021)
EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com