



## QUALITY CONTROL RANGES

### HUMAN FGF-21 ELISA CATALOG # EZHFGF21-19K

### HUMAN FGF-21 STANDARD, LYOPHILIZED CATALOG# E8019-K

STANDARD	HUMAN FGF-21 VALUE
Lot Number. 819-10K	2000 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5 mL** distilled or deionized water.

### QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #E6019-K

Reconstitute each vial with **0.5 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	FGF21-110	85 – 176	pg/mL
QC 2	FGF21-210	555 – 1154	pg/mL

Cat# E6019-K: FGF21-110.210 (26 February 2021)  
EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX)

• [www.emdmillipore.com](http://www.emdmillipore.com)