

Product Information

α -Casein–Agarose

Saline Suspension

C3905

Storage temperature 2-8 °C

Application

α -Casein is a protein that can be conjugated to agarose for use in affinity chromatography¹⁻³ and protein chromatography.⁴ Casein-agarose may be used for the purification of protein kinases⁵⁻¹⁰ and DNA-dependent RNA polymerases.⁵ This product has been specifically used in various protein binding assays.¹¹⁻¹²

Reagent

Suspension in 0.5 M NaCl containing preservative.

Preparation Instructions

Matrix: 4% beaded agarose activated using cyanogen bromide

Attachment: Amino group of casein

Spacer: 1 atom

Storage/Stability

Do not freeze. Freezing will damage the agarose bead structure.

Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

Procedures

The resin should be washed well with several aliquots of buffer (50 mL per 1 mL gel) to remove all NaCl and preservative, and any free ligand which may have leached from agarose during storage. The resin is stable in the range pH 5.0 to pH 9.0.

One suggested equilibration buffer is as follows: 10 mM Tris-HCl (pH 7.5), 10 mM MgCl₂, 22 mM NH₄Cl, 10 mM β -mercaptoethanol, and 5% glycerol.⁶

Elution of protein may be affected by adding NaCl from 100 mM up to 600 mM to this buffer.³ Removing tightly bound material (because of nonspecific binding) may require using at least 10 column volumes of 1 M NaCl for each column volume of resin. If material still remains, a 2 M urea solution (at least 10 column volumes) may be successful.

References

1. Hutadilok-Towatana, N. *et al.*, *J. Biosci. Bioeng.*, **87(5)**, 581-587 (1999).
2. Kocabiyik, S., and Özdemir, İ., *Biosci. Biotechnol. Biochem.*, **70(1)**, 126-134 (2006).
3. Yang, H. *et al.*, *J. Bacteriol.*, **189(7)**, 2599-2609 (2007).
4. Vijayaraghavan, P., and Vincent, S.G.P., *Biomed. Res. Int.*, **2014**, 725064 (2014).
5. Rahmsdorf, H.J. and Pai, S.-H., *Biochim. Biophys. Acta*, **267(2)**, 339-346 (1979).
6. *Affinity Chromatography: A Practical Approach*, P.D.G. Dean, W.S. Johnson, and F.A. Middle (eds.). IRL Press (Oxford, UK), pp. 125-131 (1985).
7. Zioncheck, T.F. *et al.*, *J. Biol. Chem.*, **261(33)**, 15637-15643 (1986).
8. Sanghera, J.S. *et al.*, *Biochem. J.*, **283**, 829-837 (1992).
9. Pan, Z.-Q., and Hurwitz, J., *J. Biol. Chem.*, **268(27)**, 20433-20442 (1993).
10. Biji, G.D. *et al.*, *3 Biotech*, **6(2)**, 231 (2016).
11. Politis, I. *et al.*, *J. Dairy Sci.*, **78(3)**, 484-490 (1995).
12. Olivares, J.E. *et al.*, *BMC Plant Biol.*, **11**, 134 (2011).

Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

Technical Assistance

Visit the tech service page at SigmaAldrich.com/techservice.

Standard Warranty

The applicable warranty for the products listed in this publication may be found at SigmaAldrich.com/terms.

Contact Information

For the location of the office nearest you, go to SigmaAldrich.com/offices.

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

C3905 Rev 01/22