



Product Information

Fas LIGAND/TNFSF6

Human, Recombinant
Expressed in CHO cells

Product Number **F 0427**

Synonyms: FasL, CD95 ligand, Apo-1 ligand, TNFSF6

Product Description

Recombinant human Fas ligand (FasL) consists of amino acid residues 134-281 of human FasL¹ fused at the amino terminus to a 6X histidine tag which is then linked at the carboxy-terminus to a CD33 (Met1 – Met17) signal peptide. The N-terminal 6X his-tagged recombinant protein is expressed in a CHO cells. Recombinant human Fas ligand, after removal of the CD33 signal peptide, is a non-covalently linked trimer. The monomer has a calculated molecular weight of approximately 18 kDa. As a result of glycosylation, recombinant Fas Ligand migrates as an approximately 26 to 28 kDa protein on SDS-PAGE under reducing and non-reducing conditions.

Native Fas ligand is a 40 kDa, type II membrane protein belonging to the TNF super family. Its specific receptor, Fas (CD95/Apo-1), is a 45 kDa protein that transduces the apoptotic signal into cells. Fas Ligand is expressed primarily by T cells and NK cells and to a lesser degree in the testes, cornea, and several malignant tumor cell types. Human Fas ligand shares 85% sequence identity with mouse Fas ligand and 78% sequence identity with rat Fas ligand.

The Fas/FasL system modulates the immune response by inducing cell apoptosis to maintain homeostasis and in the regulation of immune responses and privilege.² Cysteine-rich repeats of the Fas receptor are required for binding by the Fas ligand.³ This binding induces trimerization of Fas in the target cell membrane and activation. Activation of Fas causes the recruitment of Fas-associated protein with death domain (FADD) via interactions between the death domains of Fas and FADD.⁴ Pro-caspase 8 binds to Fas-bound FADD via interactions between the death effector domains (DED) of FADD and pro-caspase 8 leading to the activation of caspase 8. Activated caspase 8 then cleaves (activates) nine other procaspases, a process that ultimately leads to apoptosis.

B cell antigen receptor signaling inhibits Fas-mediated apoptosis via up-regulation of cellular FLICE-inhibitory protein (c-FLIP).⁵ The expression of the Fas ligand is

regulated by protein phosphatases(s) sensitive to okadaic acid.⁶ Serum withdrawal-induced apoptosis is mediated partially by the Fas/FasL interactions.⁷ Studies suggest the Fas ligand is also a potent chemotactic factor in polymorphonuclear neutrophils which may be independent of the death-domain mediated apoptosis previously described.⁸

Reagent

Recombinant human Fas ligand is supplied as approximately 10 µg of protein lyophilized from a sterile filtered phosphate-buffered saline (PBS) solution containing 50 µg bovine serum albumin per µg of cytokine.

Preparation Instructions

Reconstitute the vial contents with sterile PBS containing a minimum of 1% human or bovine serum albumin. Stock solution concentration should be no less than 10 µg/ml.

Storage/Stability

Lyophilized samples are stable for at least six months at –20 °C. Upon reconstitution, store at 2-4 °C for up to one month. For extended storage, store in working aliquots at –20 °C. Repeated freeze-thaw cycles should be avoided. Do not store in frost-free freezer.

Product Profile

Recombinant Fas Ligand activity is measured by its ability to induce apoptosis in Jurkat cells. Recombinant human Fas ligand effective range is 0.2 to 0.5 µg/ml in the presence of a cross-linking anti-6X histidine antibody. Optimal dilutions should be determined by each laboratory for each application.

Purity: >95% by SDS-PAGE, visualized by silver stain.

Endotoxin level: < 0.1 ng/µg of protein as determined by the LAL (Limulus amoebocyte lysate) method.

References

1. Mita, E., et al., Role of Fas ligand in apoptosis induced by hepatitis C virus infection. *Biochem. Biophys. Res. Commun.*, **204**, 468-474 (1994).
2. Starling, G. c., et al., Identification of amino acid residues important for ligand binding to Fas. *J. Exp. Med.*, **185**, 1487-1492 (1997).
3. Orlinick, J. R., et al., Requirement of cysteine-rich repeats of the Fas receptor for binding by the Fas ligand. *J. Biol. Chem.*, **272**, 28889-28894 (1997).
4. Juo, P., et al., FADD is required for multiple signaling events downstream of the receptor Fas. *J. Cell Growth Differ.*, **10**, 797-804 (1999).
5. Wang, J., et al., Inhibition of Fas-mediated apoptosis by the B cell antigen receptor through c-FLIP. *Eur. J. Immunol.*, **30**, 155-163 (2000).
6. Morimoto, Y., et al., Upregulation of the expression of Fas antigen and Fas ligand in a human submandibular gland ductal cell line by okadaic acid. *Arch. Oral Biol.*, **45**, 657-666 (2000).
7. Hu, C. L., Apoptosis of bovine granulosa cells after serum withdrawal is mediated by Fas antigen (CD95) and Fas ligand. *Biol. Reprod.*, **64**, 518-526 (2001).
8. Seino, K., et al., Cutting Edge: Chemotactic Activity of Soluble Fas Ligand Against Phagocytes. *J. Immunol.*, **161**, 4484-4488 (1998).

lcm 2/03

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.